# mectron

medical technology

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USE AND MAINTENANCE MANUAL

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IT


EN


ES

# combi touch

CE  
0051

mectron

medical technology

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IT

# combi touch

---

**Copyright**

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---

**Copyright**

© Mectron S.p.A. 2021. All rights reserved. No part of this document may be reproduced, in any form whatsoever, without the prior written consent of the copyright holder.

## CONTENTS

|          |   |           |
|----------|---|-----------|
| <b>1</b> | <b>Introduction</b>                                       | <b>1</b>  |
| 1.1      | Intended Use  | 1         |
| 1.1.1    | Piezoelectric Ultrasonic Scaler                           | 2         |
| 1.1.2    | Polisher for Dental Prophylaxis and Hygiene               | 2         |
| 1.2      | Description of the Appliance                              | 2         |
| 1.2.1    | Eligible Patient Group                                    | 3         |
| 1.2.2    | Patient Selection Criteria                                | 3         |
| 1.2.3    | Indications for Use                                       | 3         |
| 1.2.4    | Users   | 3         |
| 1.3      | Disclaimer  | 4         |
| 1.4      | Safety Requirements                                       | 5         |
| 1.5      | Symbols   | 8         |
| <b>2</b> | <b>Identification Data</b>                                | <b>10</b> |
| 2.1      | Appliance ID Plate  | 10        |
| 2.2      | Scaler Handpiece Identification Data                      | 10        |
| 2.3      | Inserts Identification Data                               | 11        |
| 2.4      | Air-Polishing Handpiece Identification Data               | 11        |
| <b>3</b> | <b>Delivery</b>   | <b>12</b> |
| 3.1      | List of components  | 12        |
| <b>4</b> | <b>Installation</b>                                       | <b>17</b> |
| 4.1      | First Installation  | 17        |
| 4.2      | Safety Requirements During Installation                   | 17        |
| 4.3      | Connecting the Accessories                                | 18        |
| <b>5</b> | <b>Use</b>  | <b>23</b> |
| 5.1      | Power On/Off  | 23        |
| 5.2      | Description of the Keyboard                               | 24        |
| 5.2.1    | ULTRASOUND - Scaler Part                                  | 24        |
| 5.2.2    | POLISHING - Polisher Part                                 | 28        |
| 5.2.3    | Symbols   | 29        |
| 5.3      | Safety Requirements Before and During Use                 | 29        |
| 5.3.1    | Ultrasonic Scaler   | 30        |
| 5.3.2    | Jet Polisher  | 32        |
| 5.4      | Instructions for Use - Scaler Part                        | 32        |
| 5.5      | Important information on the Inserts                      | 35        |
| 5.6      | Instructions for Use - Polisher Part                      | 36        |
| 5.6.1    | "Refill" function   | 38        |
| <b>6</b> | <b>Flush function</b>                                     | <b>40</b> |
| 6.1      | ULTRASOUND - Scaler Part                                  | 40        |
| 6.2      | POLISHING - Polisher Part                                 | 43        |
| <b>7</b> | <b>Disassembling Parts for Cleaning and Sterilisation</b> | <b>45</b> |
| <b>8</b> | <b>Cleaning and Sterilisation</b>                         | <b>49</b> |
| 8.1      | Preparation   | 50        |
| 8.1.1    | Disinfecting the External Irrigation Circuit              | 51        |
| 8.1.1.1  | Required Materials  | 51        |
| 8.1.1.2  | Procedure - Polisher Handpiece                            | 51        |
| 8.1.1.3  | Procedure - Scaler Handpiece                              | 60        |
| 8.1.2    | Disinfecting the Tank Irrigation Circuit                  | 68        |
| 8.1.2.1  | Required Materials  | 68        |
| 8.1.2.2  | Procedure   | 69        |
| 8.2      | Cleaning the Non-Sterilisable Parts                       | 74        |

|           |   |            |
|-----------|---|------------|
| 8.2.1     | Required Materials  | 74         |
| 8.2.2     | Cleaning Method   | 74         |
| 8.3       | Cleaning the Tank and Cap   | 75         |
| 8.3.1     | Preparation   | 75         |
| 8.3.2     | Required Materials  | 75         |
| 8.3.3     | Cleaning Method   | 75         |
| 8.4       | Cleaning the Sterilisable Accessories   | 77         |
| 8.4.1     | Manual Cleaning   | 77         |
| 8.4.1.1   | Required Materials  | 77         |
| 8.4.1.2   | Scaler Handpiece  | 78         |
| 8.4.1.3   | Inserts   | 82         |
| 8.4.1.4   | Inserts Tightening Wrench   | 84         |
| 8.4.1.5   | Air-polishing handpieces  | 86         |
| 8.4.1.6   | K9 Wrench   | 89         |
| 8.4.2     | Automatic Cleaning  | 91         |
| 8.4.2.1   | Required Materials  | 91         |
| 8.5       | Cleaning Check  | 92         |
| 8.5.1     | Required Materials  | 92         |
| 8.6       | Drying and Lubrication  | 94         |
| 8.6.1     | Required Materials  | 94         |
| 8.7       | Sterilisation   | 96         |
| 8.7.1     | Preparation   | 96         |
| 8.7.2     | Sterilisation Method  | 97         |
| 8.7.2.1   | Special Information   | 97         |
| <b>9</b>  | <b>Maintenance</b>  | <b>98</b>  |
| 9.1       | Replacing the Peristaltic Pump  | 98         |
| 9.2       | Cleaning and/or Replacing the Water Filter  | 100        |
| 9.3       | Eliminating the Condensate  | 101        |
| 9.4       | Cleaning the Powder Containers and Caps   | 101        |
| <b>10</b> | <b>Methods and Precautions for Disposal</b>   | <b>102</b> |
| <b>11</b> | <b>Technical Data</b>   | <b>103</b> |
| 11.1      | Ultrasound part   | 104        |
| 11.2      | Polishing part  | 104        |
| 11.3      | Electromagnetic Compatibility IEC/EN 60601-1-2  | 105        |
| 11.3.1    | Guide and Manufacturer's Declaration - Electromagnetic Emissions  | 105        |
| 11.3.2    | Accessible Parts of the Casing  | 106        |
| 11.3.3    | Guide and the Manufacturer's Declaration - Electromagnetic Immunity   | 107        |
| 11.3.3.1  | Power Connection BC Input   | 107        |
| 11.3.3.2  | Points of Contact with the Patient  | 109        |
| 11.3.3.3  | Parts Accessible to the Input / Output Signals  | 110        |
| 11.3.4    | Specifications of the tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device | 111        |
| <b>12</b> | <b>Troubleshooting</b>  | <b>112</b> |
| 12.1      | Diagnostic System and Symbols on the Keyboard   | 112        |
| 12.2      | Quick Troubleshooting   | 114        |
| 12.3      | Replacing the Fuses   | 119        |
| 12.4      | Shipping to an Authorised Mectron Service Centre  | 120        |
| <b>13</b> | <b>Warranty</b>   | <b>121</b> |

**PAGE LEFT INTENTIONALLY BLANK**

# 1 INTRODUCTION

Carefully read this manual before proceeding with installation operations, use, maintenance or other intervention operations on the appliance.

This manual must always be available to the operator.

**Important:** To avoid damages to persons or property, carefully read all the "Safety Requirements" in the manual.

Depending on the level of severity, the safety requirements are classified with the following indications:

**⚠ WARNING:** (Always referred to damage to persons)

**ⓘ CAUTION:** (Referring to possible damage to property)

The aim of this manual is to make the operators aware of the safety regulations,

installation procedures, instructions for proper use and maintenance of the appliance and its accessories.

Use of this manual for aims other than those strictly linked to the installation, use and maintenance of the appliance is prohibited.

The information and illustrations in this manual were updated on the edition date shown on the last page.

MECTRON is engaged in continuously updating its products with possible changes to the components of the apparatus.

In case you encounter discrepancies between the descriptions found in this manual and the equipment in your possession you can:

- check for any available updates in the *MANUALS* section of the *MECTRON* website<sup>1</sup>;
- ask clarifications to Your Dealer;
- contact MECTRON After Sales Service.

## 1.1 Intended Use

combi touch combines, in a single appliance, a multi-purpose piezoelectric scaler and water jet, air, and powder polisher, distributed by Mectron, intended for complete supra- and subgingival dental prophylaxis.

The Mectron air-polishing devices and piezoelectric ultrasonic scalers can be used on patients of any age and gender requiring a tooth polishing treatment. There are no contraindications for specific population segments.

**⚠ WARNING:** The appliance must be used in a dental practice or clinic, or in professional oral hygiene and preventive care centres. Do not use the device in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

**⚠ WARNING: Qualified and specialist personnel.** The appliance should only be used by qualified personnel with adequate medical culture; for the use of the device are not requested special training activities. Use of the device does not produce any side effects if used correctly. Improper use is manifested by the transfer of heat to the tissues.

**⚠ WARNING: Intended use.** Use the appliance exclusively for the use for which it is intended. Failure to observe this regulation can cause serious injury to the patient, the operator, and damage/failure of the device.

<sup>1</sup> <https://manuals.mectron.com/>

## 1.1.1 Piezoelectric Ultrasonic Scaler

With the appropriate inserts it is possible to perform the following treatments:

- **Scaling:** all procedures for the removal of bacterial plaque deposits and supragingival, subgingival, and interdental stones, and removal of stains;
- **Periodontology:** periodontal root planing and debridement therapy, including cleaning and irrigation of the periodontal pocket;
- **Cleaning treatment of the implant surface;**
- **Endodontics:** all treatments for root canal preparation, irrigation, filling, gutta-percha condensation, endodontic reprocessing and retrograde cavity preparation;
- **Restorative and prosthetics:** cavity preparation and removal of carious tissue, removal of prosthetics and excess restoration materials, amalgam condensation, finishing of prosthetic stump;
- **Extraction techniques.**

## 1.1.2 Polisher for Dental Prophylaxis and Hygiene

combi touch features a water jet, air, and powder polisher for dental use, designed for a complete supra- and subgingival prophylaxis treatment, with a specific intended use depending on the type of powder distributed by Mectron.

Air-polisher for supragingival use:

- Removal of supragingival bacterial plaque;
- Removal of stains on the surface of teeth;
- Preparation of cavity for better adhesion between enamel and filling material;
- Prophylaxis on patients undergoing orthodontic therapy.

Air-polisher for subgingival use:

- Removal of subgingival bacterial plaque;
- Removal of Biofilm in the prevention of Peri-implantitis.

**⚠ WARNING:** In the presence of composite resin surfaces the jet must be directed towards the part to be treated, on average 2-3 seconds per tooth.

## 1.2 Description of the Appliance

combi touch combines, in a single appliance, a multi-purpose piezoelectric scaler and water jet, air, and powder polisher, intended for complete supra- and subgingival dental prophylaxis.

In regards to the various possible ultrasonic treatments, combi touch can be used either connected directly to the main water supply in the dental practice, or with independent irrigation through the special liquid container, which can hold different types of medical solutions.

The appliance is equipped with an automatic tuning circuit which compensates the wear of the inserts, thus always allowing operation in conditions of maximum efficiency.

The operating principle of the polisher is based on the mechanical action obtained from

a jet of various crystal types accelerated by a flow of compressed air. The kinetic energy thus imprinted at the particles, dissipates almost completely due to impact against the surface of the enamel, producing a gentle but effective cleaning action. The action is completed by a jet of water which, using the vacuum created around the nozzle, has a bell-shape around the main flow, thus producing a double effect: to prevent much of the rebound and the leakage of the cloud of powder and perform continuous washing of the treated area, dissolving the powder.

### 1.2.1 Eligible Patient Group

This medical device is designed to be used with the following patient groups:

- Children;
- Adolescents;
- Adults;
- Seniors.

This medical device can be used on patients of any age, weight, height, gender, and nationality, where applicable.

### 1.2.2 Patient Selection Criteria

It is not recommended to use the device in the following cases.

Air-polisher for supra- and subgingival use:

1. Upper respiratory tract infections, chronic bronchitis/asthma;
2. Pregnant and breastfeeding women;
3. Patients undergoing treatment (radiotherapy, chemotherapy, antibiotics);
4. Acute infectious oral lesions.

Piezoelectric ultrasonic scaler:

1. Patients with active implantable medical devices (for example: pacemakers, hearing aids, and/or other electromagnetic prostheses) without prior authorisation from their doctor;
2. Patients with clinical conditions not suitable for treatment of the sites (for example: local anaesthesia).

It is not recommended to use the powders in the following cases:

1. Allergy to the aroma of the powder;
2. Patients on a sodium restricted diet or who suffer from serious respiratory problems, such as chronic bronchitis, asthma, emphysema, etc., unless

otherwise specified by their doctor.

All air-polishing devices and piezoelectric ultrasonic scalers are intended for professional use only. The user is therefore the only person able to decide if and how to treat their patients.

### 1.2.3 Indications for Use

Use of the device is suitable for all eligible patients (see *Chapter 1.2.1 on page 3*) whose doctor has prescribed a treatment among those described in the intended use of the device (see *Chapter 1.1.1 on page 2* for the function of the piezoelectric ultrasonic scaler and *Chapter 1.1.2 on page 2* for the function of the air-polisher).

### 1.2.4 Users

The device must be used exclusively by specialised and properly trained personnel, specifically the dental surgeon/dentist or dental hygienist, who must be able-bodied, adult, of any weight, age, height, gender, and nationality.

### 1.3 Disclaimer

The manufacturer MECTRON disclaims all liability, express or implied, and cannot be held responsible for direct or indirect personal injury and/or property damage, occurring as a result of incorrect procedures linked to the use of the appliance and its accessories.

The manufacturer MECTRON cannot be held liable, expressly or implicitly, for any type of personal injury and/or property damage, inflicted by the user of the product and its accessories and which occurs in the following cases, by way of example and not of limitation:

- Misuse or use during procedures other than those specified in the destination of use of the product;
- The environmental conditions for preservation and storage of the device are not complying with the requirements indicated in *Chapter 11 on page 103*;
- The appliance is not used in accordance with all the instructions and requirements described in this manual;
- The electrical system of the places where the equipment is used do not comply with the laws in force and the related regulations;
- Assembly operations, extensions, re-adjustments, upgrades and repairs of the device are carried out by personnel not authorised by MECTRON;
- Misuse, abuse, abnormal use, negligent use, intentional misconduct or use exceeding the limits of the appliance indicated and allowed and/or normal wear or deterioration, ill-treatment and/or incorrect interventions;
- Any attempt to tamper with or modification of the appliance under every circumstance;
- Use of non original MECTRON inserts resulting in permanent damage to the thread of the handpiece with impaired functioning and risk of injury to the patient;
- Use of non original MECTRON inserts used in accordance with the settings designed and tested on the original MECTRON inserts. Correct use of the settings is only guaranteed with original MECTRON inserts;
- Shortage of stock material (handpiece, inserts, wrenches) to be used in case of a standstill due to faults or problems;
- Incorrect/omitted maintenance compared to what is stated in *Chapter 9 on page 98* of this manual;
- Breach of the requirements and the information contained in *Chapter 5.5 on page 35* of this manual;
- Breach of the requirements and the information contained in *Chapter 8 on page 49* of this manual;
- Unauthorised repairs in accordance with the indications contained in *Chapter 12.4 on page 120* this manual.

## 1.4 Safety Requirements

### **⚠ WARNING: Contraindications.**

Do not use combi touch on patients with Pacemakers or other implantable electronic devices. This regulation also stands for the operator.

**⚠ WARNING: Contraindications.** Do not perform scaling treatments without water spray in order to avoid an overheating of the insert that can cause damage to the tooth. Treatments without water spray can only be those performed with the "Dry Work" inserts without water passage.

**⚠ CAUTION: Contraindications. Ultrasonic scaler.** Do not perform treatments on metal or ceramic prosthetics. Ultrasonic vibrations could lead to decementation of the dentures.

### **⚠ WARNING: Contraindications. Interference from other equipment.**

An electrosurgical scalpel or other electrosurgical units positioned near the combi touch appliance may interfere with the correct operation of the appliance itself.

**⚠ WARNING: Contraindications. Interference with other equipment.** Even if compliant with standard IEC 60601-1-2, combi touch may interfere with other devices in the vicinity. combi touch must not be used near or stacked with other appliances. However, if this should become necessary, it is necessary to check and monitor the correct functioning of the appliance in that configuration.

**⚠ WARNING: Risk of explosion.** The appliance cannot operate in environments where the atmosphere is saturated with flammable gases (aesthetic mixtures, oxygen, etc.).

**⚠ CAUTION:** In case the final user, operating in their own medical room or surgery, in order to comply with mandatory requirements, must periodically inspect the equipment present in the surgery, the test procedures to apply to medical electrical equipment and medical electrical systems for the safety assessment must be carried out following the standard EN 62353 'Medical electrical equipment - Periodic inspections and tests to be carried out after repair of medical electrical equipment'. The frequency of periodic inspections in the intended conditions of use described in this "Use and Maintenance" manual is once per year or every 2000 hours of use, whichever condition is satisfied first.

### **⚠ WARNING: Checking the condition of the device before treatment.**

Always check that there is no water under the appliance. Before every treatment always check the perfect operation of the appliance and efficiency of the accessories. If operation anomalies are observed do not carry out the treatment. If the faults concern the appliance, contact an Authorised MECTRON Service Centre.

**⚠ CAUTION:** The electrical system of the premises in which the appliance is installed and used must comply with the rules in force and the relevant requirements of electrical safety.

**⚠ CAUTION:** To avoid the risk of electrical shocks, this appliance has to be connected exclusively to a power supply with an earth safety device.

**⚠ WARNING: Cleaning and sterilising new or repaired instruments.** All the new or repaired appliance accessories are not sterile. At first use and after each treatment they must be cleaned and sterilised by carefully following the instructions in *Chapter 8* on page 49.

## **⚠ WARNING: Infection control.**

For maximum patient and operator safety, before using all the parts and reusable accessories, make sure they have been previously cleaned and sterilised by following the instructions in *Chapter 8 on page 49*.

**⚠ CAUTION: Contraindications.** After having autoclaved the handpiece, inserts, torque wrench or any other sterilisable accessory, wait for them to completely cool down before reusing them.

## **⚠ WARNING: Insert Breakage and Wear.**

High frequency oscillations and wear can, in rare cases, lead to the breakage of the insert. Deformed inserts or inserts that have been otherwise damaged are prone to breakage during use. Broken or worn inserts must not ever be used.

In the event of breakage check that no fragments remain in the treated part and at the same time use suction effectively to remove them.

The patient must be instructed to breathe through the nose during the treatment, or use a dental dam, so as to avoid ingesting fragments of broken inserts.

Check the state of wear of the insert and its integrity before and during each use. If a drop in performance occurs, see to its replacement.

The state of wear of the most common inserts (S1, S1-S, S2, S5, P2, P4, P10) can be checked using the supplied INSERT-CARD. To use the INSERT-CARD correctly:

- Position the insert on the INSERT-CARD so that the profile corresponds to the one printed on the card. The printed profile on the card has a red line indicating the limit of wear;
- If the insert is shorter than the limit of wear, its performance will be significantly inferior compared to that of a new insert, and should therefore be replaced.

If the layer of titanium nitride (gold-plated surface), where present, is visibly worn, the insert must be replaced. The use of a worn insert reduces its efficiency.

**Diamond Inserts:** diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.

When the nitriding is consumed, the cutting edge loses effectiveness; a possible regrinding harms the insert, therefore it is prohibited. Verify that the insert is not worn. During the operation frequently check that the insert is intact, especially in the apical part.

During the operation avoid prolonged contact with the retractor or with metal instruments in use. Do not exert excessive pressure on the inserts during use.

**⚠ WARNING:** Exclusively use original MECTRON inserts, accessories and spare parts.

**⚠ CAUTION:** It is prohibited to make any changes to this appliance.

**⚠ WARNING: Contraindications - Jet polisher.** Patients wearing contact lenses must remove them prior to receiving treatment with the jet polisher.

**⚠ WARNING: Contraindications - Supragingival powder jet polisher.** Do not aim the air, supragingival powder, or water jet onto the soft tissues or inside the gingival sulcus. Failure to comply with this prescription can cause a gingival tissue emphysema (emphysema of the mucous and/or subcutaneous). For this type of application, use only subgingival powder.

**⚠ WARNING: Contraindications - Jet polisher.** Do not use the device near areas subject to recent dental extraction and in traumatised/damaged areas (or areas nearby) due to the risk of emphysema.

 **WARNING: Temperature of the water spray - Jet polisher.** The appliance is equipped with a double safety device that controls the temperature of the water spray. Before treatment, it is however recommended to instruct the patient to inform the operator if he perceives an excessive increase in the temperature of the water.

 **WARNING:** In case of an adverse event and/or accident attributable to the device during correct use and in accordance with the intended use, a report must be made to the Competent Authority and manufacturer indicated on the product label.

## 1.5 Symbols

| Symbol  | Description   | Symbol  | Description   |
|---|---|---|---|
|    | Device compliant with Regulation (EU) 2017/745. Notified body: IMQ S.p.A.                         |    | Nemko brand Compliance with UL - CSA regulations      |
|    | Medical device  |    | Caution! Read the instructions for use                |
|    | Operating instructions  |    | Manufacturer  |
|    | Date of manufacture   |    | Serial number   |
|    | Batch number  |    | Product code  |
|    | Disposable  |    | Expiry date   |
|    | Sterilised with Ethylene Oxide (EO)   |    | Not sterile   |
|    | The sterilisable materials must be autoclaved and can resist up to a maximum temperature of 135°C |    | "B" Part applied according to regulation EN 60601-1   |
|    | Alternating current   |   | Generic warning signal <sup>a)</sup>                  |
|  | Power switch set to "on"  |  | Power switch set to "off"                             |
|  | Connection of the control pedal   |  | Biological Risk                                       |
|  | The appliance and its accessories must not be disposed of or treated as municipal solid waste     |  | Temperature limits for transport and storage          |
|  | Moisture limits for transport and storage   |  | Atmospheric pressure limits for transport and storage |

| Symbol  | Description   | Symbol  | Description                          |
|---------|---|---|--------------------------------------|
| QTY.1   | Quantity in package: 1  |  | Do not use if the package is damaged |
| Rx Only | For US market only<br><b>CAUTION:</b> US federal law restricts sale to orders by licensed dentists or dental hygienists only. |   |                                      |

**Table 1** – Symbols

- a) The symbol is represented by a yellow triangle and black graphic symbol.

**NOTE:** For other symbols, refer to *Chapter 12.1* on page 112.

## 2 IDENTIFICATION DATA


A correct description of the model and of the serial number of the appliance will allow the After Sales Service to provide fast and effective answers.

Always provide this information every time that you contact MECTRON After Sales Service.

### 2.1 Appliance ID Plate


Each appliance is equipped with an identification plate indicating the main technical characteristics and serial number. The identification plate is located under the appliance. The complete technical specifications are provided in Chapter 11 on page 103.

**NOTE:** The complete list of symbols is provided in Chapter 1.5 on page 8.

Additional symbols and characteristics of the device are provided on a separate plate. This identification plate is located under the appliance.


**NOTE:** The complete list of symbols and their description is provided in Chapter 1.5 on page 8

### 2.2 Scaler Handpiece Identification Data

The handpiece name, MECTRON logo, serial number, product code, HIBC matrix code, and a series of symbols are laser engraved on the scaler handpiece.


**NOTE:** The complete list of symbols is provided in Chapter 1.5 on page 8.

EN

## 2.3 Inserts Identification Data

The name of the insert itself (Reference 1), MECTRON logo (Reference 2), and insert batch number (Reference 3) are laser engraved on each insert.


EN

## 2.4 Air-Polishing Handpiece Identification Data

The HIBC matrix data (Reference 4), product code and handpiece batch (Reference 5), handpiece name (Reference 6), and MECTRON logo (Reference 7) are laser engraved on each air-polishing handpiece.

**NOTE:** The complete list of symbols is provided in *Chapter 1.5* on page 8.

## 3 DELIVERY

### 3.1 List of components

See the Figure on the inside cover. combi touch is supplied with the standard equipment (see Table 2 on page 12), a set of variable accessories depending on the configuration and client requests (see Table

3 on page 15), and accessories that can be ordered separately (see Table 4 on page 16).

**NOTE:** Both the items included in the standard supply and all accessories can be ordered separately by the client.

| Standard supply                        |              |  |      |
|--|--------------|--|------|
| Item                                   | Code         | Description  | Ref. |
| Device core unit                       | 04000107     | b)   | A    |
| Irrigation tank complete with cap      | 03150102     | b)   | B    |
| Pedal FS-05 with cord and plug         | 02900114     | b)   | C    |
| Tank safety cap                        | 01950014     | Grey safety cap <sup>b)</sup>  | D    |
| Tank cap                               | 03020082     | Grey tank cap <sup>b)</sup>  | E    |
| Water supply hose with quick coupling  | 02900119     | b)   | F    |
| Air supply hose with quick coupling    | 02900012     | b)   | G    |
| Water filter                           | 00420004     | b)   | H    |
| combi touch O-ring kit                 | 03020206-001 | b)   | I    |
| Cleaning needle kit                    | 02900151-001 | Includes: <ul style="list-style-type: none"> <li>• Cleaning needle Ø 0.4 mm (cod. 02900123);</li> <li>• Cleaning needle Ø 0.8 mm (cod. 02900148). <sup>b)</sup></li> </ul> | J    |
| Scaler Handpiece                       | 03120238     | MANIPOLO ABLATORE SLIM (SLIM SCALER HANDPIECE) <sup>b)</sup>   | K2   |
| Mectron water circuit disinfection kit | 02900171     | b)   | W    |

**Table 2** – Standard supply.

| Accessories available for order with standard supply               |          |   |      |
|--|----------|---|------|
| Item   | Code     | Description   | Ref. |
| MANIPOLO AIR-POLISHING<br>90° (90° AIR-POLISHING<br>HANDPIECE)     | 03120239 | 90° air-polishing handpiece <sup>b)</sup>                         | M1   |
| MANIPOLO AIR-POLISHING<br>120° (120° AIR-POLISHING<br>HANDPIECE)   | 03120240 | 120° air-polishing handpiece <sup>b)</sup>                        | M2   |
| MANIPOLO AIR-POLISHING<br>PERIO (PERIO AIR-POLISHING<br>HANDPIECE) | 03120241 | PERIO air-polishing handpiece <sup>b)</sup>                       | M3   |
| SUBGINGIVAL PERIO TIPS   | 02900120 | 40 pc set of SUBGINGIVAL PERIO<br>TIPS <sup>b)</sup>              | N    |
| K9 Wrench  | 02900122 | <sup>b)</sup>   | O    |
| Jar of specific subgingival<br>powder                              | 03140022 | Prophylaxis Powder: SENSITIVE +<br>(single package) <sup>c)</sup> | P    |
|  | 03140021 | Prophylaxis Powder: SENSITIVE +<br>(2x160g package) <sup>c)</sup> |      |
| Jar of specific supragingival<br>powder                            | 03140028 | Prophylaxis Powder: SOFT +<br>(single package) <sup>c)</sup>      | Q    |
|  | 03140027 | Prophylaxis Powder: SOFT +<br>(4x250g package) <sup>c)</sup>      |      |
|  | 03140009 | Prophylaxis Powder: SMOOTH<br>(single package) <sup>b)</sup>      |      |
|  | 03140010 | Prophylaxis Powder: SMOOTH<br>(4x250g package) <sup>b)</sup>      |      |
| Torque wrench  | 02900137 | K10 torque wrench <sup>b)</sup>                                   | R    |
|  | 02900074 | K6 torque wrench <sup>d) b)</sup>                                 |      |
|  | 02900081 | K7 torque wrench <sup>e) b)</sup>                                 |      |
| Inserts  | 0296xxxx | "S" series reusable scaler inserts<br><sup>b)</sup>               | S    |
|  | 0308xxxx | "PE" series reusable scaler inserts<br><sup>b)</sup>              |      |
|  | 0305xxxx | "R" series reusable scaler inserts<br><sup>b)</sup>               |      |
|  | 0345xxxx | "ER" series reusable scaler inserts<br><sup>b)</sup>              |      |

| Accessories available for order with standard supply |          |   |      |
|--|----------|---|------|
| Item   | Code     | Description   | Ref. |
|  | 0235xxxx | "E" series reusable scaler inserts <sup>b)</sup>    |      |
|  | 0299xxxx | "D" series reusable scaler inserts <sup>b)</sup>    |      |
|  | 0219xxxx | "CM" series reusable scaler inserts <sup>b)</sup>   |      |
|  | 0313xxxx | "ME" series reusable scaler inserts <sup>b)</sup>   |      |
|  | 03570004 | "ICS" series reusable base insert <sup>b)</sup>     |      |
|  | 03590009 | "IC1" tip <sup>b)</sup>                             |      |
|  | 02900112 | IC1 tips kit (5 pcs.) <sup>b)</sup>                 |      |
|  | 03570001 | DB1 reusable base support for inserts <sup>b)</sup> |      |
|  | 03590001 | TA12D60 diamond insert <sup>b)</sup>                |      |
|  | 03590004 | TA12D90 diamond insert <sup>b)</sup>                |      |
|  | 03590002 | TA14D60 diamond insert <sup>b)</sup>                |      |
|  | 03590005 | TA14D90 diamond insert <sup>b)</sup>                |      |
|  | 03590007 | TA14D120 diamond insert <sup>b)</sup>               |      |
|  | 03590003 | TA16D60 diamond insert <sup>b)</sup>                |      |
|  | 03590006 | TA16D90 diamond insert <sup>b)</sup>                |      |
|  | 03590008 | TA16D120 diamond insert <sup>b)</sup>               |      |
|  | 03590010 | TF12D60 diamond insert <sup>b)</sup>                |      |
|  | 03590011 | TF16D60 diamond insert <sup>b)</sup>                |      |
|  | 03590012 | TF12D90 diamond insert <sup>b)</sup>                |      |

| Accessories available for order with standard supply |          |  |      |
|--|----------|--|------|
| Item   | Code     | Description  | Ref. |
|  | 03590013 | TF16D90 diamond insert <sup>b)</sup>                   |      |
| Use and Maintenance Manual                           | 02150313 | IT/EN/ES version <sup>b)</sup>                         | T    |
|  | 02150322 | DE/FR/SV version <sup>b)</sup>                         |      |
| Information: Documentation                           | Online   | 02150650   |      |
| Power supply cable                                   | 00050030 | Black power supply cable 3X0.75 2m IT <sup>c)</sup>    | X    |
|  | 00050031 | Black power supply cable 3X0.75 2m D <sup>c)</sup>     |      |
|  | 00050034 | Black power supply cable 3X0.75 2m UK <sup>c)</sup>    |      |
|  | 00050010 | Black power supply cable 3X1 2.5m AU <sup>c)</sup>     |      |
|  | 00050032 | Black power supply cable 3X0.75 2m CH <sup>c)</sup>    |      |
|  | 00050029 | Black power supply cable 3X0.75 HG 2m JP <sup>c)</sup> |      |
|  | 00050035 | Black power supply cable 3X0.75 2m BR <sup>c)</sup>    |      |

**Table 3** – Accessories available for order with standard supply.

| Accessories available for separate order |          |   |      |
|--|----------|---|------|
| Item                                     | Code     | Description   | Ref. |
| Tank safety cap                          | 01950015 | Blue safety cap <sup>b)</sup>   | D    |
| Tank cap                                 | 03020081 | Blue tank cap <sup>b)</sup>   | E    |
| Scaler Handpiece                         | 03120224 | MANIPOLO ABLATORE LED COMBI TOUCH (COMBI TOUCH LED SCALER HANDPIECE) <sup>b)</sup>                | K1   |
| Front cone without light                 | 03020171 | Front cone for MANIPOLO ABLATORE LED COMBI TOUCH (COMBI TOUCH LED SCALER HANDPIECE) <sup>b)</sup> | L1   |

| Accessories available for separate order |          |   |      |
|--|----------|---|------|
| Item                                     | Code     | Description   | Ref. |
| Front cone for LED handpiece             | 03020205 | Front cone for MANIPOLO ABLATORE LED COMBI TOUCH (COMBI TOUCH LED SCALER HANDPIECE) <sup>b)</sup> | L2   |
| PROPHY powder container                  | 03020194 | Container only, without cap <sup>b)</sup>   | U1   |
| PERIO powder container                   | 03020200 | Container only, without cap <sup>b)</sup>   | U2   |
| Powder container cap                     | 03020195 | <sup>b)</sup>   | V    |
| Light guide kit                          | 02900146 | Front cone, light guide and decorative ring for SLIM scaler handpiece <sup>b)</sup>               | X    |

**Table 4** – Accessories available for separate order.

- b) Produced by Mectron.
- c) Distributed by Mectron.

- d) To be used only for ME type inserts.
- e) To be used only for the "Crown Preparation" kit.

The packaging of the appliance cannot undergo strong impacts as contains electronic components, therefore the transport and the storage must be carried out with particular care.

Do not stack various boxes in order to avoid pinching the underlying packaging.

All the material sent by MECTRON was controlled at the time of dispatch.

The appliance is shipped appropriately protected and packaged.

Upon receipt of the appliance, check for any possible damage caused during transport and

in case any damage and/or defects is found, complain to the transporter.

Keep the packaging in case any items need to be sent to an Authorised MECTRON Service Centre and to store the appliance during long periods of non-use.

**⚠ WARNING:** Before starting the treatment, always make sure to have a stock of material (scaler handpiece, inserts, wrenches, air-polishing handpieces, SUBGINGIVAL PERIO TIPS) to use in the event of failures or drawbacks.

## 4 INSTALLATION

### 4.1 First Installation

The appliance must be installed in a suitable place that is convenient for its use.

**⚠ WARNING:** The place where the device is installed must meet the requirements set out in *Chapter 4.2 on page 17*.

combi touch can be purchased ready for use

or may need to be enabled by typing in an activation key.

In case your appliance requires an activation key, the procedures to be followed may vary from country to country.

Always refer to your dealer for details.

### 4.2 Safety Requirements During Installation

**⚠ WARNING: Contraindications. Interference with other equipment.** Even if compliant with standard IEC 60601-1-2, combi touch may interfere with other devices in the vicinity. combi touch must not be used near or stacked with other appliances. However, if this should become necessary, it is necessary to check and monitor the correct functioning of the appliance in that configuration.

**⚠ WARNING: Contraindications. Interference from other equipment.** An electrosurgical scalpel or other electrosurgical units positioned near the combi touch appliance may interfere with the correct operation of the appliance itself.

**⚠ CAUTION:** The electrical system of the places where the appliance will be installed and used must comply with the laws in force and the relative electrical safety regulations.

**⚠ CAUTION:** To avoid the risk of electrical shocks, this appliance has to be connected exclusively to a power supply with an earth safety device.

**⚠ CAUTION:** Place the appliance in such a way as to always have the power plug easily accessible as it is considered to be a disconnection device.

**⚠ WARNING: Risk of explosion.** The appliance cannot operate in environments where there are saturated atmospheres of flammable gases (aesthetic mixtures, oxygen, etc.).

**⚠ WARNING:** Install the appliance in a safe place protected from impact or accidental water or liquid spray.

**⚠ WARNING:** Do not install the appliance above or near sources of heat. Arrange during installation for a suitable circulation of air around the appliance.

**⚠ CAUTION:** Do not expose the appliance to direct sunlight or sources of UV light.

**⚠ CAUTION:** The appliance can be transported but must be handled with care when moved. Position the pedal on the floor so that it can only be activated intentionally by the operator.

**⚠ CAUTION:** Before connecting the handpiece to its cord, check that the electrical contacts are perfectly dry, on both parts. If necessary dry them with compressed air.


**⚠ CAUTION:** Each tank can hold liquids up to a maximum of 500 ml.

### 4.3 Connecting the Accessories

Connect the pedal to the back of the appliance in the socket marked with the symbol  by inserting the pedal cord plug until you hear a 'click'.


**⚠ CAUTION:** Pay attention to the positioning of the pedal, which must be such that the pedal is only activated intentionally by the operator.

1


Drain the condensate from the compressed air system. Connect the supply hose to the compressed air circuit in the medical practice using a suitable reduction and shut-off valve (not included in the Mectron supply). Connect the quick coupling to the male attachment on the rear of the appliance;

2


Connect the supply hose to the water circuit in the medical practice using a suitable reduction and shut-off valve (not included in the Mectron supply). Connect the quick coupling to the male attachment on the rear of the appliance;

3

Insert the supply cable in its connection situated on the back of the appliance. Connect it to the socket in the wall;

4


**⚠ CAUTION:** Position the device in a way that the power plug is easily accessible at all times; this is considered a means of isolation.

EN

Unscrew the cap of the tank and fill it with the desired liquid;

**⚠ CAUTION:** Each tank can hold liquids up to a maximum of 500 ml.


5

Check that the small tube inside the cap is properly installed, then screw the cap onto the tank;

**⚠ CAUTION:** Check that the female coupling of the tank cap is clean and free from obstructions.

6



**⚠ CAUTION:** Check that the male coupling on the appliance core unit is clean and that its O-rings are not worn.

**⚠ CAUTION:** Use the support only to install the Mectron 500 ml tank and to stow the handpiece. Do not use the support for other purposes.

7


Holding the tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

**⚠ CAUTION:** Do not overturn the tank, the cap is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

Correctly insert the scaler handpiece onto its cord by matching the alignment notch on the handpiece connector with the alignment key on the cord connector. Check that the electrical contacts of both are perfectly dry and if necessary dry them by blowing with compressed air;


8

Position the handpiece on its dedicated support.

9

**⚠ CAUTION:** Use the support only to install the Mectron 500 ml tank and to stow the handpiece. Do not use the support for other purposes.


**⚠ CAUTION:** Before extracting the powder container or unscrewing the caps, check that the appliance is switched off and disconnected from the power mains, or that the "refill" function has been carried out and the corresponding LED stays lit (see [Chapter 5.2.2 on page 28](#)).

10

Unscrew the cap of the PROPHY white powder container on the left.

**⚠ CAUTION:** The PROPHY powder container is white coloured and is positioned on the left.

Pour the specific supragingival powder distributed by Mectron into the container, making sure the level stays below the diffuser located inside.



**⚠ WARNING:** Insert only specific subgingival powder distributed by Mectron in the PROPHY container.

**⚠ CAUTION:** Correct level of powder in the tank.

**Minimum level:** The level of powder in the container must not be lower than one centimetre to prevent cleaning performance from dropping.

**Maximum level:** The level of the powder in the container must remain below the diffuser (at least 5 mm).


11

EN

Close the tap on the container without tightening it excessively.


12

Unscrew the cap of the PERIO light blue powder container on the right.

13

**⚠ CAUTION:** The PERIO powder container is light blue coloured and is positioned on the right.

Pour the specific subgingival powder distributed by Mectron into the container, making sure the level stays below the diffuser located inside.



**⚠ WARNING:** Insert only specific subgingival powder distributed by Mectron in the PERIO container.

**⚠ CAUTION:** Correct level of powder in the tank.

**Minimum level:** The level of powder in the container must not be lower than one centimetre to prevent cleaning performance from dropping.


**Maximum level:** The level of the powder in the container must remain below the diffuser (at least 5 mm).

14

15

Close the tap on the container without tightening it excessively.

## 5 USE

### 5.1 Power On/Off

#### Switching the appliance on

The switch is positioned at the rear of the device, on the left.

Flick the switch into the "I" position, taking care not to press the pedal.

All the warning lights on the device will turn on and off again. The keyboard stays off for a few moments until an acoustic signal indicates the end of the diagnostic cycle. At the end of the diagnostic cycle, the system loads the preset configuration and is ready for use.

Preset configuration:


- ULTRASOUND part:
  - mode: "endo";
  - power: "1";
  - light: "Off";
  - irrigation: tank - medium flow.
- POLISHING part:
  - "Prophy" function.

**NOTE:** The settings of the scaler part and polisher part can both be modified only with the handpieces stowed in their housings. If either of the two handpieces are extracted, only the settings relative to the active function can be modified.

#### Switching the appliance off

The switch is positioned at the rear of the device, on the left.

Flick the switch into the "O" position, taking care not to press the pedal.


**⚠ CAUTION:** Position the device in a way that the power plug is easily accessible at all times; this is considered a means of isolation.

## 5.2 Description of the Keyboard

### TOUCH KEYBOARD

The user can configure the touch keyboard by simply touching the system by. Depending on the configuration, the electronic feedback system will automatically adjust the correct working frequency.

**NOTE:** A brief acoustic signal is emitted to confirm a key selection. A prolonged acoustic signal indicates that the current configuration does not allow the selection of the pressed key.


### 5.2.1 ULTRASOUND - Scaler Part

#### FUNCTIONS

(Ref. 1 inside cover)

Depending on the type of operation, it is possible to select one of 3 options from the "function" list, as follows:

- **endo:** dedicated to endodontic treatments such as root canal cleaning and the retrograde approach.
- **perio/scaler:** dedicated to all supra- and subgingival dental prophylaxis, root planing, and implant surface cleaning procedures.
- **restorative:** dedicated to preservative techniques, implants and extractions.

**NOTE:** when "restorative" is selected with power level 6, the PULSE function is activated, which is used to optimise the performance of inserts used in extraction and prosthetic techniques.

EN

**SOFT MODE****(Ref. 2 inside cover)**

The treatment can be made more delicate for sensitive patients, by activating the "soft mode" function.

"Soft mode" can only be activated in the "perio/scaler" function with power levels from 1 to 5 and "restorative" function with power levels from 1 to 4.

In "endo" mode the function is not available.

**POWER****(Ref. 3 inside cover)**




For all "mode" functions, the operating power can be adjusted by selecting the numbers on the "power" scale. There are 5 power levels, from 1 to 5. The operating power can be adjusted incrementally (1: minimum power, 5: maximum power).

**LIGHT****(Ref. 4 inside cover)**

Depending on the type of treatment to be performed, the "light" function can be activated:


- When the "on/off" key is activated, the LED light on the scaler handpiece nose cone lights up when the pedal is pressed and automatically switches off 3 seconds after the pedal is released;
- When the "on/off" key is deactivated, the LED light on the scaler handpiece nose cone stays off.

**NOTE:** Adopting the "MANIPOLO ABLATORE LED COMBI TOUCH" ("COMBI TOUCH LED SCALER HANDPIECE" - cod. 03120224) equipped with the "Front cone without light" (cod. 03020171), this function will not be available.

The position of the LED light on the nose cone of the "MANIPOLO ABLATORE LED COMBI TOUCH" ("COMBI TOUCH LED SCALER HANDPIECE" cod. 03120224) can be adjusted as follows:

1. Hold the handpiece body and slightly unscrew the metal ring at the base of the nose cone by turning it counter-clockwise;
2. Rotate the nose cone so that the LED light is in the desired and necessary position;
3. To fix the position, screw the metal ring in a clockwise direction.

## IRRIGATION



### (Ref. 5 inside cover)

The appliance, for the scaler part, allows the use of two types of irrigation:

- the normal water circuit or
- the tank circuit.

The flow rate of the two circuits can be regulated as follows:

- for the normal water circuit, the flow can be continuously regulated using the left-hand knob;
- for the tank circuit, the flow can be continuously regulated on 7 levels via touch screen using the “-” and “+” keys.


Irrigation can be excluded in “endo” and “perio/scaler” modes with power levels from 1 to 5 and in “restorative” mode with power levels from 1 to 4.

**! CAUTION:** The exclusion of irrigation with non “Dry Work” inserts can cause overheating with consequent breakage of the handpiece.

**NOTE:** On regulation of the irrigation flow rate (“-” and “+”) a prolonged acoustic signal is emitted when high full scale is reached.

**NOTE:** the irrigation flow rate can be adjusted using the “-” and “+” keys also during the treatment (holding the pedal pressed).

## “FLUSH” FUNCTION

(Ref. 6 inside cover)


The appliance has a “flush” key, which allows the irrigation circuit to be filled and rinsed.

**Filling the irrigation circuit:** Before starting the treatment, the “flush” function can be used to allow the liquid to reach the insert, so as to start the treatment with the necessary irrigation.

### Cleaning the hydraulic circuit:

The “flush” function allows the irrigation circuit rinse cycle to be performed. This function must be used at the end of the intervention and before cleaning and sterilising all parts.

**NOTE:** If both handpieces are stowed in their housings, the “flush” function cannot be activated.


## 5.2.2 POLISHING - Polisher Part

### FUNCTIONS

#### (Ref. 7 inside cover)

Depending on the type of application, one of the 2 types of powder can be selected as follows:


- “prophy”: dedicated to the clinical indications of the specific supragingival powder distributed by Mectron;
- “perio”: dedicated to the clinical indications of the specific subgingival powder distributed by Mectron.

#### “Refill” FUNCTION

#### (Ref. 8 inside cover)

The "refill" function must be used to depressurise the powder containers in order to then be able to open and remove them from the appliance.


### IRRIGATION

#### (Ref. 9 inside cover)

The appliance uses just normal water for the polisher part.

The water circuit flow rate can be regulated continuously through the right hand knob;


EN

### “Flush” function (Ref. 6 inside cover)

The “flush” function can be used when the type of powder needs to be changed, to make sure the powder circuit has been properly cleaned of the previously used powder.

**NOTE:** If both handpieces are stowed in their housings, the “flush” function cannot be activated.


EN

### 5.2.3 Symbols

#### (Ref. 10 inside cover)

combi touch is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch keyboard by means of a symbol.

To help the user identify the malfunctioning component, there are three symbols described in *Chapter 12.1* on page 112.

## 5.3 Safety Requirements Before and During Use

**⚠ WARNING: Exclusively use original Mectron inserts, accessories and replacement parts.**

**⚠ WARNING:** Before starting work always ensure that you have all replacement parts (handpiece, inserts, keys) to use should there be downtime due to a fault or problem.

**⚠ WARNING: Checking the condition of the device before treatment.** Always check that there is no water under the appliance. Before every treatment always check the perfect operation of the appliance and efficiency of the accessories. If operation anomalies are observed do not carry out the treatment. Contact an Authorised Mectron Service Centre if the faults concern the device.

**⚠ WARNING: Infection control. First use:** All parts and reusable accessories (new or returned from an Authorised MECTRON Service Centre) are delivered in NON-STERILE conditions and must be treated before each use, according to the instructions provided in *Chapter 8* on page 49. **Subsequent uses:** After each treatment, clean and sterilise all parts and reusable accessories according to the instructions in *Chapter 8* on page 49.

**⚠ WARNING: Infection control.** Do not leave liquids in the tank for long periods of time. The tank must be filled after having cleaned and sterilised all parts, and before a treatment. If the tank has been filled without having used the device, it must be emptied at the end of the day, and all parts and accessories must be cleaned and sterilised.

**⚠ WARNING:** The patient must not come into contact with the machine core unit or with the pedal.

**⚠ WARNING:** During the intervention on the patient, do not perform any maintenance tasks on the system.

## 5.3.1 Ultrasonic Scaler

**⚠ WARNING: Use of non-original inserts not produced by MECTRON:** this entails permanent damage to the thread of the handpiece and compromises correct operation with the risk of causing harm to the patient.

**ⓘ CAUTION: Contraindications.** Do not perform treatments on metal or ceramic prosthetics. Ultrasonic vibrations could lead to decementation of the dentures.

**⚠ WARNING: Contraindications.** Do not use combi touch on patients with Pacemakers or other implantable electronic devices. This regulation also stands for the operator.

**⚠ WARNING:** Do not perform scaling treatments without water spray in order to avoid an overheating of the insert that can cause damage to the tooth. Treatments without water spray can only be those performed with the "Dry Work" inserts without water passage.

**ⓘ CAUTION:** In treatments requiring irrigation, use only inserts with a liquid passage.

**⚠ WARNING: Treatments that require irrigation.** Always check the operation of the irrigation before and during use. Make sure that the liquid exits from the insert. Do not use the appliance if the irrigation does not work or if the pump is faulty.

**⚠ WARNING:** To ensure that the handpiece cools down activate it always with irrigation circuit installed and properly filled. To fill the irrigation circuit, always use the flush function.

**ⓘ CAUTION:** For correct use of the appliance, the pedal must be pressed and the appliance started with the insert not in contact with the part to be treated, so that the electronic circuit is able to recognise the best point of resonance of the insert without interference, thus allowing optimum performance.

**⚠ WARNING:** Before each treatment ensure that the insert that is suitable for the treatment is inserted on the handpiece. Use only the Mectron torque wrench to secure the insert to the handpiece.

**⚠ WARNING:** Do not change the insert whilst the handpiece is in operation in order to avoid causing injury to the operator.

**ⓘ CAUTION: Flush function.** After use with aggressive and non-aggressive solutions, the tubes and handpiece must be cleaned using the flush function (See *Chapter 6 on page 40*). If the tubes are not cleaned, the crystallisation of the salts may seriously damage the device.

**⚠ WARNING: Insert breakage and wear.** High frequency oscillations and wear can, in rare cases, lead to the breakage of the insert. Do not bend, change shape or sharpen an insert in no way. Folding an insert or levering on it can cause the insert to break. Deformed or otherwise damaged inserts, are susceptible to breakage during use. These inserts should never be used. Excessive pressure on the inserts during use can lead to breakage. In the event of breakage check that no fragments remain in the treated part and at the same time use suction effectively to remove them. The patient must be instructed to breathe through the nose during the treatment, or use a dental dam, so as to avoid ingesting fragments of broken inserts.

Check the state of wear of the insert and its integrity before and during each use. If a drop in performance occurs, see to its replacement.

The state of wear of the most common inserts (S1, S1-S, S2, S5, P2, P4, P10) can be checked using the supplied INSERT-CARD. To use the INSERT-CARD correctly:

- Position the insert on the INSERT-CARD so that the profile corresponds to the one printed on the card. The printed profile on the card has a red line indicating the limit of wear;
- If the insert is shorter than the limit of wear, its performance will be significantly inferior compared to that of a new insert, and should therefore be replaced.

If the layer of titanium nitride (gold-plated surface), where present, is visibly worn, the insert must be replaced. The use of a worn insert reduces its efficiency.

**Diamond Inserts:** diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.

When the nitriding is consumed, the cutting edge loses effectiveness; a possible regrinding harms the insert, therefore it is prohibited. Verify that the insert is not worn. During the operation, frequently check that the insert is intact, especially in the apical part.

During the operation avoid prolonged contact with the retractor or with metal instruments in use.

⚠ **CAUTION: Contraindications.** After autoclaving the handpiece, the inserts, the torque wrench, or any other sterilisable accessory, wait until it has completely cooled before reuse.

⚠ **CAUTION: The electrical contacts inside the handpiece and cord connectors must be dry.** Before connecting the handpiece to its cord, make sure that the electrical contacts of the connector on both sides are perfectly dry, especially after the sterilisation cycle in an autoclave. If necessary dry the contacts by blowing them with compressed air.

⚠ **CAUTION:** Owing to its conformation, the handpiece can rotate. When not used the handpiece must always be put back on its support stand.

### 5.3.2 Jet Polisher

**⚠ WARNING: Contraindications.** Patients on a sodium restricted diet must not receive treatment with sodium bicarbonate powders.

**⚠ WARNING: Contraindications.** Patients who suffer from serious respiratory problems such as chronic bronchitis, asthma, emphysema, etc., must not receive prophylaxis treatment unless otherwise indicated by their doctor.

**⚠ WARNING: Contraindications.** Patients wearing contact lenses or glasses must remove them prior to receiving treatment with the jet polisher.

**⚠ WARNING: Contraindications - Jet polisher.** Do not aim the air, supragingival powder, or water jet onto the soft tissues or inside the gingival sulcus. Failure to comply with this prescription can cause a gingival tissue emphysema (emphysema of the mucous and/or subcutaneous). For this type of application, use only subgingival powder.

**⚠ WARNING: Contraindications.** Do not use the device near areas subject to recent dental extraction and in traumatised/damaged areas (or areas nearby) due to the risk of emphysema.

**⚠ WARNING: Temperature of the water spray.** The appliance is equipped with a double safety device that controls the temperature of the water spray. Before the treatment, it is however recommended to instruct the patient to inform the operator if he perceives an excessive increase in the temperature of the water.

**⚠ WARNING: Infection control and cleaning the water and air circuits.** To ensure utmost patient and operator safety, after each treatment follow the indications provided in *Chapter 8 on page 49*.

**⚠ WARNING: Do not use the appliance without water.** Make sure that the appliance is connected to the water circuit and the water tap is open.

**ⓘ CAUTION:** Do not attempt to unscrew the tap of the powder container before having performed the "refill" cycle.

EN


## 5.4 Instructions for Use - Scaler Part

After having connected all the accessories as described in *Chapter 4.3 on page 18* proceed as follows:

1

Lift the scaler handpiece, with or without insert, fill the irrigation circuit by selecting "flush" on the keyboard. The "flush" function LED flashes. The appliance allows the use of two types of irrigation: the normal water circuit or the tank circuit.

**ⓘ CAUTION:** Make sure the polisher handpiece is stowed in its housing, otherwise the system will remain inactive.


2

Press the pedal for a moment to start the cycle. The "flush" function and selected irrigation type LEDs will flash. The performance of the cycle is indicated with brief sequence acoustic signals.


The cycle lasts 23 seconds, but can be interrupted by pressing the pedal as soon as the liquid starts coming out of the handpiece.

At the end of the "flush" cycle, the device becomes available again with the most recently used configuration;

3

Screw in the preselected insert on the combi touch handpiece until it comes into contact;


4

Tighten the insert using the Mectron torque wrench.

For correct use of the Mectron torque wrench, proceed as follows:

- Place the insert inside the torque wrench, as shown;


Hold the handpiece body firmly;

**⚠ CAUTION:** Do not grab the handpiece on the end and/or on the cord, but only on the body. Do not turn the handpiece, rather it must be held steady turning only the wrench.

Turn the torque wrench in a clockwise direction until the notch clicks (the external body of the torque wrench turns with respect to the body of the handpiece, emitting mechanical "CLICK" sounds).

The insert is now optimally locked;


5

Select the type of function, power, and irrigation necessary, and the light if desired, on the keyboard;


Lift the handpiece and press the pedal to start treatment;

6

At the end of treatment, stow the scaler handpiece in its housing.

7

## 5.5 Important information on the Inserts

### **WARNING:**

- Before using the sterilised insert, check the integrity of the sterile packaging and inspect the product to ensure there is no damage. The insert will no longer be sterile if the packaging is broken or damaged. If the packaging is damaged, the insert **MUST** be sterilised again before use.
  - When the layer of titanium nitride is visibly consumed the insert must be replaced. The use of an insert too worn decreases its efficiency.
  - **Diamond Inserts:** diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.
  - Do not activate the handpiece while the insert is in contact with the part to be treated so that the electronic circuit can recognise the best point of resonance of the insert allowing optimum performance.
  - Check the state of wear of the insert and its integrity before and during each use. If you experience a loss of performance provide for its replacement.
  - Use only original Mectron inserts. The use of non-original inserts will void the warranty and damage the thread of the combi touch handpiece, with the risk of no longer being able to correctly screw in the original inserts for subsequent use. In addition, the machine settings are tested and guaranteed for correct operation using only original Mectron inserts.
  - Do not change in any way the shape of the insert, by bending it or filing it. This could cause it to break.
  - Do not use an insert that has undergone deformation of any type.
- Do not attempt to sharpen an insert used.
  - Always check that the threaded parts of the insert and the handpiece are perfectly clean - See *Chapter 8 on page 49 - Cleaning and sterilisation.*
  - Excessive pressure applied to the insert can cause it to break and potentially injure the patient.
  - The Mectron inserts vibrate with longitudinal oscillation, with a back and forward movement. During the treatment, always keep the instrument in a tangential direction with respect to the surface of the tooth. Move the handpiece back and forward, applying slight side pressure.
  - Do not point the instrument directly onto the surface of the enamel or implant. Position the tip/operating part only in a tangential direction with respect to the surface of the tooth or implant.
  - The insert must be held in constant movement. If the insert is blocked it may cause overheating of the treated part. It is recommended to use continuous motions to minimise contact between the tip and the part. Do not block them against the tissues so as not to cause their overheating. It is advisable to increase the level of irrigation as the power level increases.
  - Leave the ultrasonic vibrations to work, do not exert excessive pressure on the inserts during use. Apply slight force on the insert for optimum performance.
  - When the insert is used in interproximal spaces, do not block the instrument or leverage on the operating part. The inserts must be left free to vibrate.

- In endodontic root canal therapy, never operate the files outside of the root canal to prevent them from breaking. To prevent breakages, create a smooth path with a manual endo file and plan the most direct possible path of access to minimise any bending of the insert.

Make delicate movements. Examine the file often for signs of wear. If the file breaks inside the canal, avoid contact between the instrument and the broken file to stop it from being pushed further in. Do not apply pressure on the insert in an axial direction.


## 5.6 Instructions for Use - Polisher Part

After having connected all the accessories as described in *Chapter 4.3* on page 18 proceed as follows:

1 Lift the quick coupling on the air-polishing handpiece cord. The POLISHING function is activated.

**⚠ CAUTION:** Make sure the scaler handpiece is stowed in its housing, otherwise the system will remain inactive.


1

2 Choose the type of air-polishing handpiece depending on the operation that needs to be performed;


**⚠ CAUTION:** The SUBGINGIVAL PERIO TIP can only be used with the perio air-polishing handpiece.

2

3 On the quick coupling of the cord, slide the ring and hold it in position, then insert the air-polishing handpiece until it comes into contact, and finally release the ring again;

3

Select the desired “prophy” or “perio” function.

- “Prophy” function: use the specific supragingival powder distributed by Mectron;
- “Perio” function: use the specific subgingival powder distributed by.

**⚠ WARNING:** Before performing any other operations, make sure the air-polishing handpiece is properly and completely inserted in the quick coupling (the quick coupling ring and air-polishing handpiece must be in contact) and that the ring has been returned to the initial position.

**⚠ CAUTION:** If the PERIO air-polishing handpiece has been selected, the SUBGINGIVAL PERIO TIP must be inserted on the front part.

**⚠ CAUTION:** Handle the SUBGINGIVAL PERIO TIP with care.


**⚠ CAUTION:** With the PERIO air-polishing handpiece, only specific subgingival powder distributed by Mectron must be used.

Insert the SUBGINGIVAL PERIO TIP on the PERIO air-polishing handpiece, pushing it in until it comes into contact.


**⚠ WARNING:** Make sure the SUBGINGIVAL PERIO TIP is properly inserted on the air-polishing handpiece; the two pieces must be in contact.

**⚠ CAUTION:** Use only original MECTRON accessories.

4


5

EN


Press the pedal to start treatment. The flow of water can be regulated via the right hand knob until the desired amount is reached.

6

At the end of the treatment, put the polisher handpiece back in its seat.

7

## 5.6.1 "Refill" function

The "refill" function allows the powder containers to be depressurised, enabling their opening or removal and preventing powder from escaping. This function must be used whenever a container is to be reloaded or cleaned.


**NOTE:** The powder containers are only pressurised when, after being selected, the pedal is pressed.

**⚠ CAUTION:** During the "refill" cycle, air and powder escape from the polisher handpiece.

To activate the "refill" function, proceed as follows:


Press the "refill" button.

1

Wait for completion of the cycle, during which the LED of the "refill" key will flash together with that of the pressurised powder container, emitting an acoustic signal.

2


**NOTE:** If both the powder containers are pressurised, the "refill" cycle will depressurise both.

**NOTE:** The duration of the "refill" cycle varies depending on whether one or both pressurised powder containers are present.

3


At the end of the cycle, the "refill" key LED will remain lit and it will therefore be possible to open or remove the containers.

4

At this point simply close or re-insert the containers to continue using the machine.

**NOTE:** If the powder containers are neither opened nor removed, suffice it to simply press the "refill" key to return to the last used program.

## 6 FLUSH FUNCTION

Using the "flush" function it is possible to refill or rinse the irrigation circuit.

**⚠ CAUTION: Flush function.**  
The "flush" function must be used after each treatment, before starting the cleaning and sterilisation procedures.

**⚠ CAUTION:** If the tubes are not rinsed, the crystallisation of the salts may seriously damage the appliance.


EN

### 6.1 ULTRASOUND - Scaler Part

Fill the tank with the appropriate liquid. Depending on the purpose of the FLUSH, different liquids can be used:


- If the tank irrigation circuit needs to be filled, use the solution suited to the treatment;
- If the tank irrigation circuit needs a rinse cycle, use water, preferably demineralised.

1

Connect the tank to the device;


2

Extract the scaler handpiece from its housing;


**⚠ CAUTION:** Make sure the polisher handpiece is stowed in its housing, otherwise the system will remain inactive.

3

4

Place the scaler handpiece, with or without insert, over a container to catch the liquid released during the rinse cycle;


5

To enter "flush" mode, select "flush" on the touch keyboard.

The previously selected irrigation type will remain active and the corresponding LED will flash (faded out).


All the other selectable options on the keyboard will be disabled.

**NOTE:** The "flush" mode can be deactivated at any time by pressing the "flush" key again. The keyboard will be reactivated and configured in the last setting used.

6

The irrigation circuit used to carry out the "flush" cycle can be changed by pressing the key showing a tap and tank;


7

Press the pedal once and release it to start the cycle.


The function and selected irrigation type LEDs will flash. A sequence of short acoustic signals indicates that the cycle is running;

**NOTE:** The function can be interrupted at any time by briefly pressing the pedal. The keyboard will be reactivated and configured in the last setting used.

8

The cycle lasts 23 seconds. Once ended, the keyboard will be reactivated and configured in the last setting used;


9

If the "flush" cycle is to be performed on the water circuit, repeat the steps from point 4 selecting irrigation with water.

At the end of the selected function time or at the end of the "flush" cycle, the system will be reactivated and configured in the last setting used.

**NOTE:** If running a water circuit "flush" cycle, remember to open the tap.


After rinsing the irrigation circuit, disassemble the individual parts (see *Chapter 7 on page 45*) and proceed to clean and sterilise them (see *Chapter 8 on page 49*).

10

**⚠ CAUTION:** If the handpiece is stowed back in its support, with the "flush" function selected (flashing LED) but not running, the function is deactivated;

**⚠ CAUTION:** During the "flush" cycle, i.e. when the function has been selected and activated by pressing the pedal, if the handpiece is stowed back in the support, the cycle is not interrupted.

EN


## 6.2 POLISHING - Polisher Part

Lift the air-polishing handpiece or quick coupling on the cord of the air-polishing handpiece;

1


**⚠ CAUTION:** Make sure the scaler handpiece is stowed in its housing, otherwise the system will remain inactive.

**⚠ CAUTION:** Make sure that both powder containers are inserted in the device before activating the "flush" function.


2

Position the quick coupling on the cord of the air-polishing handpiece, with or without the air-polishing handpiece, over a container or sink to catch the liquid and powder that will be released during the "flush" cycle;


To enter “flush” mode, select “flush” on the touch keyboard: both the powder containers and corresponding LEDs will flash (faded out). All the other selectable options on the keyboard will be disabled.

3

Press the pedal once and release it to start the cycle. The LEDs of the function and powder containers will flash. A sequence of short acoustic signals indicates that the cycle is running;


4

**NOTE:** The execution of the “flush” cycle cannot be interrupted. Before running the “flush” cycle, check that the right-hand tap is open.

5

The “flush” cycle is run on both powder containers simultaneously. The cycle lasts about 40 seconds. Once ended, the keyboard will be reactivated and configured in the last setting used.

EN


At the end of the irrigation circuit "flush" cycle, disassemble the individual parts (see *Chapter 7 on page 45*) and proceed to clean and sterilise them (see *Chapter 8 on page 49*).

**⚠ CAUTION:** If the handpiece is stowed back in its support, with the "flush" function selected (flashing LED) but not running, the function is deactivated;

**⚠ CAUTION:** During the "flush" cycle, i.e. when the function has been selected and activated by pressing the pedal, if the handpiece is stowed back in the support, the cycle is not interrupted.

6


EN

## 7 DISASSEMBLING PARTS FOR CLEANING AND STERILISATION

Before carrying out the cleaning procedures described in *Chapter 8 on page 49*, disconnect all accessories and components of the combi touch .

**⚠ WARNING:** Switch the appliance off. Always switch the appliance off using the switch and disconnect the power supply cable from the wall socket and machine core unit before carrying out cleaning and sterilisation tasks.

1


Disconnect the pedal from the appliance: grasp the pedal connector, press the release tab and pull the connector back;

**⚠ CAUTION:** Do not try to unscrew or turn the connector during disconnection: this might damage the connector.


**⚠ CAUTION:** While disconnecting the cable from the pedal, keep always and only the connector of the cord. Never pull the cord itself.

2


If present, unscrew the insert from the handpiece using the torque wrench;

3


Disconnect the scaler handpiece from the cord;

4

**⚠ CAUTION:** Do not try to unscrew or turn the connector when disconnecting the handpiece. The connector can be damaged.


Unscrew the nose cone from the scaler handpiece;

**⚠ CAUTION:** The nose cone has a light pipe. If the nose cone is unscrewed, the light pipe will no longer be held in place and may slide and be disconnected. Be careful not to lose the light pipe.

**NOTE:** On the "Front cone for LED handpiece" (cod. 03020205) of the MANIPOLO ABLATORE LED COMBI TOUCH (COMBI TOUCH LED SCALER HANDPIECE - cod. 03120224), the metal ring cannot be separated from the plastic cone.


5

EN


Remove the light pipe from the nose cone;

6


Disconnect the tank from the machine core unit, pulling it outwards;

7


On the quick coupling of the cord, slide the ring and hold it in position, then extract the air-polishing handpiece, and finally release the ring of the handpiece;

8

If the PERIO air-polishing handpiece was used, remove the SUBGINGIVAL PERIO TIP using the supplied K9 wrench and proceed with disposal (see Chapter 10 on page 102);


9

**NOTE:** Insert the K9 wrench exactly in the position shown in the figure.

**! CAUTION:** Before extracting the powder container or unscrewing the caps, check that the appliance is switched off and disconnected from the power mains, or run the "refill" function (see Chapter 5.6.1 on page 38).

10


Extract the powder container from the appliance, remove the cap and empty it.

Proceed with polishing. (see Chapter 9.4 on page 101)

Repeat the operation on the second powder container of the appliance if both have been used.

When the cleaning operation has been completed, re-position the containers in the appliance.

11

**! CAUTION:** Insert the powder containers in such a way that the groove of the bottom of the appliance is aligned with the convexity of the powder container.

## 8 CLEANING AND STERILISATION

This table is purely indicative. For the complete cleaning and sterilisation procedures of the individual parts, refer to the paragraphs indicated in the table.

**⚠ CAUTION:** Methods not considered in the table below must not be used.

**⚠ WARNING:** The sterile single-use SUBGINGIVAL PERIO TIPs must be used for one treatment and one patient only. The single-use SUBGINGIVAL PERIO TIPs must not be reused. Sort and dispose of each single-use SUBGINGIVAL PERIO TIP in accordance with current hospital waste disposal laws.

### Machine Core Unit, Pedal, Quick Coupling and Relative Cord of Polisher Handpiece and Scaler Handpiece

| Phase | Chapter | Procedure                                |
|-------|---------|--|
| I     | 8.1     | Preparation                              |
| II    | 8.2     | Manual cleaning with enzymatic detergent |

### Tank and Cap

| Phase | Chapter | Procedure  |
|-------|---------|--|
| III   | 8.3     | Manual cleaning with enzymatic detergent and running water |

### Accessories (Scaler Handpiece, Scaler Nose Cone, Inserts, Air-polishing Handpieces, K9 Torque, Torque Wrench)

| Phase | Chapter | Procedure              | Scaler Hand-piece                        | Scaler Nose Cone | In-serts | Air-pol-ishing hand-pieces                                | K9 Wrench | Torque wrench |                                    |
|-------|---------|------------------------|--|------------------|----------|---|-----------|---------------|------------------------------------|
| IV    | 8.4.1   | Manual cleaning        | Manual cleaning with Enzymatic Detergent |                  |          | Immersion into enzymatic detergent<br>Ultrasonic Cleaning |           |               |                                    |
|       | 8.4.2   | Automatic cleaning     |  |                  |          |   |           |               | Thermal-disinfector with detergent |
| V     | 8.5     | Check cleaning         | Visual inspection with Lens 2.5x         |                  |          |   |           |               |                                    |
| VI    | 8.6     | Drying and lubrication | Drying                                   |                  |          |   |           |               | Drying and Lu-<br>brica-<br>tion   |
| VII   | 8.7     | Sterilisa-<br>tion     | Packaging, Sterilisation and Storage     |                  |          |   |           |               |                                    |

**NOTE:** Repeated reconditioning has a minimal effect on the devices and their accessories. The end of the service life of the accessories is generally determined by wear or damage resulting from use. Mectron guarantees the integrity of its sterilisable air-polishing nozzles and scaler handpieces for up to 250 reconditioning cycles.

## 8.1 Preparation

1. Run the flush function (see *Chapter 6 on page 40*);
2. Disinfect the external irrigation circuit as described in *Chapter 8.1.1 on page 51*.
3. Disinfect the tank irrigation circuit as described in *Chapter 8.1.2 on page 68*.
4. Check that all of the following accessories have been removed/disconnected from the machine core unit (see *Chapter 7 on page 45*):
  - Power supply cable;
  - Pedal;
  - Scaler handpiece;
  - Inserts;
  - Air-polishing handpiece;
  - Water and air hose.

**ⓘ CAUTION:** The cleaning and sterilising operations described in the following sections are to be performed at first use and at all subsequent uses.

**⚠ WARNING:** Always switch off the appliance using the O/I switch and disconnect it from the power mains before carrying out the post-preparation cleaning tasks.

**ⓘ CAUTION:** Always disconnect the insert from the handpiece before cleaning and sterilising it.

**ⓘ CAUTION:** Do not immerse the handpiece in disinfecting solutions or other liquids as these could damage it.

**ⓘ CAUTION:** In case of excessive powder humidity, remove/disconnect the air-polishing handpiece from the quick coupling on the cord, free the air-polishing handpiece channel of any powder residue using the supplied Ø 0.8 mm cleaning needle, then carry out the "flush" function without the air-polishing handpiece.

## 8.1.1 Disinfecting the External Irrigation Circuit

After performing the flush function (see *Chapter 6 on page 40*) and before carrying out the subsequent cleaning procedures, proceed to disinfect the external irrigation circuit.

The procedure must be carried out on both the scaler handpiece and polisher handpiece.

To disinfect the irrigation circuit, proceed as follows:


### 8.1.1.1 Required Materials

- CIDEX® OPA disinfectant solution;
- Distilled water;
- "Mectron Water Circuit Cleaning Kit" (cod. 02900171-001) to disinfect the external irrigation circuit;
- 60ml syringe with eccentric Luer Slip cone, without needle.

### 8.1.1.2 Procedure - Polisher Handpiece


**⚠ WARNING: Switch the appliance off.** Switch the appliance off using the switch (Ref. 1), disconnect the power supply cable from the wall socket (Ref. 2) and machine core unit (Ref. 3) before carrying out the cleaning and sterilisation tasks.

1


2

Disconnect the quick connector of the external water circuit from the device (see *Chapter 4.3 on page 18*);


3

Unscrew the male coupling for the external water circuit connection (Ref. A - see *Chapter 9.2 on page 100*); Remove the water filter (Ref. B - see *Chapter 9.2 on page 100*);


Check the integrity of the filter, rinse it, dry it with compressed air, making sure to remove any residual impurities;

4


Completely immerse the male coupling for the external water circuit connection in the disinfectant solution (CIDEX® OPA). Leave to soak for 12 minutes at 20 °C ± 2 °C.

5


Remove the male coupling for the external water circuit connection from the disinfectant solution (CIDEX® OPA) and rinse it under running water; Place the male coupling for the external water circuit connection in a clean and safe place for later use.

6

Disinfect the internal surfaces of the seat where the male coupling for the external water circuit connection will be inserted using a clean, soft, low-lint cloth dampened with the disinfectant solution (CIDEX® OPA).

7


8 Screw in the male coupling supplied with the "Mectron Water Circuit Cleaning Kit" in its seat until it comes into contact;


⚠ **CAUTION:** Check the integrity of the "Mectron Water Circuit Cleaning Kit". Replace the "Mectron Water Circuit Cleaning Kit" in case of wear or damage.

8


9 If previously disconnected (for example, for cleaning), connect the hose supplied with the kit to the male water coupling;

9


10 Connect the power supply cable to the device (Ref. 1) and to a wall socket (Ref. 2). Switch the device on using the switch located at the back (Ref. 3);

10


11 Extract the quick coupling on the cord of the air-polishing handpiece and place it above a container or sink to catch the liquid released during the procedure;

11


Completely open the right-hand knob in the front part of the device (Ref. 9 inside cover);

12


Using the syringe, aspirate 60ml of disinfectant (CIDEX® OPA), making sure there is no air inside;

13

Connect the syringe to the end of the hose previously connected to the male water coupling (see point 9 of the procedure);


14

Press the "flush" button on the touch keyboard (see Chapter 6.2 on page 43);

15

**NOTE:** Make sure that both powder containers are present and correctly stowed in their seats, otherwise it will not be possible to activate the "flush" function.

Press and release the pedal. When the "flush" cycle begins, start injecting the disinfectant with the syringe.

**During the "flush" cycle, inject 50ml of disinfectant, leaving the remaining 10ml in the syringe.**

The "flush" cycle on the "Polishing" part lasts about 40 seconds;


**NOTE:** The "flush" cycle on the "Polishing" part cannot be interrupted. Before running the "flush" cycle on the "Polishing" part check that the right-hand knob in the front part of the device (Ref. 9 inside cover) is open. At the end of the "flush" cycle, stop injecting the liquid with the syringe.

**! CAUTION:** Do not press too hard on the piston of the syringe. The irrigation circuit hoses are quite small and it is normal for injection to occur slowly.

**NOTE:** If the 50ml of disinfectant are injected before the flush cycle is completed, wait for it to finish. If the cycle ends before the full 50ml have been injected, repeat the flush cycle until the 50ml have been completely injected.


Wait 12 minutes;

16


EN

17

Press the “flush” button on the touch keyboard (see Chapter 6.2 on page 43);

18

Press and release the pedal. When the “flush” cycle begins, start injecting the remaining disinfectant with the syringe.

**During the “flush” cycle, inject the 10ml of remaining disinfectant.**


The “flush” cycle on the “Polishing” part lasts about 40 seconds;

19

**NOTE:** The “flush” cycle on the “Polishing” part cannot be interrupted. Before running the “flush” cycle on the “Polishing” part check that the right-hand knob in the front part of the device (Ref. 9 inside cover) is open. At the end of the “flush” cycle, stop injecting the liquid with the syringe.


**! CAUTION:** Do not press too hard on the piston of the syringe. The irrigation circuit hoses are quite small and it is normal for injection to occur slowly.

**NOTE:** If the 10ml of disinfectant are injected before the flush cycle is completed, wait for it to finish. If the cycle ends before the full 10ml have been injected, repeat the flush cycle until the 10ml have been completely injected.


Disconnect the syringe and aspirate air;

20


Reconnect the syringe to the end of the hose from which it had been disconnected.;

21


Press the “flush” button on the touch keyboard (see Chapter 6.2 on page 43);

22


Press and release the pedal. When the “flush” cycle begins, start injecting air with the syringe. The “flush” cycle on the “Polishing” part lasts about 40 seconds;

23

Repeat points 20 to 23 another 4 times. The number of cycles run must guarantee the complete drainage of the external irrigation circuit. In any case, check that liquids are no longer coming out of the quick coupling on the cord of the air-polishing handpiece;

24


Disconnect the syringe and repeat points 13 to 24 another 2 times;

25

Disconnect the syringe and fill it with distilled water;


26

**NOTE:** Use of distilled water is recommended for long periods of non-use. If the device is used after the disinfection procedure, use common running water.


Reconnect the syringe to the end of the hose from which it had been disconnected.

27


Press the "flush" button on the touch keyboard (see Chapter 6.2 on page 43).

28

Press and release the pedal and inject the distilled water with the syringe. The "flush" cycle on the "Polishing" part lasts about 40 seconds.

29

Repeat points 26 to 29 another 4 times;

30


Repeat points 20 to 23 twice until the irrigation circuit is completely drained of the water used in the previous steps.

31

**NOTE:** It is advisable to drain the circuit during long periods of non-use. If the device is used after the disinfection procedure, and provided different liquids don't need to be used for irrigation, the water circuit does not need to be drained.


At the end of the disinfection procedure, remove the syringe (Ref. 1), the hose (Ref. 2) and unscrew the male coupling used (Ref. 3).

32

Re-insert the previously cleaned water filter (see point 4 of the procedure) in its seat - Ref. 1 -, screw in the original and previously disinfected male coupling (see points 5 and 6 of the procedure) - Ref. 2 - and reconnect the quick coupling of the external water circuit - Ref. 3.

33


EN

### 8.1.1.3 Procedure - Scaler Handpiece

**⚠ WARNING: Switch the appliance off.** Switch the appliance off using the switch (Ref. 1), disconnect the power supply cable from the wall socket (Ref. 2) and machine core unit (Ref. 3) before carrying out the cleaning and sterilisation tasks.

1


Carry out steps 2 to 10 described in Chapter 8.1.1.2 on page 51;

2

Extract the quick coupling on the cord of the scaler handpiece and place it above a container or sink to catch the liquid released during the procedure;

**⚠ CAUTION:** Before connecting and using the scaler handpiece, make sure that the electrical contacts on both the scaler handpiece and its cord are perfectly dry.


3


4 Completely open the left-hand knob in the front part of the device (Ref. 5 inside cover);

4


5 Using the syringe, aspirate 60ml of disinfectant (CIDEX® OPA), making sure there is no air inside;

5

6 Connect the syringe to the end of the hose previously connected to the male water coupling;


6

7 Press the “flush” button on the touch keyboard (see Chapter 6.2 on page 43);


7

**NOTE:** Make sure that both powder containers are present and correctly stowed in their seats, otherwise it will not be possible to activate the “flush” function.

Select the external irrigation circuit (“tap” symbol - see Chapter 6.1 on page 40);

8


Press and release the pedal. When the “flush” cycle begins, start injecting the disinfectant with the syringe.

9

**During the “flush” cycle, inject 50ml of disinfectant, leaving the remaining 10ml in the syringe.**

The “flush” cycle on the “Ultrasound” part lasts about 23 seconds;

**NOTE:** The “flush” cycle on the “Ultrasound” part can be interrupted at any time by pressing and releasing the pedal. Before running the “flush” cycle on the “Ultrasound” part check that the left-hand knob in the front part of the device (Ref. 5 inside cover) is open. At the end of the “flush” cycle, stop injecting the liquid with the syringe.


**! CAUTION:** Do not press too hard on the syringe. The irrigation circuit hoses are quite small and it is normal for injection to occur slowly.

**NOTE:** If the 50ml of disinfectant are injected before the flush cycle is completed, wait for it to finish. If the cycle ends before the full 50ml have been injected, repeat the flush cycle until the 50ml have been completely injected.

Wait 12 minutes;


10

EN


Press the “flush” button on the touch keyboard (see Chapter 6.2 on page 43);

11

Select the external irrigation circuit (“tap” symbol - see Chapter 6.1 on page 40);

12

13

Press and release the pedal. When the "flush" cycle begins, start injecting the remaining disinfectant with the syringe.


**During the "flush" cycle, inject the 10ml of remaining disinfectant.**

The "flush" cycle on the "Ultrasound" part lasts about 23 seconds;

**NOTE:** The "flush" cycle on the "Ultrasound" part can be interrupted at any time by pressing and releasing the pedal. Before running the "flush" cycle on the "Ultrasound" part check that the left-hand knob in the front part of the device (Ref. 5 inside cover) is open. At the end of the "flush" cycle, stop injecting the liquid with the syringe.


**! CAUTION:** Do not press too hard on the syringe. The irrigation circuit hoses are quite small and it is normal for injection to occur slowly.

**NOTE:** If the 10ml of disinfectant are injected before the flush cycle is completed, wait for it to finish. If the cycle ends before the full 10ml have been injected, repeat the flush cycle until the 50ml have been completely injected.


14

Disconnect the syringe and aspirate air;


Reconnect the syringe to the end of the hose from which it had been disconnected;

15


Press the “flush” button on the touch keyboard (see Chapter 6.2 on page 43);

16


Select the external irrigation circuit (“tap” symbol - see Chapter 6.1 on page 40);

17

Press and release the pedal. When the “flush” cycle begins, start injecting air with the syringe. The “flush” cycle on the “Ultrasound” part lasts about 23 seconds;

18

Repeat points 14 to 18 another 4 times. The number of cycles run must guarantee the complete drainage of the external irrigation circuit. In any case, check that liquids are no longer coming out of the quick coupling on the cord of the scaler handpiece;

**⚠ CAUTION:** Before connecting and using the scaler handpiece, make sure that the electrical contacts on both the scaler handpiece and its cord are perfectly dry.


Disconnect the syringe and repeat points 5 to 19 another 2 times;

Disconnect the syringe and fill it with distilled water;


**NOTE:** Use of distilled water is recommended for long periods of non-use. If the device is used after the disinfection procedure, use common running water.

Reconnect the syringe to the end of the hose from which it had been disconnected.

19


20


21

22


Press the “flush” button on the touch keyboard (see Chapter 6.2 on page 43);

23


Select the external irrigation circuit (“tap” symbol - see Chapter 6.1 on page 40);

24

Press and release the pedal. When the “flush” cycle begins, start injecting the distilled water with the syringe. The “flush” cycle on the “Ultrasound” part lasts about 23 seconds;

25

Repeat points 21 to 25 another 4 times;

26


Repeat points 14 to 18 twice until the irrigation circuit is completely drained of the water used in the previous steps.

27

**NOTE:** It is advisable to drain the circuit during long periods of non-use. If the device is used after the disinfection procedure, and provided different liquids don't need to be used for irrigation, the water circuit does not need to be drained.


At the end of the disinfection procedure, remove the syringe (Ref. 1), the hose (Ref. 2) and unscrew the male coupling used (Ref. 3).

28

Re-insert the previously cleaned water filter (see point 4 of the procedure) in its seat - Ref. 1 -, screw in the original and previously disinfected male coupling - Ref. 2 - and reconnect the quick coupling of the external water circuit - Ref. 3.

29

### 8.1.2 Disinfecting the Tank Irrigation Circuit

After having run the flush function (see *Chapter 6 on page 40*), disinfected the external irrigation circuit (see *Chapter 8.1.1 on page 51*), and before carrying out the subsequent cleaning procedures, proceed to disinfect the tank irrigation circuit. To disinfect the irrigation circuit, proceed as follows:

#### 8.1.2.1 Required Materials


- CIDEX® OPA disinfectant solution;
- Distilled water.

## 8.1.2.2 Procedure


Disconnect the irrigation solution tank from the machine core unit, pulling it outwards;

1


Unscrew the cap from the irrigation solution tank and empty it;  
Rinse the irrigation solution tank with running water;

2


Fill the irrigation solution tank with 250ml of disinfectant (CIDEX® OPA);

3

Holding the irrigation solution tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

4


**⚠ CAUTION:** Do not overturn the tank, the cap is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

Lift the scaler handpiece connector;

**⚠ CAUTION:** Make sure the polisher handpiece is stowed in its housing, otherwise the system will remain inactive.


5

Position the scaler handpiece connector above a container or sink to catch the liquid released during the disinfection cycle;

**⚠ CAUTION:** Before connecting and using the scaler handpiece, make sure that the electrical contacts on both the scaler handpiece and its cord are perfectly dry.

6


Select "flush" on the touch keyboard (see Chapter 6 on page 40).

The previously selected irrigation type will remain active and the corresponding LED will flash (faded out).

All the other selectable options on the keyboard will be disabled.


**NOTE:** The "flush" mode can be deactivated at any time by pressing the "flush" key again. The keyboard will be reactivated and configured in the last setting used.

7

Select the tank irrigation circuit by pressing the "irrigation" key and making sure that the corresponding tank LED lights up;

8


Press the pedal once and release it to start the cycle.

The function and selected irrigation type LEDs will flash. A sequence of short acoustic signals indicates that the cycle is running;

**NOTA:** the function can be interrupted at any time by briefly pressing the pedal. The keyboard will be reactivated and configured in the last setting used.


9

EN

The cycle lasts 23 seconds. Once ended, the keyboard will be reactivated and configured in the last setting used;

10


Repeat the "flush" cycle (see points 7 to 10 of the procedure) a second time;

11

At the end of the "flush" cycles, wait 12 minutes;

12


Repeat the "flush" cycle (see points 7 to 10 of the procedure) 3 times;

13

At the end of the "flush" cycles, wait 12 minutes;

14


Repeat the "flush" cycle (see points 7 to 10 of the procedure) 3 times;

15

At the end of the "flush" cycles, wait 12 minutes;

16


Disconnect the irrigation solution tank from the device;

Unscrew the cap from the irrigation solution tank and drain the remaining liquid;


Rinse the irrigation solution tank with running water;

17

**With the tank disconnected**, run a flush cycle (see points 7 to 10) to empty the irrigation circuit of the liquid inside;


18

Fill the irrigation solution tank with 300ml of distilled water;

**NOTE:** Use of distilled water is recommended for long periods of non-use. If the device is used after the disinfection procedure, use common running water.


19

Holding the irrigation solution tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

**! CAUTION:** Do not overturn the tank, the cap is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

20

Run a number of "flush" cycles (see points 7 to 10) sufficient to completely drain the water from the tank;


21

After cleaning the irrigation circuit, disassemble the individual parts (see Chapter 7 on page 45) and proceed to clean and sterilise them (see Chapter 8 on page 49).

**! CAUTION:** If the handpiece is stowed back in its support, with the "flush" function selected (flashing LED) but not running, the function is deactivated;

**! CAUTION:** During the "flush" cycle, i.e. when the function has been selected and activated by pressing the pedal, if the handpiece is stowed back in the support, the cycle is not interrupted.

22

## 8.2 Cleaning the Non-Sterilisable Parts

The following procedure must be performed on non-sterilisable parts of the appliance, except for the tank and its cap. The parts in question are:

- Machine core unit;
- Pedal and relative connection cable to the machine;

- Quick coupling of polisher handpiece and relative cord;
- Quick coupling of scaler handpiece and relative cord.

**⚠ WARNING:** Always switch off the appliance using the O/I switch and disconnect it from the power mains before carrying out cleaning tasks.

### 8.2.1 Required Materials


- Clean, soft, low-lint cloths;

- Cleaning solution (pH 6-9).

### 8.2.2 Cleaning Method


Clean the surface of the parts in question using a clean, dry, low-lint cloth, dampened with a detergent solution (pH 6-9);

1

Dry the parts using a dry, non-abrasive, low-lint cloth.

2

**⚠ CAUTION:** Do not sterilise the non-sterilisable parts. They may stop working and cause damage to people and/or property.

**⚠ WARNING:** The device and its parts that cannot be sterilised are not protected against the penetration of liquids. Do not spray liquids directly on the surface of the device or its parts that cannot be sterilised.

**⚠ CAUTION:** Do not use running water to clean the non-sterilisable parts.

**⚠ CAUTION:** Do not soak the non-sterilisable parts in liquids and/or various kinds of solutions.

**⚠ CAUTION: DO NOT USE** as disinfectants::

- Very alkaline products (pH > 9);
- Products containing sodium

hypochlorite;

- Products containing hydrogen peroxide;
- Products containing abrasive substances;
- Very acid products (PH < 4);
- Products containing aldehyde, amine and/or phenols
- Acetone;

- Methylethylketone;

since they may discolour and/or damage the plastic materials.

The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.

## 8.3 Cleaning the Tank and Cap

The following procedure must be carried out on the tank and cap of the device.

### 8.3.1 Preparation

1. Disconnect the tank from the machine core unit (see *Chapter 7 on page 45*);
2. Unscrew the cap from the tank.

**ⓘ CAUTION:** Do not sterilise the tank and cap in an autoclave. They could be damaged.

### 8.3.2 Required Materials

- Water;
- Cleaning solution (pH 6-9);
- Clean, soft, low-lint cloth;
- Demineralised water.

- Very acid products (PH < 4);
- Products containing aldehyde, amine and/or phenols
- Acetone;
- Methylethylketone;

**ⓘ CAUTION: DO NOT USE** as disinfectants::


- Very alkaline products (pH > 9);
- Products containing sodium hypochlorite;
- Products containing hydrogen peroxide;
- Products containing abrasive substances;

since they may discolour and/or damage the plastic materials.

The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.


### 8.3.3 Cleaning Method

Thoroughly rinse the inside and outside of both the tank and cap under running water;


Clean the external and internal surfaces of the tank and cap with a clean, soft, low-lint cloth, dampened with a mild detergent solution (pH 6-9);

2


Thoroughly rinse the inside and outside of the tank and cap under running water to eliminate all residue of the detergent solution;

3


Remove any residue on the external or internal surfaces of the tank and cap using a soft, low-lint cloth, dampened with demineralised water;

4

When cleaning operations are over, check the pieces under an appropriate light source, paying attention to any soil residue and, if necessary, repeat the cleaning cycle.

5

## 8.4 Cleaning the Sterilisable Accessories

**NOTE:** The cleaning procedures must be performed immediately after each use. Immerse the insert and/or instrument in demineralised water or in an enzymatic detergent solution immediately after use. Do not leave residue or blood deposits on the inserts and instruments, eliminate larger impurities with a disposable cloth or paper towel.

The sterilisable parts of the appliance are:

- Scaler handpiece;
- Scaler nose cone;
- Inserts;
- Inserts tightening wrench;
- Air-polishing handpiece;
- K9 Wrench.

Before proceeding with operations to check cleanliness (*Chapter 8.5 on page 92*), drying

and lubrication (*Chapter 8.6 on page 94*) and then sterilisation (*Chapter 8.7 on page 96*), one of the two possible cleaning methods explained in detail in the following chapters must be selected.

**⚠ CAUTION:** The instructions supplied below have been validated by the manufacturer of the medical device as ABLE to prepare a medical device for re-use. The process manager is responsible for ensuring that the processes repeated are effectively performed using the equipment, materials and staff in the reprocessing structure in order to obtain the desired result. This generally requires the validation and systematic monitoring of the process. Similarly, all deviations from the instructions provided by the processes manager must be adequately assessed to judge their efficiency and potential undesired consequences.

### 8.4.1 Manual Cleaning

Manual cleaning can be carried out as an alternative to the automatic cleaning described in *Chapter 8.4.2 on page 91*.

#### 8.4.1.1 Required Materials

- Enzymatic detergent at pH 6-9;
- Water;
- Container for immersion in the enzymatic liquid;
- Ultrasonic tank;
- Clean, soft, low-lint cloths;
- Brush with soft nylon bristles;
- Syringe;
- Demineralised water

**⚠ CAUTION: DO NOT USE** as disinfectants::

- Very alkaline products (pH > 9);
- Products containing sodium hypochlorite;
- Products containing hydrogen peroxide;
- Products containing abrasive substances;
- Very acid products (PH < 4);
- Products containing aldehyde, amine and/or phenols
- Acetone;
- Methyleneethylketone;

since they may discolour and/or damage the plastic materials.


The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.

Prepare a solution of enzymatic detergent<sup>a)</sup> with pH 6-9, according to the instructions of the manufacturer;

**⚠ CAUTION:** Once used, dispose of the enzymatic detergent correctly, do not recycle it.

- a) Process validated by independent bodies with Mectron ENZYMEC enzymatic detergent, 0.8% v/v.


1

## 8.4.1.2 Scaler Handpiece

Completely immerse the scaler handpiece, nose cone, and light pipe in the enzymatic solution. Leave to soak for 10 minutes at 22 °C ±2 °C;




2

3


Gently brush the entire surface of the scaler handpiece, nose cone, and light pipe for at least 20 seconds using a brush with soft nylon bristles, paying particular attention to the following areas:

- thread of the scaler handpiece;
- titanium shaft;
  
- external and internal parts of nose cone;
  
- external and internal parts of light pipe.


4

Rinse the inner channel of the scaler handpiece using a 20 ml syringe previously filled with a new enzymatic solution. Repeat twice;


5

Remove the handpiece, nose cone, and light pipe from the enzymatic solution and gently brush their surfaces using a brush with soft nylon bristles under hot ( $40\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ ) running water, for at least 10 minutes;


6

Fill the ultrasonic tank with the enzymatic detergent solution prepared according to the manufacturer's instructions. Use water at ambient temperature ( $22\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ );


7

Place the scaler handpiece, nose cone, and light pipe in the ultrasonic tank immersed in the enzymatic detergent solution at  $22\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$  and run a cycle for at least 10 minutes;


8

Remove the handpiece, nose cone, and light pipe from the enzymatic solution and gently brush their surfaces using a brush with soft nylon bristles under hot ( $40\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ ) running water, for at least 10 minutes;


Place the scaler handpiece, nose cone, and light pipe in the ultrasonic tank immersed in demineralised water at  $40\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$  and run a cycle for at least 10 minutes;

9


Rinse the inner channel of the scaler handpiece using a 20 ml syringe previously filled with demineralised water. Repeat five times;

10

Dry the surface of the scaler handpiece, nose cone, and light pipe with a clean, non-abrasive, low-lint cloth.

11

EN

Process validated by independent bodies with enzymatic detergent.

## 8.4.1.3 Inserts

Place the insert in a clean container, in a horizontal position. Add enough enzymatic detergent solution to completely cover it.

Leave the insert to soak in the enzymatic detergent solution for 10 minutes at ambient temperature ( $22\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ );

While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed.

Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.


**NOTE:** Thoroughly brush all the following parts of the insert for about 20 seconds:

- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- Internal cavities, grooves and fissures.


Remove the insert from the enzymatic detergent solution.

Thoroughly rinse and brush all the surfaces of the insert (see previous point) under running water for at least 10 minutes;


12

13


14


Use a disposable syringe to withdraw and inject the enzymatic cleaning solution into the hard-to-reach areas (through-holes/cannulae). Repeat this operation three times to ensure the effective removal of the soil from the internal surfaces of the through-hole.;

15


Rinse the internal channel of the insert with demineralised water at a pressure of 3.8 bar for at least 10 seconds so as to eliminate any soil residue.;

16


Place the insert in the ultrasonic tank immersed in the enzymatic detergent solution at 24 °C ± 2 °C and run a cycle for at least 20 minutes.;

17


Remove the insert(s) from the ultrasonic tank and rinse them in demineralised water.;

18


Brush the internal and external surfaces of the insert with a clean brush with soft nylon bristles, under running water;

19


Rinse the insert(s) in demineralised water;

20

Rinse the internal channel of the insert with demineralised water at a pressure of 3.8 bar for at least 10 seconds so as to eliminate any soil residue.


21

## 8.4.1.4 Inserts Tightening Wrench

Place the torque wrench in a clean container, in a horizontal position. Add enough enzymatic detergent solution to completely cover the torque wrench. Leave the torque wrench to soak in the enzymatic detergent solution for 10 minutes at ambient temperature (22 °C ±2 °C);

22

While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed. Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.


**NOTE:** Thoroughly brush all the following parts of the torque wrench for about 20 seconds:

- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- External metal ring;
- Internal cavities, grooves and fissures.

Remove the torque wrench from the enzymatic detergent solution. Thoroughly rinse and brush all the surfaces of the torque wrench (see previous point) under running water for at least 10 minutes;


Place the torque wrench in the ultrasonic tank immersed in the enzymatic detergent solution at  $24\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$  and run a cycle for at least 20 minutes;

23


EN

24


25

Remove the torque wrench from the ultrasonic tank and rinse it under running water so as to eliminate all detergent residue.

Brush the internal and external surfaces of the wrench with a clean brush with soft nylon bristles, under running water.

26

Process validated by independent bodies with enzymatic detergent.

### 8.4.1.5 Air-polishing handpieces


Free the channel of the air-polishing handpiece of any powder residue using the supplied  $\varnothing$  0.4 mm cleaning needle, from both sides.

Blow compressed air into the central hole of the handpiece, from both ends;

**⚠ CAUTION:** Clean the channel of the nozzle using only the  $\varnothing$  0.4 mm cleaning needle supplied with the appliance.


**⚠ CAUTION:** If the single-use SUBGINGIVAL PERIO TIP is on the air-polishing handpiece, remove and dispose of it.

27


28

Place the air-polishing handpiece in a clean container, in a horizontal position, and add enough enzymatic detergent solution to completely cover the handpiece;


Leave to soak for 10 minutes at  $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ . This procedure reduces the amount of blood, protein and mucous present on the air-polishing handpiece;

29


After 10 minutes of immersion in the enzymatic solution, gently brush the surface of the air-polishing handpiece using the brush with soft nylon bristles for at least 20 seconds, and eliminate all visible residue. Thoroughly clean the hard-to-reach areas such as the edges, recesses, and joints;

30


Rinse the air-polishing handpiece under hot running water ( $40^{\circ}\text{C} \pm 5^{\circ}\text{C}$ ) for about 10 minutes, moving it slightly to allow the water to reach the entire surface;

31


Place the air-polishing handpiece in the ultrasonic tank immersed in the enzymatic detergent solution for at least 10 minutes at  $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ;

32


Gently brush the surface of the air-polishing handpiece again using the brush with soft nylon bristles;

33


Rinse the air-polishing handpiece under hot running water (40 °C ±5 °C) for about 10 minutes, moving it slightly to allow the water to reach the entire surface;

34


Place the air-polishing handpiece in the ultrasonic tank immersed in demineralised water for at least 10 minutes at 40°C ±2°C;

35

After 10 minutes of immersion in the ultrasonic tank with demineralised water, rinse the internal channels of the air-polishing handpiece by injecting 20 ml of deionised water at ambient temperature (20°C – 25°C) using a syringe with needle, five times in a row.

36


Process validated by independent bodies with enzymatic detergent.

## 8.4.1.6 K9 Wrench


Place the K9 wrench in a clean container. Add enough enzymatic detergent solution to completely cover it. Leave the K9 wrench to soak in the enzymatic detergent solution for 10 minutes at ambient temperature ( $22\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ );

37

While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed. Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.

38


**NOTE:** Thoroughly brush all the following parts of the K9 wrench for about 20 seconds:

- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- External metal ring;
- Internal cavities, grooves and fissures.

Remove the K9 wrench from the enzymatic detergent solution. Thoroughly rinse and brush all the surfaces of the K9 wrench (see previous point) under running water for at least 10 minutes;


39

EN


Place the K9 wrench in the ultrasonic tank immersed in the enzymatic detergent solution at 24 °C  $\pm$ 2 °C and run a cycle for at least 20 minutes;

40

Remove the K9 wrench from the ultrasonic tank and rinse it under running water so as to eliminate all detergent residue. Brush the internal and external surfaces of the K9 wrench with a clean brush with soft nylon bristles, under running water.

41

## 8.4.2 Automatic Cleaning

Automatic cleaning can be carried out as an alternative to the manual cleaning described in Chapter 8.4.1 on page 77.

**NOTE:** Procedure validated with:

- Miele PG8536 washer/disinfector;
- Miele DES-VAR-TD program;

- Alkaline detergent: neodisher® FA (0.2 % v/v);
- Neutralising liquid: neodisher® Z (0.1 % v/v)

### 8.4.2.1 Required Materials

- Alkaline detergent: neodisher® FA (0.2 % v/v);
- Neutralising liquid: neodisher® Z (0.1 % v/v);
- Water;
- Metal basket;
- Adapters;
- Thermal-disinfector.

**NOTE:** Make sure that the accessories are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

**⚠ WARNING:** Avoid overloading the thermal-disinfector as this could compromise the effectiveness of cleaning.


**⚠ WARNING:** On completion of the cleaning cycle in the thermal-disinfector, the scaler handpiece remains at the wash temperature for a long time. Use appropriate precautions when extracting the scaler handpiece from the thermal-disinfector to prevent injury to the operator.

**ⓘ CAUTION:** Owing to its conformation, the scaler handpiece can rotate. When not in use, the scaler handpiece must always be placed on its support.

Position the accessories in a metal basket. Connect the relevant adapter (supplied as an optional) to the scaler handpiece connector and then to the water jet cleaning connections of the thermal-disinfector.

Repeat the same operation for the inserts by connecting them to the appropriate adapters supplied as optionals.

1

Sequence and parameters applicable to the cycle:

- 1 min, Rinse with cold water;
- 5 min, Wash with alkaline detergent at  $55^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ;
- 1 min, Neutralisation with suitable solution (1/3 cold water, 2/3 hot water);
- 1 min, Rinse with water (1/3 cold water, 2/3 hot water);
- 5 min, Thermal disinfection at  $93^{\circ}\text{C}$  with demineralised water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1, Table B.1 [4] thermal disinfection at a temperature of  $90^{\circ}\text{C}$  for 5 min determines a value A0 3000.

2

## 8.5 Cleaning Check


### 8.5.1 Required Materials

- Light source;

Once the cleaning operations have been completed, check the scaler handpiece and scaler nose cone under an adequate source of light, if necessary using a magnifying glass 2.5X, paying attention to any parts that may conceal soil residue (threading, cavities, grooves) and, if need be, repeat the cleaning cycle if soil is still visible. Finally, check the integrity of those parts and those elements that could have deteriorated during use;

- 2.5X Magnifier.





1

Repeat the checking operations for the other accessories (inserts, inserts tightening wrenches, air-polishing handpieces, K9 wrench), repeating the cleaning cycle if necessary.

2

EN

## 8.6 Drying and Lubrication

### 8.6.1 Required Materials



- Compressed air;
- Soft, low-lint cloth;
- Medical grade lubricant.

Thoroughly dry all parts of the scaler handpiece, scaler nose cone, and light pipe by blowing them with compressed air;


**⚠ CAUTION:** The scaler handpiece electrical contacts must be dry at the start and end of the sterilisation cycle, before the cord is connected to the appliance. Always make sure that the electrical contacts of the connector are perfectly dry, if necessary dry them blowing with compressed air.

**⚠ CAUTION:** Before starting the sterilisation cycle, make sure that the insert is thoroughly dry both internally and externally. To do so, blow compressed air both externally and through the internal through-hole. This will prevent the appearance of stains, streaks on the surface, or oxidation inside the insert.

1


2


**⚠ CAUTION:** Before starting the sterilisation cycle, make sure that the air-polishing handpiece is thoroughly dry both internally and externally. To do so, blow compressed air both externally and through the internal through-holes. This will prevent the appearance of stains, streaks on the surface, or oxidation inside the air-polishing handpiece.

3

Dry the inserts tightening wrench using a soft, low-lint cloth;


4

Lubricate the inserts tightening wrench with medical-grade lubricants in the indicated point.


5

**⚠ CAUTION:** Do not use oil or silicone-based lubricants.

Dry the K9 wrench using a soft, low-lint cloth.

6


## 8.7 Sterilisation

### 8.7.1 Preparation


Seal the scaler handpiece (without inserts), scaler nose cone, and light pipe individually and separately in disposable sterilisation bags.

1


Seal the inserts individually inside a disposable bag for sterilisation. Seal the air-polishing handpieces individually inside a disposable bag for sterilisation.

2


Seal the wrench individually inside a disposable bag for sterilisation.

3

Seal the K9 wrench individually inside a disposable bag for sterilisation.

4

EN

## 8.7.2 Sterilisation Method

The scaler handpiece and other sterilisable accessories are manufactured with materials resistant to a maximum temperature of 135°C for a maximum time of 20 minutes.

Once the scaler handpiece and the other sterilisable accessories have been put into bags individually, perform the steam sterilisation process in the autoclave.

**⚠ CAUTION:** Use sterilisation bags compliant with standard UNI EN ISO 11607-1.

The sterilisation process validated by Mectron S.p.A., in an autoclave with steam, guarantees a SAL  $10^{-6}$  by setting the parameters indicated below:

- **Type of cycle:** 3 times Pre-vacuum (pressure min. 60 mBar).
- **Minimum sterilisation temperature:** 132°C (interval 0°C ÷ +3°C).
- **Minimum sterilisation time:** 4 minutes.
- **Minimum drying time:** 20 minutes.

All stages of sterilisation must be carried out by the operator in compliance with the most current version of standards: UNI EN ISO

17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

**⚠ CAUTION:** Do not sterilise the scaler handpiece with the insert screwed in.

**⚠ WARNING: Infection control - Sterilisable parts.** Thoroughly remove all residual organic matter before sterilisation.

**⚠ CAUTION:** Carry out the sterilisation by using only an autoclave with steam from water. Do not use any other sterilisation procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

**⚠ CAUTION:** Do not exceed the permissible load of the steam-steriliser.

**⚠ WARNING:** On completion of the sterilisation cycle in autoclave, the scaler handpiece remains at the sterilisation temperature for a long time. Use appropriate precautions when extracting the scaler handpiece from the autoclave to prevent injury to the operator.

**⚠ CAUTION:** Wait for the scaler handpiece to cool down completely before use.

### 8.7.2.1 Special Information

Sterilisation parameters, in steam autoclave, used in Great Britain:

- Temperature: 134 °C;
- Time: 3 minutes.

## 9 MAINTENANCE

If the appliance is not used for a long time observe the following recommendations:

1. Empty the powder containers;
2. Run a complete cleaning cycle of the irrigation circuit using the "flush" function (see *Chapter 6 on page 40*) for both the Ultrasound and Polishing parts;
3. Drain the irrigation circuits of any residual water, removing the tank and running the "flush" cycle for the Ultrasound part;
4. Eliminate the condensate water from the air filter (see *Chapter 9.3 on page 101*);
5. Disconnect the appliance from the power mains and from the air and water circuits;
6. Clean and dry the water filter (see *Chapter 9.2 on page 100*);
7. In case of long periods of non-use, place the appliance back in its original packaging, in a safe place;
8. Before using the appliance again, clean and sterilise the handpiece and accessories following the instructions in *Chapter 8 on page 49*;
9. Check that the inserts are not worn, deformed or broken, with particular attention to the integrity of the tip.

**⚠ WARNING:** Periodically check the integrity of the power supply cable; if damaged, replace it with an original Mectron spare part.


### 9.1 Replacing the Peristaltic Pump

There is a plastic protection on the left side of the appliance, which covers the peristaltic pump housing. Remove this protection by pressing on the sides and pulling towards yourself;


**⚠ CAUTION:** Before carrying out operations on the peristaltic pump, make sure the appliance is unplugged and that the liquids container is not connected.

Remove the two tubes of the pump from the respective clutches positioned below it;

1


2


Extract the peristaltic pump from its seat, pulling it towards yourself. Use caution as pieces may detach;

3

Attach the new peristaltic pump onto the base until you hear it click in, and connect the two pump tubes to the respective couplings underneath it;


4

**⚠ WARNING:** Only use original spare parts.

Refit the plastic protection onto the peristaltic pump.

5

## 9.2 Cleaning and/or Replacing the Water Filter


Check and clean the water filter monthly, performing the following operations:

Disconnect the water supply hose from the male coupling.

Unscrew the knurled bush of the male coupling;


**⚠ CAUTION:** Before cleaning and/or replacing the water filter, make sure that the appliance is disconnected from the mains electricity and that the liquids container is not connected.

1

2


Extract the filter and wash it under running water to eliminate any impurities obstructing it;

3

Re-introduce the filter into its seat and screw the knurled bush back into its housing tightly until it is fully home.

**NOTE:** Replace the filter with a new one if it is damaged or washing is not effective.

### 9.3 Eliminating the Condensate


The appliance has an air filter, which intercepts any impurities and the condensate present in the pneumatic circuit.

To prevent the condensate entering into circulation in the appliance, check and empty the air filter weekly by performing the following operations:

Position an absorbent cloth underneath the appliance to collect the condensate;

**⚠ CAUTION:** This maintenance operation must be performed with the appliance on in order for the air circuit to be pressurised.


1

With the appliance on and in a perfectly horizontal position, press the air filter vent valve positioned on the bottom of the appliance, until just air escapes.

**NOTE:** However, it is recommended to use dry compressors and introduce a dehumidifier into the practice's pneumatic circuit.

2


### 9.4 Cleaning the Powder Containers and Caps

Check cleanliness of the powder container and, in particular, the cap as powder residues in the presence of moisture could solidify and make opening and closing operations difficult.

**⚠ CAUTION:** Always switch the appliance off using the I/O switch and disconnect it from the mains electricity before cleaning the powder containers and the caps.

**⚠ WARNING:** Before blowing compressed air into the powder containers, make sure that they have been emptied.

1


Blow compressed air inside the container and on the threading of both the powder containers and the caps.

When the cleaning operation has been completed, re-position the containers on the appliance.

**⚠ CAUTION:** Insert the powder containers in such a way that the groove of the bottom of the appliance is aligned with the convexity of the powder container.

2

## 10 METHODS AND PRECAUTIONS FOR DISPOSAL

**⚠ WARNING:** Hospital waste. Treat the following objects as hospital waste:

- Inserts, when worn or broken;
- Inserts tightening wrench, when worn or broken;
- Air-polishing handpiece, when worn or broken;
- Cleaning needles, when worn or broken;
- K9 wrench, when worn or broken;
- SUBGINGIVAL PERIO TIP, after each use.

Disposable materials and materials that carry biological risk must be disposed of according to local regulations in force regarding hospital waste.

**⚠ WARNING:** When handling the inserts, pay particular attention to the sharp, pointed, irregular parts to avoid any wounds or injuries.

combi touch must be disposed of and treated as waste subject to separate collection.

Failure to comply with the previous points can result in a penalty pursuant to the directive on waste electrical and electronic equipment (WEEE).

The purchaser is entitled to deliver the device to be disposed of to the dealer who supplies them with new equipment; at Mectron, instructions for proper disposal are available.

## 11 TECHNICAL DATA

|   |   |
|---|---|
| <b>Device compliant with Regulation (EU) 2017/745</b>       | Class IIa   |
| <b>Classification under the IEC/EN 60601-1</b>              | I<br>Parts applied: type B (insert)<br>IP 20 (appliance)<br>IP 22 (pedal model FS-05)   |
| <b>Essential performance</b>                                | According to the standard IEC 80601-2-60 the device has no essential performance  |
| <b>Appliance for intermittent operation</b>                 | 55sec. ON - 30sec. OFF with irrigation (ULTRASOUND function and POLISHING function)<br>30sec. ON - 120sec. OFF without irrigation (mode: "endo" and "perio/scaler" power from 1 to 5, "restorative" power from 1 to 4).                 |
| <b>Power Supply</b>   | 100-240 V~ 50/60 Hz   |
| <b>Max. Power Consumption</b>                               | 90 VA   |
| <b>Fuses</b>  | Type 5 x 20 mm, T 2AL, 250V   |
| <b>Water supply:</b>  | Working pressure from 1 to 6 bar.<br>Polisher and scaler part water circuit cleaning function - See <i>Chapter 6 on page 40</i> .<br>Connection via the supplied hose with quick coupling through an incorporated and removable filter. |
| <b>Air supply:</b>  | Inlet pressure between 4 and 8 bar.<br>Air circuit cleaning function - See <i>Chapter 9.3 on page 101</i><br>Connection via the supplied hose with quick coupling through an incorporated filter and pressure reducer.                  |
| <b>Operating Conditions</b>                                 | from 10°C to +35°C<br>Relative humidity from 30% to 75%<br>Air pressure P: 800hPa/1060hPa   |
| <b>Transport and Storage Conditions (Excluding Powders)</b> | from -10°C to +60°C<br>Relative humidity from 10% to 90%<br>Air pressure P: 500hPa/1060hPa  |
| <b>Altitude</b>   | lower than or equal to 2000 metres  |
| <b>Weights and dimensions</b>                               | 4.8Kg<br>410 x 260 x 145 mm (L x W x H) <sup>a)</sup>   |

**Table 5** – Technical Data

a) W = width ; L = length ; H = height

## 11.1 Ultrasound part

|                                       |   |
|---------------------------------------|---|
| <b>Working Frequency</b>              | Automatic Scan<br>From 24 KHz to 36 KHz   |
| <b>Power Level</b>                    | "endo"<br>"perio/scaler"<br>"restorative"<br>"soft mode"  |
| <b>Irrigation</b>                     | Water circuit with continuous regulation.<br>Peristaltic pump with continuous regulation via touch screen: <ul style="list-style-type: none"> <li>• 7 flow rate levels, from 1 (approx. 5 ml/min) to 7 (approx. 30 ml/min).</li> </ul> Possibility to exclude irrigation in "endo", "perio/scaler" and "restorative" powers, with power levels from 1 to 4. |
| <b>Handpiece LED system:</b>          | ON/OFF function enabled:<br>The handpiece LED lights up as soon as the machine starts working, and turns off 3 seconds after the pedal is released.<br>ON/OFF function disabled:<br>The handpiece LED is off.<br>White LED light power risk free according to standard IEC/EN 62471.  |
| <b>Protections of the APC Circuit</b> | Missing handpiece;<br>Breaking wire cord;<br>Insert not correctly tightened or broken;  |

## 11.2 Polishing part

|                           |   |
|---------------------------|---|
| <b>Polisher functions</b> | Can be selected via touch screen:<br>"Prophy" function - "Perio" function |
| <b>Irrigation</b>         | Continuous regulation via knob.<br>Water heating via heating element.     |

## 11.3 Electromagnetic Compatibility IEC/EN 60601-1-2

**⚠ WARNING: Contraindications. Interference with other equipment**

Even if compliant with standard IEC/EN 60601-1-2, combi touch may interfere with other devices in the vicinity. combi touch must not be used near or stacked with other appliances. However, if this is necessary, the correct operation of the equipment in that configuration, and of all equipment, must be checked and monitored before starting the intervention.

**⚠ WARNING:** Portable and mobile radio communication equipment may influence the correct operation of the appliance.

**⚠ WARNING: Contraindications. Interference from other equipment**

An electrosurgical scalpel or other electrosurgical units positioned near the combi touch appliance may interfere with the correct operation of the appliance itself.

**⚠ WARNING:** The appliance requires particular EMC precautions and must be installed and put into service according to the EMC information provided in this chapter.

**⚠ WARNING:** The use of cables and accessories not supplied by MECTRON may adversely affect the EMC performances.

### 11.3.1 Guide and Manufacturer's Declaration - Electromagnetic Emissions

combi touch is designed to operate in the electromagnetic environment specified below. The purchaser or user of combi touch should ensure that it is used in such an environment.

| Emissions Test  | Compliance | Electromagnetic Environment Guidance  |
|---|------------|---|
| RF Emissions<br>CISPR 11                                      | Group 1    | combi touch uses RF energy only for its internal operation. Therefore, its RF emissions are very low and probably do not cause any interference with nearby electronic devices.                                 |
| RF Emissions<br>CISPR 11                                      | Class B    | combi touch is suitable for use in all buildings, including domestic buildings, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2                           | Class A    |   |
| Emissions of fluctuations<br>voltage/flicker<br>IEC 61000-3-3 | Compliant  |   |

### 11.3.2 Accessible Parts of the Casing

combi touch is designed to operate in the electromagnetic environment specified below. The purchaser or user of combi touch should ensure that it is used in such an environment.

| Phenomenon   | Essential EMC standard or test method | Immunity test values   | Electromagnetic Environment Guidance  |
|--|---------------------------------------|--|---|
| Electrostatic discharge (ESD)                          | IEC 61000-4-2                         | ±8 kV on contact<br>±2 kV, ±4 kV, ±8 kV,<br>±15 kV in air                              | The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.   |
| Radiated RF EM fields <sup>a)</sup>                    | IEC 61000-4-3                         | 3 V/m <sup>f)</sup><br>80 MHz - 2,7 GHz <sup>b)</sup><br>80% AM at 1 kHz <sup>c)</sup> | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter. |
| Magnetic field at the mains frequency <sup>d) e)</sup> | IEC 61000-4-8                         | 30 A/m <sup>g)</sup><br>50 Hz or 60 Hz   | Power supply frequency magnetic fields should have levels characteristic of a typical location in a commercial or hospital environment.   |

- a) If used, the interface between the PATIENT'S physiological signal simulation and combi touch , must be positioned within a 0.1 m radius of the vertical plane of the uniform field area in the same direction as combi touch .
- b) combi touch, which intentionally receives RF electromagnetic energy for its operation, must be tested at the receiving frequency. The test can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test evaluates the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an environmental signal is in the passband. It is understood that the receiver may not receive normally during the test.
- c) The test can be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

- d) Only applicable to equipment and systems with magnetically sensitive components or circuits.
- e) During the tests, combi touch can be supplied with any NOMINAL input voltage, but with the same frequency as the test signal.
- f) Before applying the modulation.
- g) This test value assumes a minimum distance between combi touch and the magnetic field sources with a supply frequency of at least 15 cm. If the RISK ANALYSIS indicates that combi touch will be used less than 15 cm from the sources of the magnetic field with power frequency, the immunity test value should be adjusted according to the minimum expected distance.

EN

### 11.3.3 Guide and the Manufacturer's Declaration - Electromagnetic Immunity

#### 11.3.3.1 Power Connection BC Input

combi touch is designed to operate in the electromagnetic environment specified below. The purchaser or user of combi touch should ensure that it is used in such an environment.

| Phenomenon   | Essential EMC standard or test method | Immunity test values  | Electromagnetic Environment Guidance  |
|--|---------------------------------------|---|---|
| Electrical fast transient/burst<br>a) l) o)              | IEC 61000-4-4                         | ±2 kV on contact<br>100 KHz repetition frequency  | The quality of the network voltage should be that of a typical commercial or hospital environment.  |
| Pulses Differential mode<br>a) b) j) o)                  | IEC 61000-4-5                         | ± 0.5 kV, ± 1 kV  | The quality of the network voltage should be that of a typical commercial or hospital environment.  |
| Pulses common mode<br>a) b) j) k) o)                     | IEC 61000-4-5                         | ± 0.5 kV, ± 1 kV, ± 2kV   | The quality of the network voltage should be that of a typical commercial or hospital environment.  |
| Conductive disturbances induced by RF fields<br>c) d) o) | IEC 61000-4-6                         | 3 V <sup>m)</sup><br>0.15 MHz - 80 MHz<br>6 V <sup>m)</sup> in the ISM bands between 0.15 MHz and 80 MHz <sup>n)</sup><br>80 % AM at 1 KHz <sup>e)</sup>  | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter. |
| Voltage dips<br>f) p) r)                                 | IEC 61000-4-11                        | 0% UT; 0,5 cycle g)<br>at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % UT; 1 cycle e<br>70 % UT; 25/30 cycle <sup>h)</sup><br>Single phase: at 0° | The quality of the network voltage should be that of a typical commercial or hospital environment.  |
| Voltage interruptions<br>f) i) o) r)                     | IEC 61000-4-11                        | 0 % UT; 250/300 cycle <sup>h)</sup>   | The quality of the network voltage should be that of a typical commercial or hospital environment.  |

- a) The test can be performed at any supply voltage within the combi touch range of NOMINAL voltages. If combi touch is tested at one supply voltage value, it need not be retested at other voltage values.
- b) During the test, all combi touch cables must be connected.
- c) The calibration of the current injection terminals must be performed in a system at 150  $\Omega$ .
- d) If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) The test can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) Devices and systems with a Direct Current (DC) input power supply using AC to DC converters must be tested with a converter that complies with the MANUFACTURER'S specifications. Immunity test levels are applied to the AC power input of the converter.
- g) Only applicable to devices and systems connected to a single-phase Alternating Current (AC) power supply.
- h) For example, 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) Devices and systems with nominal input current above 16 A / phase must be disconnected from the power supply once every 250/300 cycles at any angle and from all phases simultaneously (if applicable). Devices and systems with battery backup, after the test, must resume operation using the power supply line. For devices and systems with a nominal input current lower than 16 A, all phases must be disconnected simultaneously.
- j) Devices and systems that do not have a surge protection device in the primary power circuit can only be tested at  $\pm 2$  kV between the line(s) and the ground (common mode) and at  $\pm 1$  kV between line(s) and line(s) (differential mode).
- k) Not applicable to CLASS II devices and systems.
- l) It is necessary to use direct coupling.
- m) R.M.S. , applied before modulation.
- n) The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.
- o) Applicable to devices and systems with NOMINAL input current less than or equal to 16 A / phase and devices and systems with NOMINAL input current greater than 16 A / phase.
- p) Applicable to devices and systems with NOMINAL input current less than or equal to 16 A / phase.
- q) At some phase angles, the application of this test to devices with a transformer on the input power supply may cause the opening of an overcurrent protection device. This can occur due to the saturation of the magnetic flow of the transformer core after the voltage drop. If this happens, the device must ensure BASIC SAFETY during and after the test.
- r) For devices and systems that have multiple voltage settings or self-regulating voltage capabilities, the test must be performed at the minimum and at the maximum NOMINAL input voltage. Devices and systems with a NOMINAL input voltage range lower than 25% of the highest NOMINAL input voltage have to be tested with a NOMINAL input voltage within the range.

**11.3.3.2 Points of Contact with the Patient**

combi touch is designed to operate in the electromagnetic environment specified below. The purchaser or user of combi touch should ensure that it is used in such an environment.

| Phenomenon   | Essential EMC standard or test method | Immunity test values   | Electromagnetic Environment Guidance   |
|--|---------------------------------------|--|--|
| Electrostatic discharges (ESD) <sup>c)</sup>               | IEC 61000-4-2                         | ±8 kV on contact<br>±2 kV, ±4 kV, ±8 kV,<br>±15 kV in air  | The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Conductive disturbances induced by RF fields <sup>a)</sup> | IEC 61000-4-6                         | 3 V <sup>b)</sup><br>0.15 MHz - 80 MHz<br>6 V <sup>b)</sup> in ISM bands<br>between 0.15 MHz and<br>80 MHz<br>80 % AM at 1 KHz | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter. |

a) The following applies:

- All connection cables with the patient must be tested, either individually or grouped together.
- The connection cables with the patient must be tested using a current clamp unless the current clamp is not suitable. If a current clamp is not suitable, an EM clamp must be used.
- In any case, no intentional decoupling device should be used between the injection site and POINT OF CONNECTION TO THE PATIENT.
- The tests can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- The tubes that are intentionally filled with conductive liquids and intended to be placed in contact with the PATIENT must be considered connection cables with the patient.

- If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

b) R.M.S., applied before modulation.

c) The discharges must be applied without connection to an artificial hand and without connection to the PATIENT simulation. The PATIENT simulation can be connected after the test, if necessary, to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

## 11.3.3.3 Parts Accessible to the Input / Output Signals

combi touch is designed to operate in the electromagnetic environment specified below. The purchaser or user of combi touch should ensure that it is used in such an environment.

| Phenomenon   | Essential EMC standard or test method | Immunity test values  | Electromagnetic Environment Guidance   |
|--|---------------------------------------|---|--|
| Electrostatic discharges (ESD) <sup>e)</sup>                     | IEC 61000-4-2                         | ±8 kV on contact<br>±2 kV, ±4 kV, ±8 kV,<br>±15 kV in air   | The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Electrical fast transient/burst <sup>b) f)</sup>                 | IEC 61000-4-4                         | ±1 kV on contact<br>100 KHz repetition frequency  | The quality of the network voltage should be that of a typical commercial or hospital environment.   |
| Pulses Common mode <sup>a)</sup>                                 | IEC 61000-4-5                         | ± 2kV   | The quality of the network voltage should be that of a typical commercial or hospital environment.   |
| Conductive disturbances induced by RF fields <sup>b) d) g)</sup> | IEC 61000-4-6                         | 3 V <sup>h)</sup><br>0.15 MHz - 80 MHz<br>6 V <sup>h)</sup> in the ISM bands between 0.15 MHz and 80 MHz <sup>i)</sup><br>80% AM at 1 KHz <sup>c)</sup> | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter. |

- a) This test applies only to output lines connected directly to the external cables.
- b) SIP/SOPS with maximum cable length less than 3 m are excluded.
- c) The tests can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) The calibration of the current injection terminals must be performed in a system at 150 Ω.
- e) The connectors have to be tested in accordance with paragraph 8.3.2 and Table 4 of the standard IEC 61000-4-2:2008. For insulated connector housings, perform the air discharge test on the connector housing and pins using the probe with the rounded tip of the ESD generator, with the exception that only the connector pins that are tested are those that can be reached or touched, under the conditions of INTENDED USE, by the standard probe shown in Figure 6 of the general standard, applied in a bent or straight position.
- f) It is necessary to use the capacitive coupling.
- g) If an ISM or amateur radio band is not present among

the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

- h) R.M.S., applied before modulation.
- i) The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

### 11.3.4 Specifications of the tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device

combi touch is designed to operate in an electromagnetic environment in which radiated RF disturbances are under control. The purchaser or operator of combi touch can help prevent electromagnetic interferences by guaranteeing a minimum distance between the mobile and portable RF communication devices (transmitters) and combi touch, as recommended below, in relation to the maximum output power of the radio communication devices.

| Test Freq. (MHz) | Band <sup>a)</sup> (MHz) | Service <sup>a)</sup>  | Modulation <sup>b)</sup>                            | Max power (W) | Distance (m) | Immunity test value (V/m) |
|------------------|--------------------------|--|---|---------------|--------------|---------------------------|
| 385              | 380 - 390                | TETRA 400  | Pulse modulation <sup>b)</sup><br>18 Hz             | 1.8           | 0.3          | 27                        |
| 450              | 430 - 470                | GMRS 460<br>FRS 460  | FM <sup>c)</sup><br>± 5 kHz deviation<br>1 kHz sine | 2             | 0.3          | 28                        |
| 710              | 704 - 787                | LTE band 13,<br>17   | Pulse modulation <sup>b)</sup><br>217 Hz            | 0.2           | 0.3          | 9                         |
| 745              |                          |  |   |               |              |                           |
| 780              |                          |  |   |               |              |                           |
| 810              | 800 - 960                | GSM 800/900<br>TETRA 800<br>iDEN 820<br>CDMA 850<br>Band LTE 5               | Pulse modulation <sup>b)</sup><br>18 Hz             | 2             | 0.3          | 28                        |
| 870              |                          |  |   |               |              |                           |
| 930              |                          |  |   |               |              |                           |
| 1720             | 1700 - 1990              | GSM 1800<br>CDMA 1900<br>GSM 1900<br>DECT<br>LTE Band 1, 3,<br>4, 25<br>UMTS | Pulse modulation <sup>b)</sup><br>217 Hz            | 2             | 0.3          | 28                        |
| 1845             |                          |  |   |               |              |                           |
| 1970             |                          |  |   |               |              |                           |
| 2450             | 2400 - 2750              | Bluetooth<br>WLAN<br>802.11 b/g/n<br>RFID 2450<br>Band LTE 7                 | Pulse modulation <sup>b)</sup><br>217 Hz            | 2             | 0.3          | 28                        |
| 5420             | 5100 - 5800              | WLAN<br>802.11 a/n   | Pulse modulation <sup>b)</sup><br>217 Hz            | 0.2           | 0.3          | 9                         |
| 5500             |                          |  |   |               |              |                           |
| 5785             |                          |  |   |               |              |                           |

a) For some services, only uplink frequencies are included.

b) The carrier must be modulated using a square wave signal with a duty cycle of 50%.

c) As an alternative to FM modulation, 18 Hz pulse modulation at 50% can be used because, although it does not represent true modulation, it represents the worst case.

**NOTE:** If necessary to reach the immunity test level, the distance between the transmitter antenna and combi touch can be reduced to 1 m. The test distance of 1 m is allowed by IEC 61000-4-3.

**⚠ WARNING:** Portable RF communication equipment (including peripheral devices such as antenna cables and external antennas) must not be used closer than 30 cm to any part of the combi touch device, including the cables specified by the manufacturer. Otherwise, there may be a performance degradation of these appliances.

## 12 TROUBLESHOOTING

### 12.1 Diagnostic System and Symbols on the Keyboard

combi touch is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch keyboard by means of a symbol. Users, by using the following table, are guided to identifying and the possible resolution of the malfunction detected.

| Symbol on keyboard   | Possible cause   | Solution  |
|--|--|---|
|   | Handpiece/cord electrical contacts wet                               | Thoroughly dry the contacts with compressed air (see <i>Chapter 8.6 on page 94</i> ).                           |
|  | Scaler handpiece combi touch not connected to the appliance          | Connect the scaler handpiece (see <i>Chapter 4.3 on page 18</i> ).  |
|  | Scaler handpiece failure   | Replace the scaler handpiece  |
|  | Tuning circuit malfunction   | Contact an Authorised Mectron Service Centre.   |
|  | Insert not present or not properly tightened on the scaler handpiece | Unscrew the insert and screw it back in correctly using the torque wrench (See <i>Chapter 5.4 on page 32</i> ). |
|  | Insert broken, worn or deformed                                      | Replace the insert.   |
|  | Cord/handpiece electrical contacts wet                               | Thoroughly dry the contacts with compressed air (see <i>Chapter 8.6 on page 94</i> ).                           |

| Symbol on keyboard  | Possible cause   | Solution   |
|---|--|--|
|  | Incorrect start-up procedure: the appliance was started with the pedal pressed.                    | Check that the pedal is not pressed down. If the problem persists, disconnect the pedal and potentially contact an Authorised Mectron Service Centre.                                    |
|   | Peristaltic pump malfunction   | Check that there are no impediments to the rotation of the peristaltic pump. Check that the peristaltic pump and two tubes are correctly installed (see <i>Chapter 9.1 on page 98</i> ). |
|   | The appliance has been switched off and on again without waiting 5 seconds                         | Turn off and wait for 5 seconds before switching on the appliance again.   |
|   | Faults on the electrical network, excessive electrostatic discharges or internal faults            | Turn off and wait 5 seconds before switching the appliance on again<br>If the warning persists, contact an Authorised Mectron Service Centre.  |
|   | A powder container has been opened without performing the "refill" cycle.                          | Before opening one of the powder containers, the "refill" cycle must be run (See <i>Chapter 5.6.1 on page 38</i> ).  |
|   | One of the powder containers has been removed from its seat without performing the "refill" cycle. | Before removing one of the powder containers, the "refill" cycle must be run (See <i>Chapter 5.6.1 on page 38</i> ).   |
|   | The selected powder container is not correctly inserted in place.                                  | Correctly insert the powder containers until they are in contact.  |

**Table 6** – Error Messages

**NOTE:** For diagnostic warnings not included in this list, contact the technical service desk.

## 12.2 Quick Troubleshooting

| Problem   | Possible Cause   | Solution  |
|---|--|---|
| The appliance does not start after the switch has been moved to the "I" position.   | The terminal of the electrical power supply cable has not been properly inserted in the rear socket of the appliance | Check that power supply cable is firmly connected   |
|   | The power supply cable is faulty   | Check that the supply socket is working.<br>Replace the power supply cable  |
|   | The fuses are out of order   | Replace the fuses (See <i>Chapter 12.3 on page 119</i> )  |
| The appliance is on but is not working. There are no anomalies signalled on the keyboard.   | The pedal plug is not correctly inserted into the appliance socket   | Properly insert the plug of the pedal in the connector on the back of the appliance (see <i>Chapter 4.3 on page 18</i> ). |
|   | The pedal does not work properly   | Contact an Authorised Mectron Service Centre  |
| The appliance is on but is not working.<br>One of the following symbols appears on the screen:<br> | See <i>Chapter 12.1 on page 112</i> for the possible cause depending on the symbol                                   | See <i>Chapter 12.1 on page 112</i> for the action to be taken depending on the symbol                                    |
| During operation a faint whistling noise can be heard coming from the scaler handpiece.   | The insert is not correctly tightened on the handpiece   | Unscrew the insert and screw it back in correctly using the Mectron torque wrench (See <i>Chapter 5.4 on page 32</i> )    |
|   | The irrigation circuit has not been completely filled  | Fill the irrigation circuit using the "flush" function (See <i>Chapter 5.4 on page 32</i> )                               |
| The pump rotates properly, but liquid leaks from the handpiece when it is stopped.  | The peristaltic pump is worn.  | Replace the peristaltic pump (See <i>Chapter 9.1 on page 98</i> ).  |
| By pressing the pedal, a prolonged signal is emitted and the ULTRASOUND and POLISHING function LEDs flash.  | The pedal has been pressed with both handpieces stowed in their housings.  | Lift the handpiece to be used before pressing the pedal.  |

| <b>Problem</b>   | <b>Possible Cause</b>                                    | <b>Solution</b>  |
|------------------|--|--|
| Poor performance | The insert is not correctly tightened on the handpiece   | Unscrew the insert and screw it back in correctly using the Mectron torque wrench (See <i>Chapter 5.4 on page 32</i> ) |
|                  | Insert broken, worn or deformed                          | Replace the insert with a new one  |
|                  | Insufficient or excessive powder level in the container. | Restore the correct level of powder in the container (See <i>Chapter 4 on page 17</i> ).                               |

| Problem   | Possible Cause  | Solution  |
|---|---|---|
| No liquid flows out from the insert or air-polishing handpiece during operation | Appliance not connected to the water circuit                      | Check the connection to the water circuit (See <i>Chapter 4.3 on page 18</i> ).   |
|   | The insert is the "Dry Work" type without a liquid passage        | Use an insert type with a liquid passage  |
|   | The insert is clogged   | Unscrew the insert/air-polishing handpiece from the handpiece and free the water passage of the insert/air-polishing handpiece by blowing compressed air through it.<br>If the problem persists, replace the insert/air-polishing handpiece with a new one. |
|   | The quick coupling on the air-polishing handpiece cord is clogged | Contact an Authorised Mectron Service Centre  |
|   | Irrigation is disabled  | Activate irrigation and adjust the irrigation level (See <i>Chapter 5.2.1 on page 24</i> ).   |
|   | The water tap on the appliance is closed                          | Regulate the flow of water using the knob dedicated to the function used.   |
|   | Water filter clogged  | See <i>Chapter 9.2 on page 100</i> .  |
|   | The liquid tank is empty  | Fill the tank   |
|   | The tank is not properly installed                                | Properly connect the tank to the machine core unit  |
|   | The silicone tubes of the pump are not correctly installed        | Check the tube connections (see <i>Chapter 9.1 on page 98</i> ).  |
| The peristaltic pump is worn  | Replace the peristaltic pump (See <i>Chapter 9.1 on page 98</i> ) |   |

| <b>Problem</b>   | <b>Possible Cause</b>   | <b>Solution</b>   |
|--|---|---|
| No powder flows out from the air-polishing handpiece during operation. | Appliance not connected to the air circuit  | Check the connection to the air circuit (See Chapter 4.3 on page 18).   |
|  | Air-polishing handpiece clogged due to an excessive amount of moisture in the powder or insufficient cleaning/maintenance | See Chapter 9 on page 98  |
|  | Polisher handpiece channel clogged due to an excessive amount of moisture in the powder or insufficient cleaning.         | See Chapter 8.4 on page 77  |
|  | The amount of powder in the container exceeds the maximum level   | Remove the powder from the container and clean it using a dry cloth. Restore the correct level of powder in the container (See Chapter 4 on page 17). |
|  | Unsuitable powder   | Use the correct powder for correct appliance operation.   |
| Powder leak through the powder container cap                           | Cap not correctly tightened   | Screw the cap correctly.  |
|  | Powder residue in the thread  | Clean the powder container thread (See Chapter 9.4 on page 101).  |
| Poor cleaning performance  | Insufficient pressure in the air supply circuit   | Check air supply circuit pressure (4-8 bar max).  |
|  | Insufficient or excessive powder level in the container   | Correct level of powder in the container.   |
|  | Unsuitable powder   | Check the correct powder for correct appliance operation.   |
|  | Air-polishing handpiece clogged due to an excessive amount of moisture in the powder or insufficient cleaning/maintenance | Remove the powder from the container and clean it using a dry cloth.  |

| Problem   | Possible Cause   | Solution   |
|---|--|--|
| One of the powder container caps does not unscrew | The appliance is on and the powder container is pressurised  | Run the "refill" cycle before opening one of the powder containers (See <i>Chapter 5.6.1 on page 38</i> ). |
|   | The "refill" cycle has been performed but the powder containers have remained pressurised because the air-polishing handpiece is clogged | Read the section on cleaning the air-polishing handpiece (See <i>Chapter 8 on page 49</i> ).               |
|   | The "refill" cycle has been performed but the powder containers have remained pressurised because the polisher cord is clogged           | Contact an Authorised Mectron Service Centre.  |


**Table 7** – Quick Troubleshooting

## 12.3 Replacing the Fuses

**⚠ WARNING: Switch the appliance off.**  
Always switch the appliance off using the main switch and disconnect it from the power supply socket before carrying out the next intervention.

Pry with a flat screwdriver, inserting the tip into the seat of the fuse holder compartment located under the power socket;


1

Extract the fuse holder compartment;


**⚠ WARNING: Replace the fuses in respect of the characteristics indicated in Chapter 8 on page 49.**

2

Reinsert the compartment in its place.

3

## 12.4 Shipping to an Authorised Mectron Service Centre

If technical assistance is required on the machine, contact an Authorised Mectron Service Centre or your Dealer. Do not try to repair or modify the appliance and its accessories.

Clean and sterilise all parts that need to be sent to an Authorised Mectron Service Centre, following the instructions in *Chapter 8 on page 49*.

Leave the sterilised parts in the bag which certifies the sterilisation process.

The cleaning and sterilisation demands comply with the mandatory requirements for workplace health and safety protection pursuant to Legislative Decree 81/08 and Italian laws.

If the client fails to comply with the requirements, Mectron reserves the right to charge them the cost of cleaning and sterilisation or to refuse goods that arrive in unsuitable conditions, returning them at the client's own expense for proper cleaning and sterilisation.

The appliance must be returned suitably packaged and accompanied by all the accessories and a form including:

- Details of owner and contact number;
- Product name;
- Serial number and/or batch number;
- Reason for return / description of failure;
- Photocopy of the receipt or invoice for the purchase of the appliance.

### ⚠ CAUTION: Packaging

Pack the appliance in its original packaging to prevent damage during transport.

Once the material is received at the Authorised Mectron Service Centre, the qualified technical personnel will evaluate the problem. The repair will be made only upon acceptance by the end client. For further details contact your nearest Authorised Mectron Service Centre or your dealer.

Unauthorised repairs can damage the system and void the guarantee and furthermore will disclaim Mectron from any liability for direct or indirect damage to persons or property.

## 13 WARRANTY

All Mectron devices, before being marketed, are subjected to a careful final inspection that verifies their full functionality.

Mectron provides a warranty for combi touch , purchased new from a Mectron dealer or importer, which covers defects in material and workmanship for a period of:

- 2 (TWO) YEARS for the appliance from the date of purchase;
- 1 (ONE) YEAR for the handpiece from the date of purchase.

The other accessories are not covered by the warranty.

During the warranty period, Mectron undertakes to repair (or, at its discretion, to replace) the parts of products free of charge, which, according to its judgement, are proven to be defective.

Full replacement of Mectron products is not covered by the warranty.

The manufacturer's warranty and the approval of the appliance are not valid in the following cases:

- The appliance is not used in accordance with the intended use.
- The appliance is not used in accordance with all the instructions and requirements described in this manual.
- The electrical system of the places where the appliance is used does not comply with the laws in force and relative requirements.
- Assembly operations, extensions, adjustments, updates and repairs are carried out by personnel not authorised by Mectron.
- The environmental conditions for preservation and storage of the device do not comply with the requirements indicated in *Chapter 8 on page 49*.

- Use of non-original Mectron inserts, accessories, and spare parts that may compromise the correct operation of the appliance and cause injury to the patient.
- Accidental breakage during transport.
- Damage due to incorrect use or carelessness, or due to connection to a voltage other than that envisaged.
- Expired warranty.

The expected service life of the device is minimum 5 years.

The service life / duration does not establish a limit of use; the service life of the device defines the period of time after installation and/or commissioning, during which the original performance, or in any case performance suited to the intended use, is guaranteed without there being any degradation such as to compromise functionality and reliability.

The service life is a minimum qualitative objective of the design, therefore, individual parts or components may guarantee superior performance and reliability with respect to that declared by the manufacturer.

The service life assumes compliance with the maintenance schedules set out in this manual, does not include components normally subject to "wear", and is not linked to the warranty period: the service life does not establish any implicit or explicit extension of the warranty period.

### **CAUTION**

The warranty starts from the date of purchase of the device, which evidence is given by the delivery note/purchase invoice issued by the Dealer / Importer or, in case of device with activation code, from the date of activation of the same.

In order to avail of the warranty service, the client must return the appliance, at their own expense, to the MECTRON dealer/importer from which they purchased the product.

The appliance must be returned together with the original packaging, accompanied by all the accessories and by a form containing:

- The data of the owner and telephone number;
- The data of the Dealer / Importer;
- Photocopy of the delivery note/purchase invoice of the appliance by the owner where are reported the date, the name of the appliance and the serial number;
- Description of the failure.

The transport and the damage caused by transport are not covered by the warranty.

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© Mectron S.p.A. 2021. Tutti i diritti riservati. Nessuna parte di questo documento può essere riprodotta, in qualsiasi forma, senza il consenso scritto del proprietario del copyright.