USE AND MAINTENANCE MANUAL

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EN



FR

ES

multipiezo

CE
0051

USE AND MAINTENANCE MANUAL

EN

multipiezo

multipiezo

EN

Copyright

© Mectron S.p.A. 2021. All rights reserved. No part of this document may be reproduced, in any form whatsoever, without the prior written consent of the copyright holder.

EN

CONTENTS

1	Introduction	1
1.1	Intended Use	1
1.2	Description of the Appliance	2
1.2.1	Eligible Patient Group	2
1.2.2	Patient Selection Criteria	2
1.2.3	Indications for Use	2
1.2.4	Users	3
1.3	Disclaimer	3
1.4	Safety Requirements	4
1.5	Symbols	6
2	Identification Data	7
2.1	Appliance ID Plate	7
2.2	Handpiece Identification Data	8
2.3	Inserts Identification Data	8
3	Delivery	9
3.1	List of Components	9
4	Installation	13
4.1	First Installation	13
4.2	Safety Requirements During Installation	13
4.3	Connecting the Accessories	14
5	Use	16
5.1	Power On/Off	16
5.2	Description of the Keyboard	17
5.3	Safety Requirements Before and During Use	21
5.4	Instructions for Use	23
5.5	Important information on the Inserts	26
6	Pump/Flush Function	27
7	Disassembling Parts for Cleaning and Sterilisation	28
8	Cleaning and Sterilisation	31
8.1	Preparation	32
8.1.1	Disinfecting the Tank Irrigation Circuit	32
8.1.1.1	Required Materials	32
8.1.1.2	Procedure	32
8.2	Cleaning the Non-Sterilisable Parts	38
8.2.1	Required Materials	38
8.2.2	Cleaning Method	38
8.3	Cleaning the Tank and Cap	39
8.3.1	Preparation	39
8.3.2	Required Materials	39
8.3.3	Cleaning Method	40
8.4	Cleaning the Sterilisable Accessories	41
8.4.1	Manual Cleaning	42
8.4.1.1	Required Materials	42
8.4.1.2	Scaler Handpiece	43
8.4.1.3	Inserts	46
8.4.1.4	Inserts Tightening Wrench	48
8.4.2	Automatic Cleaning	50
8.4.2.1	Required Materials	50
8.5	Cleaning Check	52
8.5.1	Required Materials	52

8.6	Drying and Lubrication	53
8.6.1	Required Materials	53
8.7	Sterilisation	54
8.7.1	Preparation	54
8.7.2	Sterilisation Method	55
8.7.2.1	Special Information	56
9	Maintenance	56
9.1	Replacing the Peristaltic Pump	56
10	Methods and Precautions for Disposal	58
11	Technical Data	58
11.1	Electromagnetic Compatibility IEC/EN 60601-1-2	60
11.1.1	Guide and Manufacturer's Declaration - Electromagnetic Emissions	60
11.1.2	Accessible Parts of the Casing	61
11.1.3	Guide and the Manufacturer's Declaration - Electromagnetic Immunity	62
11.1.3.1	Power Connection BC Input	62
11.1.3.2	Points of Contact with the Patient	64
11.1.3.3	Parts Accessible to the Input / Output Signals	65
11.1.4	Specifications of the Tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device	66
12	Troubleshooting	67
12.1	Diagnostic System and Symbols on the Keyboard	67
12.2	Quick Troubleshooting	69
12.3	Replacing the Fuses	71
12.4	Shipping to an Authorised Mectron Service Centre	72
13	Warranty	73

PAGE LEFT INTENTIONALLY BLANK

1 INTRODUCTION

Carefully read this manual before proceeding with installation operations, use, maintenance or other intervention operations on the appliance.

This manual must always be available to the operator.

Important: To avoid damages to persons or property, carefully read all the "Safety Requirements" in the manual.

Depending on the level of severity, the safety requirements are classified with the following indications:

 **WARNING:** (Always referred to damage to persons)

 **CAUTION:** (Referring to possible damage to property)

The aim of this manual is to make the operators aware of the safety regulations,

installation procedures, instructions for proper use and maintenance of the appliance and its accessories.

Use of this manual for aims other than those strictly linked to the installation, use and maintenance of the appliance is prohibited.

The information and illustrations in this manual were updated on the edition date shown on the last page.

MECTRON is engaged in continuously updating its products with possible changes to the components of the apparatus.

In case you encounter discrepancies between the descriptions found in this manual and the equipment in your possession you can:

- check for any available updates in the MANUALS section of the MECTRON website¹;
- ask clarifications to Your Dealer;
- contact MECTRON After Sales Service.

1.1 Intended Use

With the appropriate inserts it is possible to perform the following treatments:

- **Scaling:** all the procedures for the removal of deposits of bacterial plaque and supragingival, subgingival and inter dental stones, and removal of stains;
- **Periodontics:** periodontal therapy for scaling and root-planing/debridement, including cleaning and irrigation of the periodontal pocket;
- **Cleaning treatment of the implant surface:**
- **Endodontics:** all treatments for canal preparation, irrigation, filling, gutta-percha condensation, endodontic reprocessing and retrograde cavity preparation;
- **Restorative and prosthetics:** cavity preparation and removal of carious tissue, removal of prosthetics and excess restoration materials, amalgam condensation, finishing of prosthetic stump;
- **Extraction techniques.**

The Mectron piezoelectric ultrasonic scaler devices can be used on patients of any age and gender requiring a tooth polishing treatment. There are no contraindications for specific population segments.

¹ <https://manuals.mectron.com/>

⚠️ WARNING: The appliance must be used in a dental practice or clinic, or in professional oral hygiene and preventive care centres. Do not use the device in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

⚠️ WARNING: Qualified and specialised staff.

The appliance is intended exclusively for use by a dental surgeon, dentist, or dental hygienist, with suitable medical/health training.

⚠️ WARNING: Intended use. Use the appliance exclusively for the use for which it is intended. Failure to observe this regulation can cause serious injury to the patient, the operator, and damage/failure of the device.

1.2 Description of the Appliance

multipiezo is a multi-purpose piezoelectric ultrasonic scaler designed to equip the operator with a product featuring an innovative design and unique technical characteristics, and to ensure utmost patient comfort during treatment.

The user interface is optimised to make all functions immediately available, by integrating them into the touch keyboard. The handpiece

is fitted with a 360° adjustable LED light and can be autoclaved at a maximum temperature of 135°C for a maximum period of 20 minutes. The appliance is equipped with an automatic tuning circuit that optimises power and frequency for each insert available, in order to always operate in conditions of maximum efficiency.

1.2.1 Eligible Patient Group

This medical device is designed to be used with the following patient groups:

- Children;
- Adolescents;
- Adults;
- Seniors.

This medical device can be used on patients of any age, weight, height, gender, and nationality, where applicable.

1.2.2 Patient Selection Criteria

It is not recommended to use the device in the following cases:

1. Patients with active implantable medical devices (for example: pacemakers, hearing aids, and/or other electromagnetic prostheses) without prior authorisation from their doctor;

2. Patients with clinical conditions not suitable for treatment of the sites (for example: local anaesthesia).

All piezoelectric ultrasonic scaler device models are intended for professional use only. The user is therefore the only person able to decide if and how to treat their patients.

1.2.3 Indications for Use

Use of the device is suitable for all eligible patients (see Chapter 1.2.1 on page 2) whose doctor has prescribed professional scaling, periodontology, an extraction technique, endodontic cleaning, implant site surface cleaning, a restorative or prosthetic technique, as part of the intended use of the device (see Chapter 1.1 on page 1).

1.2.4 Users

The device must be used exclusively by specialised and properly trained personnel, specifically the dental surgeon/dentist or dental hygienist, who must be able-bodied, adult, of any weight, age, height, gender, and nationality.

1.3 Disclaimer

The manufacturer MECTRON disclaims all liability, express or implied, and cannot be held responsible for direct or indirect personal injury and/or property damage, occurring as a result of incorrect procedures linked to the use of the appliance and its accessories.

The manufacturer MECTRON cannot be held liable, expressly or implicitly, for any type of personal injury and/or property damage, inflicted by the user of the product and its accessories and which occurs in the following cases, by way of example and not of limitation:

- Misuse or use during procedures other than those specified in the destination of use of the product;
- The environmental conditions for preservation and storage of the device are not complying with the requirements indicated in *Chapter 11 on page 58*;
- The appliance is not used in accordance with all the instructions and requirements described in this manual;
- The electrical system of the places where the equipment is used do not comply with the laws in force and the related regulations;
- Assembly operations, extensions, re-adjustments, upgrades and repairs of the device are carried out by personnel not authorised by MECTRON;
- Misuse, abuse, abnormal use, negligent use, intentional misconduct or use exceeding the limits of the appliance indicated and allowed and/or normal wear or deterioration, ill-treatment and/or incorrect interventions;
- Any attempt to tamper with or modification of the appliance under every circumstance;
- Use of non original MECTRON inserts resulting in permanent damage to the thread of the handpiece with impaired functioning and risk of injury to the patient;
- Use of non original MECTRON inserts used in accordance with the settings designed and tested on the original MECTRON inserts. Correct use of the settings is only guaranteed with original MECTRON inserts;
- Shortage of stock material (handpiece, inserts, wrenches) to be used in case of a standstill due to faults or problems;
- Incorrect/omitted maintenance compared to what is stated in *Chapter 9 on page 56* of this manual;
- Breach of the requirements and the information contained in *Chapter 5.5 on page 26* of this manual;
- Breach of the requirements and the information contained in *Chapter 8 on page 31* of this manual;
- Unauthorised repairs in accordance with the indications contained in *Chapter 12.4 on page 72* this manual.

1.4 Safety Requirements

EN

⚠ WARNING: Contraindications.

Do not use the appliance on patients with Pacemakers or other implantable electronic devices. This regulation also stands for the operator.

⚠ WARNING: Contraindications. Do not perform scaling treatments without water spray in order to avoid an overheating of the insert that can cause damage to the tooth. Treatments without water spray can only be those performed with the "Dry Work" inserts without water passage.

❗ CAUTION: Contraindications. Ultrasonic scaler. Do not perform treatments on metal or ceramic prosthetics. Ultrasonic vibrations could lead to decementation of the dentures.

⚠ WARNING: Contraindications.

Interference from other equipment.

An electrosurgical scalpel or other electrosurgical units positioned near the multipiezo appliance may interfere with the correct operation of the appliance itself.

⚠ WARNING: Contraindications.

Interference with other equipment. Even if compliant with standard IEC 60601-1-2, multipiezo may interfere with other devices in the vicinity. multipiezo must not be used near or stacked with other appliances. However, if this should become necessary, it is necessary to check and monitor the correct functioning of the appliance in that configuration.

⚠ WARNING: Risk of explosion. The appliance cannot operate in environments where the atmosphere is saturated with flammable gases (aesthetic mixtures, oxygen, etc.).

❗ CAUTION: In case the final user, operating in their own medical room or surgery, in order to comply with mandatory requirements, must periodically inspect the equipment present in the surgery, the test procedures to apply to medical electrical equipment and medical electrical systems for the safety assessment must be carried out following the standard EN 62353 'Medical electrical equipment - Periodic inspections and tests to be carried out after repair of medical electrical equipment'. The frequency of periodic inspections in the intended conditions of use described in this "Use and Maintenance" manual is once per year or every 2000 hours of use, whichever condition is satisfied first.

⚠ WARNING: Checking the condition of the device before treatment.

Always check that there is no water under the appliance. Before every treatment always check the perfect operation of the appliance and efficiency of the accessories. If operation anomalies are observed do not carry out the treatment. If the faults concern the appliance, contact an Authorised MECTRON Service Centre.

❗ CAUTION: The electrical system of the premises in which the appliance is installed and used must comply with the rules in force and the relevant requirements of electrical safety.

❗ CAUTION: To avoid the risk of electrical shocks, this appliance has to be connected exclusively to a power supply with an earth safety device.

⚠ WARNING: Cleaning and sterilising new or repaired instruments. All the new or repaired appliance accessories are not sterile. At first use and after each treatment they must be cleaned and sterilised by carefully following the instructions in Chapter 8 on page 31.

⚠ WARNING: Infection control.

For maximum patient and operator safety, before using all the parts and reusable accessories, make sure they have been previously cleaned and sterilised by following the instructions in *Chapter 8* on page 31.

❗ CAUTION: Contraindications. After having autoclaved the handpiece, inserts, torque wrench or any other sterilisable accessory, wait for them to completely cool down before reusing them.

⚠ WARNING: Insert Breakage and Wear.

High frequency oscillations and wear can, in rare cases, lead to the breakage of the insert. Deformed or damaged inserts are susceptible to breakage during use. Broken or worn inserts must not ever be used.

In the event of breakage check that no fragments remain in the treated part and at the same time use suction effectively to remove them.

It is necessary to instruct the patient to breathe through the nose during the treatment, or use a dental dam, so as to avoid to ingest fragments of broken inserts.

Check the state of wear of the insert and its integrity before and during each use.

If a drop in performance occurs, see to its replacement.

The state of wear of the most common inserts (S1, S1-S, S2, S5, P2, P4, P10) can be checked using the supplied INSERT-CARD. To use the INSERT-CARD correctly:

- Position the insert on the INSERT-CARD so that the profile corresponds to the one printed on the card. The printed profile on the card has a red line

indicating the limit of wear;

- If the insert is shorter than the limit of wear, its performance will be significantly inferior compared to that of a new insert, and should therefore be replaced.

If the layer of titanium nitride (gold-plated surface), where present, is visibly worn, the insert must be replaced. The use of a worn insert reduces its efficiency.

Diamond Inserts: diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.

When the nitriding is consumed, the cutting edge loses effectiveness; a possible regrounding harms the insert, therefore it is prohibited. Verify that the insert is not worn. During the operation frequently check that the insert is intact, especially in the apical part.

During the operation avoid prolonged contact with the retractor or with metal instruments in use. Do not exert excessive pressure on the inserts during use.

⚠ WARNING: Exclusively use original MECTRON inserts, accessories and spare parts.

❗ CAUTION: No modification is allowed to this appliance.

⚠ WARNING: In case of an adverse event and/or accident attributable to the device during correct use and in accordance with the intended use, a report must be made to the Competent Authority and manufacturer indicated on the product label.

1.5 Symbols

EN

Symbol	Description	Symbol	Description
	Device compliant with Regulation (EU) 2017/745. Notified body: IMQ S.p.A.		Nemko brand Electrical Safety ES 60601-1 Compliance with UL - CSA regulations
	Medical device		Caution! Read the instructions for use
	Operating instructions		Manufacturer
	Date of manufacture		Serial number
	Batch number		Product code
	Not sterile		The sterilisable materials must be autoclaved and can resist up to a maximum temperature of 135°C
	"B" Part applied according to regulation EN 60601-1		Alternating current
	Connection of the control pedal		Generic warning signal a)
I	Power switch set to "on"	0	Power switch set to "off"
	The appliance and its accessories must not be disposed of or treated as municipal solid waste		Temperature limits for transport and storage
	Moisture limits for transport and storage		Atmospheric pressure limits for transport and storage

Symbol	Description	Symbol	Description
QTY.1	Quantity in package: 1	Rx only	For US market only CAUTION: US federal law restricts sale to orders by licensed dentists or dental hygienists only.

Table 1 – Symbols

- a) The symbol is represented by a yellow triangle and black graphic symbol.

NOTE: For other symbols, refer to Chapter 12.1 on page 67.

2 IDENTIFICATION DATA

A correct description of the model and of the serial number of the appliance will allow the After Sales Service to provide fast and

effective answers.

Always provide this information every time that you contact MECTRON After Sales Service.



2.1 Appliance ID Plate

Each appliance is equipped with an identification plate indicating the main technical characteristics and serial number. The identification plate is located under the appliance. The complete technical specifications are provided in Chapter 11 on page 58.

NOTE: The complete list of symbols is provided in Chapter 1.5 on page 6.

Additional symbols and characteristics of the device are provided on a separate plate. This identification plate is located under the appliance.

NOTE: The complete list of symbols is provided in Chapter 1.5 on page 6.


2.2 Handpiece Identification Data

The handpiece name, MECTRON logo, serial number, product code, HIBC matrix code, and a series of symbols are laser engraved on the scaler handpiece.


NOTE: The complete list of symbols is provided in Chapter 1.5 on page 6.

EN

2.3 Inserts Identification Data

The name of the insert itself (Reference 1), MECTRON logo (Reference 2), and insert batch number (Reference 3) are laser engraved on each insert.

3 DELIVERY

3.1 List of Components

See the Figure on the inside cover.
multipiezo is supplied with the standard equipment (see Table 2 on page 9), a set of variable accessories depending on the configuration and client requests (see Table

3 on page 11), and accessories that can be ordered separately (see Table 4 on page 12).

NOTE: Both the items included in the standard supply and all accessories can be ordered separately by the client.

Standard supply			
Item	Code	Description	Ref.
Device core unit	04000100		A
Irrigation tank complete with cap	03150102		B
Pedal FS-05 with cord and plug	02900114		C
Tank safety cap	01950014	Grey safety cap	D
	01950015	Blue safety cap	
Tank cap	03020082	Grey tank cap	E
	03020081	Blue tank cap	
Scaler Handpiece	03120238	MANIPOLO ABLATORE SLIM (SLIM SCALER HANDPIECE)	F
multipiezo O-ring kit	03020175		H

Table 2 – Standard supply.

Accessories available for order with standard supply			
Item	Code	Description	Ref.
Torque wrench	02900137	K10 torque wrench ^{b)}	I
	02900074	K6 torque wrench ^{d) b)}	
	02900081	K7 torque wrench ^{e) b)}	
Inserts	0296xxxx	"S" series reusable scaler inserts ^{b)}	J
	0308xxxx	"PE" series reusable scaler inserts ^{b)}	
	0305xxxx	"R" series reusable scaler inserts ^{b)}	
	0345xxxx	"ER" series reusable scaler inserts ^{b)}	
	0235xxxx	"E" series reusable scaler inserts ^{b)}	
	0299xxxx	"D" series reusable scaler inserts ^{b)}	
	0219xxxx	"CM" series reusable scaler inserts ^{b)}	
	0313xxxx	"ME" series reusable scaler inserts ^{b)}	
	03570004	"ICS" series reusable base insert ^{b)}	
	03590009	"IC1" series tip ^{b)}	
	02900112	IC1 tips kit (5 pcs.) ^{b)}	
	03570001	DB1 reusable base support for inserts ^{b)}	
	03590001	TA12D60 diamond insert ^{b)}	
	03590004	TA12D90 diamond insert ^{b)}	
	03590002	TA14D60 diamond insert ^{b)}	
	03590005	TA14D90 diamond insert ^{b)}	

Accessories available for order with standard supply			
Item	Code	Description	Ref.
	03590007	TA14D120 diamond insert ^{b)}	
	03590003	TA16D60 diamond insert ^{b)}	
	03590006	TA16D90 diamond insert ^{b)}	
	03590008	TA16D120 diamond insert ^{b)}	
	03590010	TF12D60 diamond insert ^{b)}	
	03590011	TF16D60 diamond insert ^{b)}	
	03590012	TF12D90 diamond insert ^{b)}	
	03590013	TF16D90 diamond insert ^{b)}	
Use and Maintenance Manual	02150273	IT/DE/SV version ^{b)}	K
	02150274	EN/FR/ES version ^{b)}	
Information: Documentation	Online 02150650		
Power supply cable	00050030	Black power supply cable 3X0.75 2m IT ^{c)}	L
	00050031	Black power supply cable 3X0.75 2m D ^{c)}	
	00050034	Black power supply cable 3X0.75 2m UK ^{c)}	
	00050010	Black power supply cable 3X1 2.5m AU ^{c)}	
	00050032	Black power supply cable 3X0.75 2m CH ^{c)}	
	00050029	Black power supply cable 3X0.75 HG 2m JP ^{c)}	
	00050028	Black power supply cable 3X0.75 HG 2m USA ^{c)}	

Table 3 – Accessories available for order with standard supply.

Accessories available for separate order			
Item	Code	Description	Ref.
Scaler Handpiece	03120142	MANIPOLO ABLATORE LED (LED SCALER HANDPIECE) ^{b)}	M
Front cone for LED handpiece	03020158	Front cone for LED SCALER HANDPIECE ^{b)}	N
Light guide kit	02900146	Front cone, light guide and decorative ring for SLIM scaler handpiece ^{b)}	P
Front cone without light	03020171	Front cone for LED SCALER HANDPIECE ^{b)}	G
Peristaltic pump rotor assembly	03020226-001	c)	O

Table 4 – Accessories available for separate order.

b) Produced by Mectron.

c) Distributed by Mectron.

d) To be used only for ME type inserts.

e) To be used only for the "Crown Preparation" kit.

The packaging of the appliance cannot undergo strong impacts as contains electronic components, therefore the transport and the storage must be carried out with particular care.

Do not stack various boxes in order to avoid pinching the underlying packaging.

All the material sent by MECTRON was controlled at the time of dispatch.

The appliance is shipped appropriately protected and packaged.

Upon receipt of the appliance, check for any

possible damage caused during transport and in case any damage and/or defects is found, complain to the transporter.

Keep the packaging in case any items need to be sent to an Authorised MECTRON Service Centre and to store the appliance during long periods of non-use.

⚠ WARNING: Before starting the treatment, always make sure to have a stock of material (handpiece, inserts, wrenches) to use in the event of failures or drawbacks.

4 INSTALLATION

4.1 First Installation

The appliance must be installed in a suitable place that is convenient for its use.

⚠ WARNING: The place where the device is installed must meet the requirements set out in *Chapter 4.2 on page 13*.

multipiezo can be purchased ready for use

or may need to be enabled by typing in an activation key.

In case your appliance requires an activation key, the procedures to be followed may vary from country to country.

Always refer to your dealer for details.

4.2 Safety Requirements During Installation

⚠ WARNING: Contraindications.

Interference with other equipment. Even if compliant with standard IEC 60601-1-2, multipiezo may interfere with other devices in the vicinity. multipiezo must not be used near or stacked with other appliances. However, if this should become necessary, it is necessary to check and monitor the correct functioning of the appliance in that configuration.

⚠ WARNING: Contraindications.

Interference from other equipment. An electrosurgical scalpel or other electrosurgical units positioned near the multipiezo appliance may interfere with the correct operation of the appliance itself.

① CAUTION: The electrical system of the places where the appliance will be installed and used must comply with the laws in force and the relative electrical safety regulations.

① CAUTION: To avoid the risk of electrical shocks, this appliance has to be connected exclusively to a power supply with an earth safety device.

① CAUTION: Place the appliance in such a way as to always have the power plug easily accessible as it is considered to be a disconnection device.

⚠ WARNING: Risk of explosion.

The appliance cannot operate in environments where there are saturated atmospheres of flammable gases (aesthetic mixtures, oxygen, etc.).

⚠ WARNING: Install the appliance in a safe place protected from impact or accidental water or liquid spray.

⚠ WARNING: Do not install the appliance above or near sources of heat. Arrange during installation for a suitable circulation of air around the appliance.

① CAUTION: Do not expose the appliance to direct sunlight or sources of UV light.

① CAUTION: The appliance can be transported but must be handled with care when moved. Position the pedal on the floor so that it can only be activated intentionally by the operator.

① CAUTION: Before connecting the handpiece to its cord, check that the electrical contacts are perfectly dry, on both parts. If necessary dry them with compressed air.

① CAUTION: Each tank can hold liquids up to a maximum of 500 ml.


4.3 Connecting the Accessories

EN

Connect the pedal to the back of the appliance, on the right-hand side, in the socket marked with the symbol by inserting the pedal cord plug until you hear a 'click'.

① CAUTION: Pay attention to the positioning of the pedal, which must be such that the pedal is only activated intentionally by the operator.


1

Insert the supply cable in its connection situated on the back of the appliance. Connect it to the socket in the wall;

① CAUTION: Position the device in a way that the power plug is easily accessible at all times; this is considered a means of isolation.


2

Unscrew the cap of the tank and fill it with the desired liquid;

① CAUTION: Each tank can hold liquids up to a maximum of 500 ml.


3

Check that the small tube inside the cap is properly installed, then screw the cap onto the tank;


! CAUTION: Check that the female coupling of the tank cap is clean and free from obstructions.

4

! CAUTION: Check that the male coupling on the appliance core unit is clean and that its O-rings are not worn.


! CAUTION: Use the support only to install the Mectron 500 ml tank and to stow the handpiece. Do not use the support for other purposes.

Holding the tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;


! CAUTION: Do not overturn the tank, the cap is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

5

Correctly insert the scaler handpiece onto its cord by matching the alignment notch on the handpiece connector with the alignment key on the cord connector. Check that the electrical contacts of both are perfectly dry and if necessary dry them by blowing with compressed air;


6

Place the handpiece on its respective support.

! CAUTION: Use the support only to install the Mectron 500 ml tank and to stow the handpiece. Do not use the support for other purposes.

7

5 USE

5.1 Power On/Off

Switching the appliance on

Switching the appliance on

The switch is located on the left-hand side of the machine core unit.

Flick the switch into the "I" position, taking care not to press the pedal.

On start-up, 4 symbols appear (Ref. 7 inside cover) which will progressively turn off.


The system will then show the preset configuration and be ready for use.

Switching the appliance off

The switch is located on the left-hand side of the machine core unit.

Flick the switch into the "O" position, taking care not to press the pedal.

The system will turn off.

! CAUTION: Position the device in a way that the power plug is easily accessible at all times; this is considered a means of isolation.

NOTE: each time the appliance is started up, the default setting is applied:


- "power": 1;
- "light": AUTO;
- "function": ENDO;
- "irrigation": 3.

NOTE: When the appliance is turned on, the bottle lights up.

5.2 Description of the Keyboard

TOUCH KEYBOARD

The user can configure the machine by simply touching the touch keyboard. Depending on the configured setting, the electronic feedback system will automatically adjust the correct working frequency.

POWER

(Ref. 1 inside cover)

The operating power can be regulated by selecting the numbers on the touch keyboard in the "power" column. There are 7 power levels:

- 0: the handpiece has no power: the insert emits only irrigation;
- From 1 to 6: the operating power can be regulated incrementally.

The possibility to select the power levels is available for all functions, as follows:

- **ENDO:** 7 possible levels: from 0 to 6;
- **PERIO/SCALER:** 7 possible levels: from 0 to 6;
- **RESTORATIVE:** 7 possible levels: from 0 to 6.

FILLING AND CLEANING THE IRRIGATION CIRCUIT

(Ref. 2 inside cover)

The appliance includes the "pump/flush" key, which depending on the operating mode allows the PUMP function or FLUSH function to be run.

The PUMP function can be used at the start of the treatment to allow the liquid to reach the insert, so as to start the surgery with the necessary irrigation (See Chapter 5.4 on page 23).

The FLUSH function (see Chapter 6 on page 27) is used to run an irrigation circuit cleaning cycle.

This function must be used after each time the appliance is used and before cleaning, disinfecting and sterilising all its parts (See Chapter 8 on page 31).

LIGHT

(Ref. 3 inside cover)



Depending on the type of treatment that needs to be carried out, 3 possible options can be selected from the "light" list:

- If the **AUTO** option is selected, the LED light on the handpiece nose cone lights up when the pedal is pressed and automatically switches off 3 seconds after the pedal is released;
- If the **ON** option is selected, the LED light on the handpiece nose cone stays lit, regardless of the pedal. The light turns off 100 seconds after the last time the pedal is pressed, and the option switches from ON to AUTO;
- If the **OFF** option is selected, the LED light on the handpiece nose cone stays off.

NOTE: Adopting the "MANIPOLO ABLATORE LED" ("LED SCALER HANDPIECE" - cod. 03120142) equipped with the "Front cone without light" (cod. 03020171), this function will not be available.

The position of the LED light on the nose cone of the "MANIPOLO ABLATORE LED" ("LED SCALER HANDPIECE" - cod. 03120142) can be adjusted as follows:

1. Hold the handpiece body and slightly unscrew the metal ring at the base of the nose cone by turning it counter-clockwise;
2. Rotate the nose cone so that the LED light is in the desired and necessary position;
3. To fix the position, screw the metal ring in a clockwise direction.

FUNCTIONS

(Ref. 4 inside cover)

Depending on the type of operation, it is possible to select one of 3 options from the "function" list, as follows:

- **ENDO:** dedicated to endodontic treatments such as root canal cleaning and the retrograde approach.
- **PERIO/SCALER:** dedicated to all supra- and subgingival dental prophylaxis, root planing, and implant surface cleaning procedures.
- **RESTORATIVE:** dedicated to preservative techniques, implants and extractions.

NOTE: when RESTORATIVE is selected with power level 6, the PULSE function is activated, which is used to optimise the performance of inserts used in extraction and prosthetic techniques.


SOFT MODE

(Ref. 5 inside cover)

The treatment can be made more delicate for sensitive patients, by activating the "soft mode" function.

"Soft mode" can only be activated in the "perio/scaler" functions with powers from 1 to 6 and "restorative" with powers from 1 to 5.

In "endo" mode the function is not available.

IRRIGATION

(Ref. 6 inside cover)


The peristaltic pump flow rate can be adjusted by selecting the numbers on the touch keyboard, in the "irrigation" column. There are 7 flow rate levels:

- **0:** the pump function is off, no irrigation flows out from the insert.
- **From 1 to 6:** the pump flow rate ranges from 6 ml/min to 28 ml/min.

The possibility to choose the irrigation flow rate is linked to the selected function, as follows:

- **ENDO:** 6 flow rate levels: from 1 to 6. Irrigation can be set to 0;
- **PERIO/SCALER:** 6 flow rate levels: from 1 to 6. Irrigation can be set to 0;
- **RESTORATIVE:** 6 flow rate levels: from 1 to 6. Irrigation can be set to 0 only with powers 1 to 5. For RESTORATIVE power level 6, irrigation is necessary.

SYMBOLS

(Ref. 7 inside cover)

multipiezo is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch keyboard by means of a symbol.

To help the user identifying the part that is non-functioning there are four symbols described in Chapter 12.1 on page 67.

5.3 Safety Requirements Before and During Use

⚠ WARNING: Before starting work always ensure that you have all replacement parts (handpiece, inserts, keys) to use should there be downtime due to a fault or problem.

⚠ WARNING: Exclusively use original Mectron inserts, accessories and replacement parts.

⚠ WARNING: Use of non-original inserts not produced by MECTRON: this entails permanent damage to the thread of the handpiece and compromises correct operation with the risk of causing harm to the patient.

① CAUTION: Contraindications -

Ultrasonic Scaler. Do not perform treatments on metal or ceramic prosthetics. Ultrasonic vibrations could lead to decementation of the dentures.

⚠ WARNING: Contraindications. Do not use multipiezo on patients with Pacemakers or other implantable electronic devices. This regulation also stands for the operator.

⚠ WARNING: Do not perform scaling treatments without water spray in order to avoid an overheating of the insert that can cause damage to the tooth. Treatments without water spray can only be those performed with the "Dry Work" inserts without water passage.

① CAUTION: In treatments requiring irrigation, use only inserts with a liquid passage.

⚠ WARNING: Treatments that require irrigation. Always check the operation of the irrigation before and during use. Make sure that the liquid exits from the insert. Do not use the appliance if the irrigation does not work or if the pump is faulty.

⚠ WARNING: Checking the condition of the device before treatment. Always check that there is no water under the appliance. Before every treatment always check the perfect operation of the appliance and efficiency of the accessories. If operation anomalies are observed do not carry out the treatment. Contact an Authorised Mectron Service Centre if the faults concern the device.

⚠ WARNING: Infection control. First use: All parts and reusable accessories (new or returned from an Authorised MECTRON Service Centre) are delivered in NON-STERILE conditions and must be treated before each use, according to the instructions provided in Chapter 8 on page 31. **Subsequent uses:** After each treatment, clean and sterilise all parts and reusable accessories according to the instructions in Chapter 8 on page 31.

⚠ WARNING: Infection control. Do not leave liquids in the tank for long periods of time. The tank must be filled after having cleaned and sterilised all parts, and before a treatment. If the tank has been filled without having used the device, it must be emptied at the end of the day, and all parts and accessories must be cleaned and sterilised.

⚠ WARNING: To ensure that the handpiece cools down activate it always with irrigation circuit installed and properly filled. To fill the irrigation circuit always use the FLUSH function.

① CAUTION: For correct use of the appliance, the pedal must be pressed and the appliance started with the insert not in contact with the part to be treated, so that the electronic circuit is able to recognise the best point of resonance of the insert without interference, thus allowing optimum performance.

⚠ WARNING: Before each treatment ensure that the insert that is suitable for the treatment is inserted on the handpiece. Use only the Mectron torque wrench to secure the insert to the handpiece.

⚠ WARNING: The patient must not come into contact with the machine core unit or with the pedal.

⚠ WARNING: Do not change the insert whilst the handpiece is in operation in order to avoid causing injury to the operator.

① CAUTION: Pump/flush function. After use with aggressive and non-aggressive solutions, the tubes and handpiece must be cleaned using the FLUSH function (See Chapter 6 on page 27). If the tubes are not cleaned, the crystallisation of the salts may seriously damage the device.

① CAUTION: Pump/flush function. The "FLUSH" function must be used after each treatment, before starting the cleaning and sterilisation procedures.

⚠ WARNING: Do not use the device in case of breaks or cracks on the casing.

⚠ WARNING: Insert Breakage and Wear. High frequency oscillations and wear can, in rare cases, lead to the breakage of the insert. Deformed or damaged inserts are susceptible to breakage during use. Broken or worn inserts must not ever be used.

In the event of breakage check that no fragments remain in the treated part and at the same time use suction effectively to remove them.

It is necessary to instruct the patient to breathe through the nose during the treatment, or use a dental dam, so as to avoid to ingest fragments of broken inserts.

Check the state of wear of the insert and its integrity before and during each use. If a drop in performance occurs, see to its replacement.

The state of wear of the most common inserts (S1, S1-S, S2, S5, P2, P4, P10) can be checked using the supplied INSERT-CARD. To use the INSERT-CARD correctly:

- Position the insert on the INSERT-CARD

- so that the profile corresponds to the one printed on the card. The printed profile on the card has a red line indicating the limit of wear;
- If the insert is shorter than the limit of wear, its performance will be significantly inferior compared to that of a new insert, and should therefore be replaced.

If the layer of titanium nitride (gold-plated surface), where present, is visibly worn, the insert must be replaced. The use of a worn insert reduces its efficiency.

Diamond Inserts: diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.

When the nitriding is consumed, the cutting edge loses effectiveness; a possible regrounding harms the insert, therefore it is prohibited. Verify that the insert is not worn. During the operation frequently check that the insert is intact, especially in the apical part.

During the operation avoid prolonged contact with the retractor or with metal instruments in use. Do not exert excessive pressure on the inserts during use.

① CAUTION: Contraindications. After autoclaving the handpiece, the inserts, the torque wrench, or any other sterilisable accessory, wait until it has completely cooled before reuse.

① CAUTION: The electrical contacts inside the connectors on the handpiece and cord must be dry. Before connecting the handpiece to its cord, make sure that the electrical contacts of the connector on both sides are perfectly dry, especially after the sterilisation cycle in an autoclave. If necessary dry the contacts by blowing them with compressed air.

① CAUTION: Owing to its conformation, the handpiece can rotate. When not used the handpiece must always be put back on its support stand.

⚠ WARNING: During the intervention on the patient, do not perform any maintenance tasks on the system.

5.4 Instructions for Use

After having connected all the accessories as described in Chapter 4.3 on page 14 proceed as follows:

To fill the irrigation circuit, use the PUMP function by selecting PUMP/FLUSH on the touch keyboard: all the other selections will be deactivated and the writing PUMP/FLUSH will flash;


1


While the writing PUMP/FLUSH is flashing, press the pedal once and then release it, PUMP/FLUSH will stop flashing and the irrigation circuit will start to fill;

EN

2


3


4

The cycle can be stopped as soon as liquid is seen flowing out from the multipiezo, handpiece by pressing PUMP/FLUSH or alternatively by pressing on the pedal. The PUMP function is deactivated, and the keyboard becomes active again in the last setting used;

5

Screw in the preselected insert on the multipiezo handpiece until it comes into contact;


Tighten the insert using the Mectron torque wrench.

For correct use of the Mectron torque wrench, proceed as follows:

- Place the insert inside the torque wrench, as shown;

6


Hold the handpiece body firmly;

7

! CAUTION: Do not grab the handpiece on the end and/or on the cord, but only on the body. Do not turn the handpiece, rather it must be held steady turning only the wrench.

Turn the torque wrench in a clockwise direction until the notch clicks (the external body of the torque wrench turns with respect to the body of the handpiece, emitting mechanical "CLICK" sounds).


The insert is now optimally locked;

Select the type of function, power and irrigation necessary and the light if desired, on the keyboard;


Lift the handpiece and press the pedal to start treatment;

8

At the end of treatment stow the scaler handpiece on its housing.

9

5.5 Important information on the Inserts

EN

⚠ WARNING:

- Before using the sterilised insert, check the integrity of the sterile packaging and inspect the product to ensure there is no damage. The insert will no longer be sterile if the packaging is broken or damaged. If the packaging is damaged, the insert MUST be sterilised again before use.
- Before starting to work, properly tighten the insert to the handpiece using the torque wrench.
- When the layer of titanium nitride is visibly consumed the insert must be replaced. The use of an insert too worn decreases its efficiency.
- **Diamond Inserts:** diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.
- Do not activate the handpiece while the insert is in contact with the part to be treated so that the electronic circuit can recognise the best point of resonance of the insert allowing optimum performance.
- Check the state of wear of the insert and its integrity before and during each use. If you experience a loss of performance provide for its replacement.
- Use only original Mectron inserts. The use of non-original inserts will void the warranty and damage the thread of the multipiezo handpiece, with the risk of no longer being able to correctly screw in the original inserts for subsequent use. In addition, the machine settings are tested and guaranteed for correct operation using only original Mectron inserts.
- Do not change in any way the shape of the insert, by bending it or filing it. This could cause it to break.
- Do not use an insert that has undergone deformation of any type.

- Do not attempt to sharpen an insert used.
- Always check that the threaded parts of the insert and the handpiece are perfectly clean - See Chapter 8 on page 31 - Cleaning and sterilisation.
- Excessive pressure applied to the insert can cause it to break and potentially injure the patient.
- The Mectron inserts vibrate with longitudinal oscillation, with a back and forward movement. During the treatment, always keep the instrument in a tangential direction with respect to the surface of the tooth. Move the handpiece back and forward, applying slight side pressure.
- Do not point the instrument directly onto the surface of the enamel or implant. Position the tip/operating part only in a tangential direction with respect to the surface of the tooth or implant.
- The insert must be held in constant movement. If the insert is blocked it may cause overheating of the treated part. It is recommended to use continuous motions to minimise contact between the tip and the part. Do not block them against the tissues so as not to cause their overheating. It is advisable to increase the level of irrigation as the power level increases.
- Leave the ultrasonic vibrations to work, do not exert excessive pressure on the inserts during use. Apply slight force on the insert for optimum performance.
- When the insert is used in interproximal spaces, do not block the instrument or leverage on the operating part. The inserts must be left free to vibrate.
- In endodontic root canal therapy, never operate the files outside of the root canal to prevent them from breaking. To prevent breakages, create a smooth path with a manual endo file and plan

the most direct possible path of access to minimise any bending of the insert. Make delicate movements. Examine the file often for signs of wear. If the file breaks inside the canal, avoid contact

between the instrument and the broken file to stop it from being pushed further in. Do not apply pressure on the insert in an axial direction.

6 PUMP/FLUSH FUNCTION

The pump/flush function allows the irrigation circuit rinse cycle to be performed.


! CAUTION: Pump/flush function.

The "pump/flush" function must be used after each treatment, before starting the cleaning and sterilisation procedures.

! CAUTION: If the tubes are not rinsed, the crystallisation of the salts may seriously damage the appliance.


Fill the tank with water, preferably demineralised water. Connect the tank to the appliance;

1

Place the handpiece, with or without insert, above a container for the liquid that will be released during the cleaning cycle;

2

To enter cleaning mode, select "pump/flush" on the touch keyboard: all the other selections will be deactivated;

3

NOTE: The "pump/flush" mode can be deactivated at any time by pressing the "pump/flush" key again. The keyboard will be reactivated and configured in the last setting used.


EN

To start the "pump/flush" function while the writing "pump/flush" is flashing, press the pedal once and release it: the keyboard will stop flashing and the cleaning cycle will start;


NOTE: The function can be stopped at any time by pressing "pump/flush" again or alternatively by pressing the pedal briefly. The keyboard will be reactivated and configured in the last setting used.

4

As soon as the peristaltic pump starts working, the entire scale of values in the "irrigation" section lights up, and during the passage of the liquid, the irrigation value switches from 6 to 0. The cycle lasts 20 seconds. Once ended, the keyboard will be reactivated and configured in the last setting used.

5


After cleaning the irrigation circuit, disassemble the individual parts (see Chapter 7 on page 28) and proceed to clean and sterilise them (see Chapter 8 on page 31).

7 DISASSEMBLING PARTS FOR CLEANING AND STERILISATION

Before carrying out the cleaning procedures described in Chapter 8 on page 31, disconnect all accessories and components of the multipiezo.

⚠️ WARNING: Switch the appliance off. Always switch the appliance off using the switch and disconnect the power supply cable from the wall socket and machine core unit before carrying out cleaning and sterilisation tasks.

1


DISASSEMBLING PARTS FOR CLEANING AND STERILISATION

EN


If present, remove the insert from the handpiece using the torque wrench;

2


Disconnect the handpiece from the relative cord;

3

① CAUTION: Do not try to unscrew or turn the connector when disconnecting the handpiece.
The connector can be damaged.


Unscrew the nose cone from the handpiece;

4


① CAUTION: The nose cone has a light pipe. If the nose cone is unscrewed, the light pipe will no longer be held in place and may slide and be disconnected. Be careful not to lose the light pipe.

NOTE: On the "Front cone for LED handpiece" (cod. 03020158) of the MANIPOLO ABLATORE LED (LED SCALER HANDPIECE - cod. 03120142), the metal ring cannot be separated from the plastic cone.


Remove the light pipe from the nose cone.

5

Disconnect the tank from the machine core unit, pulling it outwards;

5


Disconnect the pedal from the appliance: grasp the pedal connector, press the release tab and pull the connector back;

6

! CAUTION: Do not try to unscrew or turn the connector during disconnection: this might damage the connector.

! CAUTION: While disconnecting the cable from the pedal, keep always and only the connector of the cord. Never pull the cord itself.

8 CLEANING AND STERILISATION

This table is purely indicative. For the complete cleaning and sterilisation procedures of the individual parts, refer to the paragraphs indicated in the table.

 **CAUTION:** Methods not considered in the table below must not be used.

EN

Machine Core Unit, Pedal, End Part of Cord							
Phase	Chapter	Procedure					
I	8.1	Preparation					
II	8.2	Manual cleaning with enzymatic detergent					
Tank and Cap							
Phase	Chapter	Procedure					
III	8.3	Manual cleaning with enzymatic detergent and running water					
Accessories (Scaler Handpiece, Scaler Nose Cone, Inserts, Torque Wrench)							
Phase	Chapter	Procedure	Scaler Handpiece	Scaler Nose Cone	Inserts	Torque wrench	
IV	8.4.1	Manual cleaning	Manual cleaning with Enzymatic Detergent	Immersion into enzymatic detergent		Ultrasonic Cleaning	
	8.4.2	Automatic cleaning		Thermal-disinfector with detergent			
V	8.5	Check cleaning	Visual inspection with Lens 2.5x				
VI	8.6	Drying and lubrication	Drying		Drying and Lubrication		
VII	8.7	Sterilisation	Packaging, Sterilisation and Storage				

NOTE: Repeated reconditioning has a minimal effect on the devices and their accessories. The end of the service life of the accessories is generally determined by wear or damage resulting from use. Mectron guarantees the integrity of its sterilisable scaler handpieces for up to 250 reconditioning cycles.

8.1 Preparation

1. Run the flush function (see Chapter 6 on page 27);
2. Disinfect the tank irrigation circuit as described in Chapter 8.1.1 on page 32.
3. Check that all of the following accessories have been removed/disconnected from the machine core unit (see Chapter 7 on page 28):
 - Power supply cable;
 - Pedal;
 - Scaler handpiece;
 - Inserts.

⚠ WARNING: Always switch off the appliance using the O/I switch and disconnect it from the power mains before carrying out cleaning tasks.

❗ CAUTION: Always disconnect the insert from the handpiece before cleaning and sterilising it.

❗ CAUTION: Do not immerse the handpiece in disinfecting solutions or other liquids as these could damage it.

❗ CAUTION: Do not immerse the handpiece in the ultrasonic tank.

❗ CAUTION: The cleaning and sterilising operations described in the following sections are to be performed at first use and at all subsequent uses.

8.1.1 Disinfecting the Tank Irrigation Circuit

After performing the pump/flush function (see Chapter 6 on page 27) and before carrying out the subsequent cleaning procedures, proceed to disinfect the tank irrigation circuit. To disinfect the irrigation circuit, proceed as follows:


8.1.1.1 Required Materials

- CIDEX® OPA disinfectant solution;
- Distilled water.

8.1.1.2 Procedure


Disconnect the irrigation solution tank from the machine core unit, pulling it outwards;

1


Unscrew the cap from the irrigation solution tank and empty it;
Rinse the irrigation solution tank with running water;

2

Fill the irrigation solution tank with 250ml of disinfectant (CIDEX® OPA);


3

Holding the irrigation solution tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;


4

! CAUTION: Do not overturn the tank, the cap is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

Lift the scaler handpiece connector;


5

Position the scaler handpiece connector above a container or sink to catch the liquid released during the disinfection cycle;

! CAUTION: Before connecting and using the scaler handpiece, make sure that the electrical contacts on both the scaler handpiece and its cord are perfectly dry.

6

Select "pump/flush" on the touch keyboard (see Chapter 6 on page 27). All the other selectable options on the keyboard will be disabled.

NOTE: The "pump/flush" mode can be deactivated at any time by pressing the "pump/flush" key again. The keyboard will be reactivated and configured in the last setting used.

7


Press the pedal once and release it to start the cycle.

The keyboard will stop flashing and the cleaning cycle will start;


NOTA: the function can be stopped at any time by pressing "pump/flush" again or alternatively by pressing the pedal briefly. The keyboard will be reactivated and configured in the last setting used.

8

As soon as the peristaltic pump starts working, the entire scale of values in the "irrigation" section lights up, and during the passage of the liquid, the irrigation value switches from 6 to 0. The cycle lasts 20 seconds. Once ended, the keyboard will be reactivated and configured in the last setting used;

9

Repeat the "flush" cycle (see points 7 to 9 of the procedure) a second time;

10

At the end of the "flush" cycles, wait 12 minutes;

11


12 min.

Repeat the "flush" cycle (see points 7 to 9 of the procedure) 3 times;

12

At the end of the "flush" cycles, wait 12 minutes;

13


12 min.

Repeat the "flush" cycle (see points 7 to 9 of the procedure) 3 times;

14

At the end of the "flush" cycles, wait 12 minutes;

15


12 min.

Disconnect the irrigation solution tank from the device;

Unscrew the cap from the irrigation solution tank and drain the remaining liquid;


Rinse the irrigation solution tank with running water;

16

With the tank disconnected, run a flush cycle (see points 7 to 10) to empty the irrigation circuit of the liquid inside;


17

Fill the irrigation solution tank with 300ml of distilled water;

NOTE: Use of distilled water is recommended for long periods of non-use. If the device is used after the disinfection procedure, use common running water.


18

Holding the irrigation solution tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

! CAUTION: Do not overturn the tank, the cap is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

19

Run a number of "flush" cycles (see points 7 to 9) sufficient to completely drain the water from the tank;


20

After cleaning the irrigation circuit, disassemble the individual parts (see Chapter 7 on page 28) and proceed to clean and sterilise them (see Chapter 8 on page 31).

21

! CAUTION: If the handpiece is stowed back in its support, with the "pump/flush" function selected (flashing LED) but not running, the function is deactivated;

! CAUTION: During the "pump/flush" cycle, i.e. when the function has been selected and activated by pressing the pedal, if the handpiece is stowed back in the support, the cycle is not interrupted.

8.2 Cleaning the Non-Sterilisable Parts

The following procedure must be performed on non-sterilisable parts of the appliance, except for the tank and its cap. The parts in question are:

- Machine core unit;

- Pedal and relative connection cable to the machine;
- End part of the scaler handpiece and relative cord.

EN

8.2.1 Required Materials


- Clean, soft, low-lint cloths;
- Cleaning solution (pH 6-9).

8.2.2 Cleaning Method

Clean the surface of the parts in question using a clean, dry, low-lint cloth, dampened with a detergent solution (pH 6-9);

Dry the parts using a dry, non-abrasive, low-lint cloth.

1

2

! **CAUTION:** Do not sterilise the parts in question. They may stop working and cause damage to people and/or property.

⚠ WARNING: Always switch off the appliance using the O/I switch and disconnect it from the power mains before carrying out cleaning tasks.

⚠ WARNING: The device and its parts that cannot be sterilised are not protected against the penetration of liquids. Do not spray liquids directly on the surface of the device or its parts that cannot be sterilised.

! **CAUTION:** Do not use running water to clean the parts in question.

! **CAUTION:** Do not soak these parts in liquid and/or various kinds of solutions.

! **CAUTION:** If disinfection is required, water-based disinfectant solutions must be used with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated, because they can fade the colour and/or damage the plastic materials. This is also valid for chemical products such as acetone and alcohol.

8.3 Cleaning the Tank and Cap

The following procedure must be carried out on the tank and cap of the device.

8.3.1 Preparation

1. Disconnect the tank from the machine core unit (see *Chapter 7 on page 28*);
2. Unscrew the cap from the tank.

! **CAUTION:** Do not sterilise the tank and cap in an autoclave. They could be damaged.

8.3.2 Required Materials

- Water;
 - Cleaning solution (pH 6-9);
 - Clean, soft, low-lint cloth;
 - Demineralised water.
- Very acid products (PH < 4);
 - Products containing aldehyde, amine and/or phenols
 - Acetone;
 - Methylethylketone;

! **CAUTION: DO NOT USE** as disinfectants::

- Very alkaline products (pH > 9);
- Products containing sodium hypochlorite;
- Products containing hydrogen peroxide;
- Products containing abrasive substances;

since they may discolour and/or damage the plastic materials.


The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.

8.3.3 Cleaning Method

Thoroughly rinse the inside and outside of both the tank and cap under running water;

EN


1


2

3

4


Thoroughly rinse the inside and outside of the tank and cap under running water to eliminate all residue of the detergent solution;

Remove any residue on the external or internal surfaces of the tank and cap using a soft, low-lint cloth, dampened with demineralised water;

When cleaning operations are over, check the pieces under an appropriate light source, paying attention to any soil residue and, if necessary, repeat the cleaning cycle.

5

8.4 Cleaning the Sterilisable Accessories

NOTE: The cleaning procedures must be performed immediately after each use. Immerse the insert and/or instrument in demineralised water or in an enzymatic detergent solution immediately after use. Do not leave residue or blood deposits on the inserts and instruments, eliminate larger impurities with a disposable cloth or paper towel.

The sterilisable parts of the appliance are:

- Scaler handpiece;
- Scaler nose cone;
- Inserts;
- Inserts tightening wrench.

Before proceeding with operations to check cleanliness (Chapter 8.5 on page 52), drying and lubrication (Chapter 8.6 on page 53) and then sterilisation (Chapter 8.7 on page 54), one of the three possible cleaning methods explained in detail in the following sub-sections must be selected.

! CAUTION: The instructions supplied below have been validated by the manufacturer of the medical device as ABLE to prepare a medical device for re-use. The process manager is responsible for ensuring that the processes repeated are effectively performed using the equipment, materials and staff in the reprocessing structure in order to obtain the desired result. This generally requires the validation and systematic monitoring of the process. Similarly, all deviations from the instructions provided by the processes manager must be adequately assessed to judge their efficiency and potential undesired consequences.

8.4.1 Manual Cleaning

Manual cleaning can be carried out as an alternative to the automatic cleaning described in Chapter 8.4.2 on page 50.

8.4.1.1 Required Materials

- Enzymatic detergent at pH 6-9;
 - Water;
 - Container for immersion in the enzymatic liquid;
 - Ultrasonic tank;
 - Clean, soft, low-lint cloths;
 - Brush with soft nylon bristles;
 - Syringe;
 - Demineralised water
- hypochlorite;
- Products containing hydrogen peroxide;
 - Products containing abrasive substances;
 - Very acid products (PH < 4);
 - Products containing aldehyde, amine and/or phenols
 - Acetone;
 - Methylethylketone;

since they may discolour and/or damage the plastic materials.


The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.

Prepare a solution of enzymatic detergent^{a)} with pH 6-9, according to the instructions of the manufacturer;

! CAUTION: Once used, dispose of the enzymatic detergent correctly, do not recycle it.

a) Process validated by independent bodies with Mectron ENZYMEC enzymatic detergent, 0.8% v/v.


1

8.4.1.2 Scaler Handpiece

Completely immerse the scaler handpiece, nose cone, and light pipe in the enzymatic solution. Leave to soak for 10 minutes at $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$;

2




3

Gently brush the entire surface of the scaler handpiece, nose cone, and light pipe for at least 20 seconds using a brush with soft nylon bristles, paying particular attention to the following areas:

- thread of the scaler handpiece;
- titanium shaft;

- external and internal parts of nose cone;


- external and internal parts of light pipe.

EN


Rinse the inner channel of the scaler handpiece using a 20 ml syringe previously filled with a new enzymatic solution. Repeat twice;

4


Remove the handpiece, nose cone, and light pipe from the enzymatic solution and gently brush their surfaces using a brush with soft nylon bristles under hot ($40^{\circ}\text{C} \pm 5^{\circ}\text{C}$) running water, for at least 10 minutes;

5


Fill the ultrasonic tank with the enzymatic detergent solution prepared according to the manufacturer's instructions. Use water at ambient temperature ($22^{\circ}\text{C} \pm 2^{\circ}\text{C}$);

6


Place the scaler handpiece, nose cone, and light pipe in the ultrasonic tank immersed in the enzymatic detergent solution at $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and run a cycle for at least 10 minutes;

7


Remove the handpiece, nose cone, and light pipe from the enzymatic solution and gently brush their surfaces using a brush with soft nylon bristles under hot ($40\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$) running water, for at least 10 minutes;

8


Place the scaler handpiece, nose cone, and light pipe in the ultrasonic tank immersed in demineralised water at $40\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ and run a cycle for at least 10 minutes;

9


Rinse the inner channel of the scaler handpiece using a 20 ml syringe previously filled with demineralised water. Repeat five times;

10

Dry the surface of the scaler handpiece, nose cone, and light pipe with a clean, non-abrasive, low-lint cloth.

11


8.4.1.3 Inserts

EN

Place the insert in a clean container, in a horizontal position. Add enough enzymatic detergent solution to completely cover it.

Leave the insert to soak in the enzymatic detergent solution for 10 minutes at ambient temperature ($22^{\circ}\text{C} \pm 2^{\circ}\text{C}$);

12


13

While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed.

Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.

NOTE: Thoroughly brush all the following parts of the insert for about 20 seconds:


- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- Internal cavities, grooves and fissures.

14


Remove the insert from the enzymatic detergent solution.

Thoroughly rinse and brush all the surfaces of the insert (see previous point) under running water for at least 10 minutes;


Use a disposable syringe to withdraw and inject the enzymatic cleaning solution into the hard-to-reach areas (through-holes/cannulae). Repeat this operation three times to ensure the effective removal of the soil from the internal surfaces of the through-hole.;

15


Rinse the internal channel of the insert with demineralised water at a pressure of 3.8 bar for at least 10 seconds so as to eliminate any soil residue.;

16


Place the insert in the ultrasonic tank immersed in the enzymatic detergent solution at $24\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ and run a cycle for at least 20 minutes;

17

Remove the insert(s) from the ultrasonic tank and rinse them in demineralised water;


18

EN


Brush the internal and external surfaces of the insert with a clean brush with soft nylon bristles, under running water;

19


Rinse the insert(s) in demineralised water;

20

Rinse the internal channel of the insert with demineralised water at a pressure of 3.8 bar for at least 10 seconds so as to eliminate any soil residue.

21


8.4.1.4 Inserts Tightening Wrench

Place the torque wrench in a clean container, in a horizontal position. Add enough enzymatic detergent solution to completely cover the torque wrench.

Leave the torque wrench to soak in the enzymatic detergent solution for 10 minutes at ambient temperature ($22^{\circ}\text{C} \pm 2^{\circ}\text{C}$). This procedure reduces the amount of blood, protein, and mucus on the instrument(s);

22


While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed.

Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.

NOTE: Thoroughly brush all the following parts of the torque wrench for about 20 seconds:


- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- External metal ring;
- Internal cavities, grooves and fissures.

Remove the torque wrench from the enzymatic detergent solution.


Thoroughly rinse and brush all the surfaces of the torque wrench (see previous point) under running water for at least 10 minutes;

24

Place the torque wrench in the ultrasonic tank immersed in the enzymatic detergent solution at $24^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and run a cycle for at least 20 minutes;


25

Remove the torque wrench from the ultrasonic tank and rinse it under running water so as to eliminate all detergent residue.

Brush the internal and external surfaces of the wrench with a clean brush with soft nylon bristles, under running water.

26

8.4.2 Automatic Cleaning

Automatic cleaning can be carried out as an alternative to the manual cleaning described in Chapter 8.4.1 on page 42.

8.4.2.1 Required Materials

- Alkaline detergent: neodisher® FA (0.2 % v/v);
- Neutralising liquid: neodisher® Z (0.1 % v/v);
- Water;
- Metal basket;
- Adapters;
- Thermal-disinfecter.

NOTE: Make sure that the accessories are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

Position the accessories in a metal basket. Connect the relevant adapter (supplied as an optional) to the scaler handpiece connector and then to the water jet cleaning connections of the thermal-disinfecter.


Repeat the same operation for the inserts by connecting them to the appropriate adapters supplied as optionals.

⚠ WARNING: Avoid overloading the thermal-disinfecter as this could compromise the effectiveness of cleaning.

⚠ WARNING: On completion of the cleaning cycle in the thermal-disinfecter, the scaler handpiece remains at the wash temperature for a long time. Use appropriate precautions when extracting the scaler handpiece from the thermal-disinfecter to prevent injury to the operator.

① CAUTION: Owing to its conformation, the scaler handpiece can rotate. When not in use, the scaler handpiece must always be placed on its support.

1

2

Sequence and parameters applicable to the cycle:

- 1 min, Rinse with cold water;
- 5 min, Wash with alkaline detergent at $55^{\circ}\text{C} \pm 2^{\circ}\text{C}$;
- 1 min, Neutralisation with suitable solution (1/3 cold water, 2/3 hot water);
- 1 min, Rinse with water (1/3 cold water, 2/3 hot water);
- 5 min, Thermal disinfection at 93°C with demineralised water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1, Table B.1 [4] thermal disinfection at a temperature of 90°C for 5 min determines a value of A0 3000.

8.5 Cleaning Check

8.5.1 Required Materials


- Light source;
- 2.5X Magnifier.

EN



Once the cleaning operations have been completed, check the scaler handpiece and scaler nose cone under an adequate source of light, if necessary using a magnifying glass 2.5X, paying attention to any parts that may conceal soil residue (threading, cavities, grooves) and, if need be, repeat the cleaning cycle if soil is still visible. Finally, check the integrity of those parts and those elements that could have deteriorated during use;

Repeat the control operations for the other accessories (inserts, tightening wrenches), repeating if necessary the cleaning cycle.

1

2


8.6 Drying and Lubrication

8.6.1 Required Materials


- Compressed air;
- Soft, low-lint cloth;
- Medical grade lubricant.

Thoroughly dry all parts of the scaler handpiece, scaler nose cone, and light pipe by blowing them with compressed air;

1


① CAUTION: The scaler handpiece electrical contacts must be dry at the start and end of the sterilisation cycle, before the cord is connected to the appliance. Always make sure that the electrical contacts of the connector are perfectly dry, if necessary dry them blowing with compressed air.


2

① CAUTION: Before starting the sterilisation cycle, make sure that the insert is thoroughly dry both internally and externally. To do this, blow compressed air both externally and through the internal passage hole. This will prevent the appearance of stains, streaks on the surface or oxidation inside the insert.

Dry the inserts tightening wrench using a soft, low-lint cloth;


3

Lubricate the inserts tightening wrench with medical-grade lubricants in the indicated point.

4

① CAUTION: Do not use oil or silicone-based lubricants.


8.7 Sterilisation

8.7.1 Preparation


Seal the scaler handpiece (without inserts), scaler nose cone, and light pipe individually and separately in disposable sterilisation bags.

1


Seal the inserts individually inside a disposable bag for sterilisation.

2

Seal the wrench individually inside a disposable bag for sterilisation.

3

8.7.2 Sterilisation Method

The scaler handpiece and other sterilisable accessories are manufactured with materials that resist a maximum temperature of 135°C for a maximum time of 20 minutes.

Once the scaler handpiece and the other sterilisable accessories have been put into bags individually, perform the steam sterilisation process in the autoclave.

① CAUTION: Use sterilisation bags compliant with standard UNI EN ISO 11607-1.

The sterilisation process validated by Mectron S.p.A., in an autoclave with steam, guarantees a SAL 10⁻⁶ by setting the parameters indicated below:

- **Type of cycle:** 3 times Pre-vacuum (pressure min. 60 mBar).
- **Minimum sterilisation temperature:** 132°C (interval 0°C ÷ +3°C).
- **Minimum sterilisation time:** 4 minutes.
- **Minimum drying time:** 20 minutes.

All stages of sterilisation must be carried out by the operator in compliance with the most

current version of standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

① CAUTION: Do not sterilise the scaler handpiece with the insert screwed in.

⚠ WARNING: Infection control - Sterilisable parts. Thoroughly remove all residual organic matter before sterilisation.

① CAUTION: Carry out the sterilisation by using only an autoclave with steam from water. Do not use any other sterilisation procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

① CAUTION: Do not exceed the permissible load of the steam-steriliser.

⚠ WARNING: On completion of the sterilisation cycle in autoclave, the scaler handpiece remains at the sterilisation temperature for a long time. Use appropriate precautions when extracting the scaler handpiece from the autoclave to prevent injury to the operator.

! **CAUTION:** Wait for the scaler handpiece to cool down completely before use.

8.7.2.1 Special Information

Sterilisation parameters, in steam autoclave, used in Great Britain:

- Temperature: 134° C;
- Time: 3 minutes.

EN

9 MAINTENANCE

If the appliance is not used for a long time observe the following recommendations:

1. Run a complete cleaning cycle on the irrigation circuit using the "pump/flush" function (see Chapter 6 on page 27);
2. Empty the circuits of any residual water, removing the tank and making the handpiece work for a few seconds;
3. Disconnect the appliance from the mains;
4. In case of long periods of non-use, place the appliance back in its original

packaging, in a safe place;

5. Before using the appliance again, clean and sterilise the handpiece, inserts, and wrench following the instructions in Chapter 8 on page 31;
6. Check that the inserts are not worn, deformed or broken, with particular attention to the integrity of the tip.


⚠ WARNING: Periodically check the integrity of the power supply cable; if damaged, replace it with an original Mectron spare part.

9.1 Replacing the Peristaltic Pump

On the back of the appliance there is the plastic cover which covers the housing of the peristaltic pump. Remove this protection by pressing on the sides and pulling towards yourself;


! **CAUTION:** Before carrying out operations on the peristaltic pump, make sure the appliance is unplugged and that the liquids container is not connected.

1


Remove the two tubes of the pump from the respective clutches positioned below it;

2


Extract the peristaltic pump from its base, pulling it towards yourself;

3


Attach the new peristaltic pump onto the base until you hear it click in, and connect the two pump tubes to the respective couplings underneath it;

4

Refit the plastic protection onto the peristaltic pump.

5

10 METHODS AND PRECAUTIONS FOR DISPOSAL

EN

⚠ WARNING: Hospital waste. Treat the following objects as hospital waste:

- Inserts, when worn or broken;
- Inserts tightening wrench, when worn or broken.

Disposable materials and materials that carry biological risk must be disposed of according to local regulations in force regarding hospital waste.

⚠ WARNING: When handling the inserts, pay particular attention to the sharp, pointed, irregular parts to avoid any wounds or injuries.

multipiezo must be disposed of and treated as waste subject to separate collection.

Failure to comply with the previous points can result in a penalty pursuant to the directive on waste electrical and electronic equipment (WEEE).

The purchaser is entitled to deliver the device to be disposed of to the dealer who supplies them with new equipment; at Mectron, instructions for proper disposal are available.

11 TECHNICAL DATA

Device compliant with Regulation (EU) 2017/745	Class IIa
Classification under the IEC/EN 60601-1	I Parts applied: type B (insert) IP 20 (appliance) IP 22 (pedal model FS-05)
Essential performance	According to the standard IEC 80601-2-60 the device has no essential performance
Appliance for intermittent operation	55sec. ON - 30sec. OFF with irrigation 30sec. ON - 120sec. OFF without irrigation
Power Supply	100-240 V~ 50/60 Hz
Max. Power Consumption	90 VA
Fuses	Type 5 x 20 mm, T 2AL, 250V
Working Frequency	Automatic Scan From 24 KHz to 36 KHz
Powers	"endo" "perio/scaler" "restorative" "soft mode"

Peristaltic pump flow rate	Adjustable using the touch screen 7 flow rate levels: from 0 (0 ml/min) to 6 (approx. 28 ml/min) (see Chapter 5.2 on page 17) Tank capacity: 500 ml. Tank lighting system: Blue LED light power risk free according to standard IEC/EN 62471.
Handpiece LED system:	Light function set to AUTO: The handpiece LED lights up as soon as the machine starts working, and turns off 3 seconds after the pedal is released Light function set to ON: The handpiece LED is always lit; after 100 seconds of non-use of the pedal, it turns off by itself and the light function switches to AUTO The handpiece LED stays off when both the options (ON and AUTO) are off. White LED light power risk free according to standard IEC/EN 62471
Protections of the APC Circuit	Missing handpiece; Breaking wire cord; Insert not correctly tightened or broken;
Operating Conditions	from 10°C to 35°C Relative humidity from 30% to 75% Air pressure P: 800hPa/1060hPa
Transport and Storage Conditions	from -10 °C to 60 °C Relative humidity from 10% to 90% Air pressure P: 500hPa/1060hPa
Altitude	lower than or equal to 2000 metres
Weights and dimensions	2.4Kg 320 x 230 x 145 mm (L x W x H) ^{a)}

Table 5 – Technical Data

a) W = width ; L = length ; H = height

11.1 Electromagnetic Compatibility IEC/EN 60601-1-2

EN

⚠ WARNING: Contraindications.
Interference with other equipment
 Even if compliant with standard IEC/EN 60601-1-2, multipiezo may interfere with other devices in the vicinity. multipiezo must not be used near or stacked with other appliances. Install multipiezo far from subsistence equipment. However, if this is necessary, the correct operation of the equipment in that configuration, and of all equipment, must be checked and monitored before starting the intervention.

⚠ WARNING: Portable and mobile radio communication equipment may influence the correct operation of the appliance.

⚠ WARNING: Contraindications.
Interference from other equipment
 An electrosurgical scalpel or other electrosurgical units positioned near the multipiezo appliance may interfere with the correct operation of the appliance itself.

⚠ WARNING: The appliance requires particular EMC precautions and must be installed and put into service according to the EMC information provided in this chapter.

⚠ WARNING: The use of cables and accessories not supplied by MECTRON may adversely affect the EMC performances.

11.1.1 Guide and Manufacturer's Declaration - Electromagnetic Emissions

multipiezo is designed to operate in the electromagnetic environment specified below. The purchaser or user of multipiezo should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1	multipiezo uses RF energy only for its internal operation. Therefore, its RF emissions are very low and probably do not cause any interference with nearby electronic devices.
RF Emissions CISPR 11	Class B	multipiezo is suitable for use in all buildings, including domestic buildings, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Emissions of fluctuations voltage/flicker IEC 61000-3-3	Compliant	

11.1.2 Accessible Parts of the Casing

multipiezo is designed to operate in the electromagnetic environment specified below. The purchaser or user of multipiezo should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields ^{a)}	IEC 61000-4-3	3 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80% AM at 1 kHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter.
Magnetic field at the mains frequency ^{d) e)}	IEC 61000-4-8	30 A/m ^{g)} 50 Hz or 60 Hz	Power supply frequency magnetic fields should have levels characteristic of a typical location in a commercial or hospital environment.

- a) If used, the interface between the PATIENT'S physiological signal simulation and multipiezo, must be positioned within a 0.1 m radius of the vertical plane of the uniform field area in the same direction as multipiezo.
- b) multipiezo, which intentionally receives RF electromagnetic energy for its operation, must be tested at the receiving frequency. The test can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test evaluates the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an environmental signal is in the passband. It is understood that the receiver may not receive normally during the test.
- c) The test can be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) Only applicable to equipment and systems with magnetically sensitive components or circuits.
- e) During the tests, multipiezo can be supplied with any NOMINAL input voltage, but with the same frequency as the test signal.
- f) Before applying the modulation.
- g) This test value assumes a minimum distance between multipiezo and the magnetic field sources with a supply frequency of at least 15 cm. If the RISK ANALYSIS indicates that multipiezo will be used less than 15 cm from the sources of the magnetic field with power frequency, the immunity test value should be adjusted according to the minimum expected distance.

11.1.3 Guide and the Manufacturer's Declaration - Electromagnetic Immunity

11.1.3.1 Power Connection BC Input

multipiezo is designed to operate in the electromagnetic environment specified below. The purchaser or user of multipiezo should ensure that it is used in such an environment.

EN

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrical fast transient/burst a) i) o)	IEC 61000-4-4	± 2 kV on contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses Differential mode a) b) j) o)	IEC 61000-4-5	± 0.5 kV, ± 1 kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses common mode a) b) j) k) o)	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conductive disturbances induced by RF fields c) d) o)	IEC 61000-4-6	3 V ^{m)} 0.15 MHz - 80 MHz 6 V ^{m)} in the ISM bands between 0.15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 KHz ^{e)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter.
Voltage dips f) p) r)	IEC 61000-4-11	0% UT; 0,5 cycle g) At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle e 70 % UT; 25/30 cycle ^{h)} Single phase: at 0°	The quality of the network voltage should be that of a typical commercial or hospital environment.
Voltage interruptions f) i) o) r)	IEC 61000-4-11	0% UT; 250/300 cycle ^{h)}	The quality of the network voltage should be that of a typical commercial or hospital environment.

- a) The test can be performed at any supply voltage within the multipiezo range of NOMINAL voltages. If multipiezo is tested at one supply voltage value, it need not be retested at other voltage values.
- b) During the test, all multipiezo cables must be connected.
- c) The calibration of the current injection terminals must be performed in a system at $150\ \Omega$.
- d) If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) The test can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) Devices and systems with a Direct Current (DC) input power supply using AC to DC converters must be tested with a converter that complies with the MANUFACTURER'S specifications. Immunity test levels are applied to the AC power input of the converter.
- g) Only applicable to devices and systems connected to a single-phase Alternating Current (AC) power supply.
- h) For example, 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) Devices and systems with nominal input current above 16 A / phase must be disconnected from the power supply once every 250/300 cycles at any angle and from all phases simultaneously (if applicable). Devices and systems with battery backup, after the test, must resume operation using the power supply line. For devices and systems with a nominal input current lower than 16 A, all phases must be disconnected simultaneously.
- j) Devices and systems that do not have a surge protection device in the primary power circuit can only be tested at $\pm 2\text{ kV}$ between the line(s) and the ground (common mode) and at $\pm 1\text{ kV}$ between line(s) and line(s) (differential mode).
- k) Not applicable to CLASS II devices and systems.
- l) It is necessary to use direct coupling.
- m) R.M.S., applied before modulation.
- n) The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.
- o) Applicable to devices and systems with NOMINAL input current less than or equal to 16 A / phase and devices and systems with NOMINAL input current greater than 16 A / phase.
- p) Applicable to devices and systems with NOMINAL input current less than or equal to 16 A / phase.
- q) At some phase angles, the application of this test to devices with a transformer on the input power supply may cause the opening of an overcurrent protection device. This can occur due to the saturation of the magnetic flow of the transformer core after the voltage drop. If this happens, the device must ensure BASIC SAFETY during and after the test.
- r) For devices and systems that have multiple voltage settings or self-regulating voltage capabilities, the test must be performed at the minimum and at the maximum NOMINAL input voltage. Devices and systems with a NOMINAL input voltage range lower than 25% of the highest NOMINAL input voltage have to be tested with a NOMINAL input voltage within the range.

11.1.3.2 Points of Contact with the Patient

multipiezo is designed to operate in the electromagnetic environment specified below. The purchaser or user of multipiezo should ensure that it is used in such an environment.

EN

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{c)}	IEC 61000-4-2	$\pm 8 \text{ kV}$ on contact $\pm 2 \text{ kV}$, $\pm 4 \text{ kV}$, $\pm 8 \text{ kV}$, $\pm 15 \text{ kV}$ in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Conductive disturbances induced by RF fields ^{a)}	IEC 61000-4-6	3 V ^{b)} 0.15 MHz - 80 MHz 6 V ^{b)} in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 KHz	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

a) The following applies:

- All connection cables with the patient must be tested, either individually or grouped together.
- The connection cables with the patient must be tested using a current clamp unless the current clamp is not suitable. If a current clamp is not suitable, an EM clamp must be used.
- In any case, no intentional decoupling device should be used between the injection site and POINT OF CONNECTION TO THE PATIENT.
- The tests can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- The tubes that are intentionally filled with conductive liquids and intended to be placed in contact with the PATIENT must be considered connection cables with the patient.

• If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

• The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

b) R.M.S., applied before modulation.

c) The discharges must be applied without connection to an artificial hand and without connection to the PATIENT simulation. The PATIENT simulation can be connected after the test, if necessary, to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

11.1.3.3 Parts Accessible to the Input / Output Signals

multipiezo is designed to operate in the electromagnetic environment specified below. The purchaser or user of multipiezo should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{e)}	IEC 61000-4-2	± 8 kV on contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst ^{b) f)}	IEC 61000-4-4	± 1 kV on contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses Common mode ^{a)}	IEC 61000-4-5	± 2 kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conductive disturbances induced by RF fields ^{b) d) g)}	IEC 61000-4-6	3 V ^{h)} 0.15 MHz - 80 MHz 6 V ^{h)} in the ISM bands between 0.15 MHz and 80 MHz ⁱ⁾ 80% AM at 1 KHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

- a) This test applies only to output lines connected directly to the external cables.
- b) SIP/SOPs with maximum cable length less than 3 m are excluded.
- c) The tests can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) The calibration of the current injection terminals must be performed in a system at 150 Ω.
- e) The connectors have to be tested in accordance with paragraph 8.3.2 and Table 4 of the standard IEC 61000-4-2:2008. For insulated connector housings, perform the air discharge test on the connector housing and pins using the probe with the rounded tip of the ESD generator, with the exception that only the connector pins that are tested are those that can be reached or touched, under the conditions of INTENDED USE, by the standard probe shown in Figure 6 of the general standard, applied in a bent or straight position.
- f) It is necessary to use the capacitive coupling.
- g) If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- h) R.M.S., applied before modulation.
- i) The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

11.1.4 Specifications of the Tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device

multipiezo is designed to operate in an electromagnetic environment in which radiated RF disturbances are under control. The purchaser or operator of multipiezo can help prevent electromagnetic interferences by guaranteeing a minimum distance between the mobile and portable RF communication devices (transmitters) and multipiezo, as recommended below, in relation to the maximum output power of the radio communication devices.

EN

Test Freq. (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Max power (W)	Distance (m)	Immunity test value (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460 FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 Band LTE 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2750	Bluetooth WLAN 802.11 b/g/n RFID 2450 Band LTE 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5420	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						

a) For some services, only uplink frequencies are included.

b) The carrier must be modulated using a square wave signal with a duty cycle of 50%.

c) As an alternative to FM modulation, 18 Hz pulse modulation at 50% can be used because, although it does not represent true modulation, it represents the worst case.

NOTE: If necessary to reach the immunity test level, the distance between the transmitter antenna and multipiezo can be reduced to 1 m. The test distance of 1 m is allowed by IEC 61000-4-3.

WARNING: Portable RF communication equipment (including peripheral devices such as antenna cables and external antennas) must not be used closer than 30 cm to any part of the multipiezo device, including the cables specified by the manufacturer. Otherwise, there may be a performance degradation of these appliances.

12 TROUBLESHOOTING

12.1 Diagnostic System and Symbols on the Keyboard

multipiezo is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch keyboard by means of a symbol. Users, by using the following table, are guided to identifying and the possible resolution of the malfunction detected.

Symbol on keyboard	Possible cause	Solution
	Cord electrical contacts wet	Thoroughly dry the contacts with compressed air.
	multipiezo handpiece not connected to the appliance	Connect the handpiece.
	Faulty handpiece	Replace the handpiece
	Tuning circuit malfunction	Contact an Authorised Mectron Service Centre.
	Insert not present or not properly tightened on the handpiece	Unscrew the insert and screw it back in correctly using the torque wrench (See Chapter 5.4 on page 23).
	Insert broken, worn or deformed	Replace the insert.
	Cord electrical contacts wet	Thoroughly dry the contacts with compressed air.

Symbol on keyboard	Possible cause	Solution
	Peristaltic pump malfunction	Check that there are no impediments for the rotation of the pump. Check that the pump and the two pipes are correctly installed.
	The appliance has been switched off and on again without waiting 5 seconds	Turn off and wait for 5 seconds before switching on the appliance again.
	Faults on the electrical network or excessive electrostatic discharges or internal faults	Turn off and wait 5 seconds before switching the appliance on again If the warning persists, contact an Authorised Mectron Service Centre.
	Incorrect start up procedure: the appliance was started with the pedal pressed	Check that the pedal is not pressed down. If the problem persists, disconnect the pedal and potentially contact an Authorised Mectron Service Centre.

Table 6 – Error Messages

NOTE: For diagnostic warnings not included in this list, contact the technical service desk.

12.2 Quick Troubleshooting

Problem	Possible Cause	Solution
The appliance does not start after the switch has been moved to the "I" position.	The terminal of the power supply cable is inserted incorrectly into the rear plug of the appliance	Check that power supply cable is firmly connected
	The power supply cable is faulty	Check that the supply socket is working. Replace the power supply cable
	The fuses are out of order	Replace the fuses (See Chapter 12.3 on page 71)
The appliance is on but is not working. The display reports no errors.	The pedal plug is not correctly inserted into the appliance socket	Correctly insert the pedal plug in the socket on the back of the appliance
	The pedal does not work properly	Contact an Authorised Mectron Service Centre
The appliance is on but is not working. One of the following symbols appears on the screen: 	See Chapter 12.1 on page 67 for the possible cause depending on the symbol	See Chapter 12.1 on page 67 for the action to be taken depending on the symbol
During operation a faint whistling noise can be heard coming from the multipiezo handpiece.	The insert is not correctly tightened on the handpiece	Unscrew the insert and screw it back in correctly using the Mectron torque wrench (See Chapter 5.4 on page 23)
	The irrigation circuit has not been completely filled	Fill the irrigation circuit using the PUMP function (See Chapter 5.4 on page 23)

Problem	Possible Cause	Solution
No liquid flows out from the insert during operation	The insert is the type without a liquid passage	Use an insert type with a liquid passage
	The insert is clogged	Unscrew the insert from the handpiece and release the water passage of the insert by blowing compressed air through it. If the problem persists, replace the insert with a new one
	The handpiece is clogged	Contact an Authorised Mectron Service Centre
	The irrigation level on the screen is set to "0"	Adjust the irrigation level
	The liquid tank is empty	Fill the tank
	The tank is not properly installed	Properly connect the tank to the machine core unit
	The silicone tubes of the pump are not correctly installed	Check the connections of the tubes
Poor performance	The peristaltic pump is worn	Replace the peristaltic pump (See Chapter 9.1 on page 56)
	The insert is not correctly tightened on the handpiece	Unscrew the insert and screw it back in correctly using the Mectron torque wrench (See Chapter 5.4 on page 23)
	Insert broken, worn or deformed	Replace the insert with a new one

Table 7 – Quick Troubleshooting


12.3 Replacing the Fuses

⚠ WARNING: Switch the appliance off.

Always switch the appliance off using the main switch and disconnect it from the power supply socket before carrying out the next intervention.

Pry with a flat screwdriver, inserting the tip into the seat of the fuse holder compartment located under the power socket;


1

Extract the fuse holder compartment;


⚠ WARNING: Replace the fuses in respect of the characteristics indicated in Chapter 8 on page 31.

2

Reinsert the compartment in its place.

3

12.4 Shipping to an Authorised Mectron Service Centre

If technical assistance is required on the machine, contact an Authorised Mectron Service Centre or your Dealer. Do not try to repair or modify the appliance and its accessories.

Clean and sterilise all parts that need to be sent to an Authorised Mectron Service Centre, following the instructions in Chapter 8 on page 31.

Leave the sterilised parts in the bag which certifies the sterilisation process.

The cleaning and sterilisation demands comply with the mandatory requirements for workplace health and safety protection pursuant to Legislative Decree 81/08 and Italian laws.

If the client fails to comply with the requirements, Mectron reserves the right to charge them the cost of cleaning and sterilisation or to refuse goods that arrive in unsuitable conditions, returning them at the client's own expense for proper cleaning and sterilisation.

The appliance must be returned suitably packaged and accompanied by all the accessories and a form including:

- Details of owner and contact number;
- Product name;
- Serial number and/or batch number;
- Reason for return / description of failure;
- Photocopy of the receipt or invoice for the purchase of the appliance.

CAUTION: Packaging

Pack the appliance in its original packaging to prevent damage during transport.

Once the material is received at the Authorised Mectron Service Centre, the qualified technical personnel will evaluate the problem. The repair will be made only upon acceptance by the end client. For further details contact your nearest Authorised Mectron Service Centre or your dealer.

Unauthorised repairs can damage the system and void the guarantee and furthermore will disclaim Mectron from any liability for direct or indirect damage to persons or property.

13 WARRANTY

All Mectron devices, before being marketed, are subjected to a careful final inspection that verifies their full functionality.

Mectron provides a warranty for multipiezo, purchased new from a Mectron dealer or importer, which covers defects in material and workmanship for a period of:

- 2 (TWO) YEARS for the appliance from the date of purchase;
- 1 (ONE) YEAR for the handpiece complete with cord from the date of purchase.

The other accessories are not covered by the warranty.

During the warranty period, Mectron undertakes to repair (or, at its discretion, to replace) the parts of products free of charge, which, according to its judgement, are proven to be defective.

Full replacement of Mectron products is not covered by the warranty.

The manufacturer's warranty and the approval of the appliance are not valid in the following cases:

- The appliance is not used in accordance with the intended use.
- The appliance is not used in accordance with all the instructions and requirements described in this manual.
- The electrical system of the places where the appliance is used does not comply with the laws in force and relative requirements.
- Assembly operations, extensions, adjustments, updates and repairs are carried out by personnel not authorised by Mectron.

- The environmental conditions for preservation and storage of the device do not comply with the requirements indicated in Chapter 8 on page 31.
- Use of non-original Mectron inserts, accessories, and spare parts that may compromise the correct operation of the appliance and cause injury to the patient.
- Accidental breakage during transport.
- Damage due to incorrect use or carelessness, or due to connection to a voltage other than that envisaged.
- Expired warranty.

The expected service life of the device is minimum 5 years.

The service life / duration does not establish a limit of use; the service life of the device defines the period of time after installation and/or commissioning, during which the original performance, or in any case performance suited to the intended use, is guaranteed without there being any degradation such as to compromise functionality and reliability.

The service life is a minimum qualitative objective of the design, therefore, individual parts or components may guarantee superior performance and reliability with respect to that declared by the manufacturer.

The service life assumes compliance with the maintenance schedules set out in this manual, does not include components normally subject to "wear", and is not linked to the warranty period: the service life does not establish any implicit or explicit extension of the warranty period.

CAUTION

The warranty starts from the date of purchase of the device, which evidence is given by the delivery note/purchase invoice issued by the Dealer / Importer or, in case of device with activation code, from the date of activation of the same.

In order to avail of the warranty service, the client must return the appliance, at their own expense, to the MECTRON dealer/importer from which they purchased the product.

The appliance must be returned together with the original packaging, accompanied by all the accessories and by a form containing:

- The data of the owner and telephone number;
- The data of the Dealer / Importer;
- Photocopy of the delivery note/purchase invoice of the appliance by the owner where are reported the date, the name of the appliance and the serial number;
- Description of the failure.

The transport and the damage caused by transport are not covered by the warranty.

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