

TriAuto mini incl. accessories and spare parts

Battery operation

Flexible speed range from 50 – 1,000 rpm

6 memory locations for treatment parameters

Modular extension with Root ZX mini
For all standard NiTi files

Display: LCD

Weight: 78g

Motor: 50 – 1,000 rpm

Battery: 3.7 Volt Lithium-Ion

Dimensions (mm): diameter 28mm, height 150mm



Complete unit consisting of:

- Handpiece
- Angle-piece
- Tabletop battery charger
- Mains cable
- Rechargeable battery
- Instructions for Use

Art. **6970**

No AR oil!

	Angle-piece "TriAuto mini" TR-CM	Art. 6970-001
	Li-Ion rechargeable battery identical to PenCure battery	Art. 6940-820
	Mains cable	Art. 6970-800
	AR spray nozzle identical to DentPortZX	Art. 6905-011
	AR oil spray without nozzle	Art. N5010216 (400ml) Art. N5010218 (6 x 400ml)
	Datalink cable	Art. 6970-006
	File electrode „TriAuto mini“ TR-CM	Art. 6970-004
	Battery cover „TriAuto mini“ TR-CM	Art. N8450188
	Handpiece "TriAuto mini" TR-CM without Akku	Art. 6970-003
	Tabletop battery charger "TriAuto mini" TR-CM	Art. 6970-002