


n. Code	Descrizione / Description	Description / Beschreibung	Descripción
1 99951680	CORPO DISTRIBUTORE DISTRIBUTION BODY	CORPS DISTRIBUTEUR VERTEILERKÖRPER	CUERPO DISTRIBUIDOR
2 99951623	LAMIERA DI CENTRAGGIO CENTERING SHEET METAL	TOLE DE CENTRAGE ZENTRIERBLECH	CHAPA DE CENTRAJE
3 99951624	SUPPORTO LATERALE LATERAL SUPPORT	SUPPOR LATERAL SEITLICHER HALTER	SOPORTE LATERAL
4 51232539	VITE TCEI M3 X 10 INOX TCEI M3 X 10 INOX SCREW	VIS TCEI M3 X 10 INOX TCEI M3 X 10 INOX SCHRAUBE	TORN. TCEI M3 X 10 INOX
5 51262080	VITE T/BOMB EI M5X12 INOX ROUNDED EI M5X12 INOX SCREW	VIS T/BOMB EI M5X12 INOX BOMB EI M5X12 INOX SCHRAUBE	TORN. T/BOMB EI M5X12 INOX
6 97960060	RACC TAPPO M G1/4 CON OR VITON PLUG FITTING M G1/4 W/OR	RACC. BOUCHON M G1/4 AVEC J.T. STOPFENANSCHL M G1/4 M O-RING	UNION TAPON M G1/4 CON OR
7 99951620	CANNOTTO GUIDA GUIDE QUILL	CANOT GUIDE FÜHRUNGSPINOLE	BOTE GUIA
8 97815298	MOLLA DE=34,56 DF=2,16 L0=50,8 DE=34,56 DF=2,16 L0=50,8 SPRIN	RESS.DE=34,56 DF=2,16 L0=50,8 DE=34,56 DF=2,16 L0=50,8 FEDE	MUELLEDE=34,56 DF=2,16 L0=50,8
9 99951621	FLANGIA SEDE MOLLA SPRING SECT. FLANGE	BRIDE RESSORT FLANSCH FEDERSITZ	BRIDA MUELLE
10 51232162	VITE TCEI M4X16 INOX SOCKET HEAD CAP SCR M4X16 S.S.	VIS TCHC M4X16 INOX ZYL.SCHR INN.SECHSK M4X16 INOX	TORN C/CIL ALLEN M4X16 INOX
11 97360059	OR D2,62 D23,47 NBR 70SH OR D2,62 D23,47 NBR 70SH	OR D2,62 D23,47 NBR 72SH OR D2,62 D23,47 NBR 70SH	OR D2,62 D23,47 NBR 72SH
12 99951618	CONVOGLIATORE SUPERIORE UPPER CONVEYOR	CONVOYEUR SUPERIEUR OBENE FOERDERANLAGE	TRANSPORTADOR SUPERIOR
13 55121046	ANELLO SEEGER X ALB D 24 INOX SEEGER RING FOR D=24 SHAFT	ANNEAU SEEGER POUR TIGE D=24 SEEGER RING FUER D=24 WELLE	ANILLO SEEGER PARA ARB D=24
14 97360058	OR D5,34 D47 NBR 70SH OR D5,34 D47 NBR 70SH	OR D5,34 D47 NBR 70SH OR D5,34 D47 NBR 70SH	OR D5,34 D47 NBR 70SH



n. Code	Descrizione / Description	Description / Beschreibung	Descripción
15 99951619	RALLA FISSAGGIO OR OR FIXING SLEWING RING	COURONNE DE DIRECTION OR BEFESTIGUNG SCHEIBE	SLEWING ANILLO
16 51242071	VITE TPSEI M3X8 10.9 TPSEI M3X8 10.9 SCREW	VIS TPSEI M3X8 10.9 TPSEI M3X8 10.9 SCHRAUBE	TORNILLO TPSEI M3X8 10.9
17 99951681	TAPPO SEDE MANIPOLO HANDPIECE SPOT CAP	BOUCHON SIEGE PIECE A MAIN VERSCHL. HANDSTÜCKSITZ	TAPA SEDE PIEZA DE MANO
19 99951628	GHIERA SEDE MANIPOLO HANDPIECE SPOT FERRULE	BAGUE SIEGE PIECE A MAIN NUTMUTTER HANDSTÜCKSITZ	TUERCA SEDE PIEZA DE MANO
20 97200389	ANELLO ADATTATORE D20 MM D20 MM ADAPTER RING	ANNEAU ADAPT. D20 MM D20 MM ADAPTERRING	ANILLO ADAPTADOR D20 MM
21 97440097	ADATTATORE PUNTE ABLATORE SCALER TYP'S ADAPTER	DOUILLE ADAPT. POINTES ZEG SPITZEN ADAPTER	CASQUILLO ADAPT. PUNTAS
22 97200388	ANELLO ADATTATORE D16 TT D16 TT ADAPTER RING	ANNEAU ADAPT. D16 TT D16 TT ADAPTERRING	ANILLO ADAPTADOR D16 TT
23 99951633	FILTRO MICROPOROSO MICROPOROUS FILTER	FILTRE MICROPOREUX MIKROPOROES FILTER	FILTRO MICROPOROSO
24 99951625	SEDE MANIPOLO HANDPIECE SITE	SIEGE PIECE A MAIN HANDSTUECK SITZ	SEDE PIEZA DE MANO
25 97360057	OR D2,62 D20,24 NBR 70SH OR D2,62 D20,24 NBR 70SH	OR D2,62 D20,24 NBR 70SH OR D2,62 D20,24 NBR 70SH	OR D2,62 D20,24 NBR 70SH
26 97150170	UGELLO ROTANTE LECHLER ROTATING NOZZLE	BUSE ROTATIVE ROTIERENDE DUESE	BOQUILLA ROTATIVA LECHLER
27 97920809	SCATOLA ACCESSORIO TETHYS TETHYS TOOLBOX	BOITE ACCESS. THETHYS TETHYS WERKZEUGKASTEN	CAJA ACCES. TETHYS